

AREBOS

Aluminium Rollator

AR-HE-RL42012

User's Manual



Please follow all security measures in this user's manual to ensure a secure use.

Please read and save these instructions. Read through this user's manual carefully before using product. Protect yourself and others by observing all safety information, warnings and cautions. Failure to comply with instructions could result in personal injury and/ or damage to product or property. Please retain instructions for future reference.

Product specifications

AR-HE-RL42012 Aluminum Rollator is designed to assist the user to walk, rest and carry small items of shopping for both indoor and outdoor use on flat surfaces by a single individual with limited mobility with a maximum weight as shown below.

It is important that the user reads and understands the safety advice given within this manual before attempting to operate the product to avoid possible injury to the user or bystanders and avoid damage to the rollator or surrounding property.

SUMMARIZATION

1. INTRODUCTION

The rollators are manufactured from high quality materials, aluminum frame, PP and nylon seat, Polyester bag and locking brake. The handles are adjustable in height and the rollator folds easily for transport or storage.

2. INDICATIONS

AR-HE-RL42012 is designed to assist the user to walk, rest and carry small items of shopping for both indoor and outdoor use on flat surfaces by a single individual with limited mobility with a maximum weight as shown below.

3. EXECUTIVE STANDARD

Type tested according to BS ENISO 11199-2:2005

4. PRODUCT STRUCTURE

Made by hand grips, adjustable parking/running brake, nylon backrest, armrest tube, seat area, height adjustable tube, nylon shopping bag, front wheel set and rear wheels.

APPEARANCE STRUCTURE DRAWING AS BELOW

1. hand grips
2. adjustable parking/running brake,
3. nylon backrest
4. armrest tube
5. seat area
6. height adjustable tube
7. nylon shopping bag
8. front wheel set
9. rear wheel



5. USE

To assist the user to walk, rest and carry small items of shopping for both indoor and outdoor.

6. WORK ENVIRONMENT

Environment temperature: -20°C~60°C

Relative humidity Range: ≤85%

Atmospheric pressure range: 850hPa≤1050hPa

Key performance specifications

Main parts of product	size
Armrest inner width	45cm
Seat height	54cm
Armrest lowest height	78.5cm
Armrest highest height	92cm
Overall width	64.5cm
Overall length	68.5cm
Seat depth	34cm
Folded width	27cm

Shopping bag weight capacity	5.5KG
Rollator weight capacity	136KG
Unit weight	7.6KG
Master carton size	66.0*29.0*37.0 CM

Fittings List

- Product one unit,
- Four knob screws,
- A manual,
- A shopping bag.

Installation & Use

1. Remove all packaging, put all parts out. Insert knob and turn it to adjust the rollator until stable.



Before inserting knob screws in the holes, please make sure all holes are on the same line. Be sure knob screws can be seen through the holes at the opposite side.

2. To open your rollator, push down on either side of the cross bars as shown below. You should feel both sides click into position.

To fold the rollator, gently pull the strap in the middle of the seat as shown below.



Open Rollator



Fold Rollator

3. The rollator's handles can be adjusted in height to suit different users. Simply depress the spring

button and set to the desired height. Repeat for the opposite side and make sure that both handles are secure and set at the same height before use.

After having depressed the spring button and set to the desired height, it should be located in the same position on either side.

(Warning: These holes' diameter is 8-8.2mm, be carefully when using to avoid hurting finger, be sure to keep children away.)



4. String accessories: Put a stick or umbrella in the holder and fix it with the bandage as shown below.



Bandage to fix umbrella or cane

Cleaning and Maintenance

1. Do not use strong corrosive action cleaner to clean the product surface.
2. Do not use acid and alkaline solution or polishing agent to wipe the surface.
3. Product cannot be put in the rain outside
4. Pay attention to safety during assembly: To avoid injury, do not put fingers into the spin position.
5. While resting, hold the break to prevent rollator from sliding.
6. Do not hang anything on the left or right side handle to avoid the rollator becoming unstable and slanted.
7. Please fold before storage.

Transportation and storage

1. AR-HE-RL42012 can be carried by any transport vehicle, make sure to prevent damage by chemical medicines and the erosion.
2. Storage environment: $-40^{\circ}\text{C} \sim 60^{\circ}\text{C}$, relative humidity: $\leq 93\%$. Atmospheric pressure: $500\text{hPa} \sim 1050\text{hPa}$
3. Storage requirement: No corrosive gases and good ventilation.

After Service

When packing or storing this rollator, please adhere to the conditions described in this manual.

Within 6 months of purchase, exchange remains free of charge. 24 months of warranty are granted.

Disposal Instruction

Disposal of the packaging

- Please make reference to the guidelines and standards for appropriate disposal of the packaging valid in your region. In part, the package may consist of plastic bags - watch this respect, with special care to ensure that this is not out of the reach of children. There is a risk of suffocation!

Disposal of waste equipment

- Equipment must be disposed of in accordance with the rules and regulations of the local waste disposal.



Meaning of the "dustbin"

Protect our environment; electrical appliances do not belong in household waste. Use the provided for the disposal of electrical equipment collection points and enter your electrical and electronic equipment that you no longer use. They help ensure that the potential effects of incorrect disposal on the environment and human health to be avoided. So do your part to recycle, recycling and other forms of recovery of waste electrical and electronic equipment. Information on where the devices are disposed of, please contact your local authorities or local Governments.

EU Declaration of Conformity

We,

Canbolat Vertriebs GmbH, Gneisenaustraße 10-11, 97074 Würzburg, Germany,

Hereby declare that the product named below, seen its design and construction as well as according to our sales, has been complied with the relevant and basic health and safety EU-requirements.

Name of the product: Arebos Aluminium Rollator

Model Nr.: AR-HE-RL42012

Art. Nr.: 4260199756577

If the product has any modification not allowed by us, this declaration loses its validity.

Tested acc. to:

BS ENISO 11199-2:2005

Date/ Manufacturer Signature: / Location:

Würzburg, 26.01.2018

Identification of the signatory:



Korhan Canbolat, head of the company

Authorised representative for the technical documentation:

Korhan Canbolat

Office address:

Canbolat Vertriebs GmbH
Gneisenaustraße 10 – 11
D-97074 Würzburg

Warehouse address:

Berenbrinker GmbH
Heiligenhäuschenweg 66
D-33378 Rheda-Wiedenbrück

Umsatzsteuer-Identifikationsnummer: DE 263752326

Gericht der Handelsregistereintragung ist Würzburg, HRB 10082

WEEE-Reg.-Nr. DE 61617071