# ease cup

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### haal

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## EN - INSTRUCTIONS

### haal

We're so pleased you're reading this! Carefully read the instructions before starting to use the EaseCup. If you still have questions after reading this, feel free to contact us using [info@easecup.nl](mailto:info@easecup.nl). We'll happily help!

### The EaseCup - what is so special about it?

- The EaseCup collects your menstrual fluid instead of absorbing it as a tampon, sponge or sanitary pad.
- The EaseCup gives you freedom of movement during your period.
- You can also use the EaseCup when you sleep, do exercise, swim or make love.
- You only need to replace the EaseCup every 12 hours when your menstrual fluid is average to slightly heavy. You can dispose of the EaseCup with your household rubbish.
- EaseCup has the shape of a diaphragm/cervical cap. It is a kind of ring with a strong ring that is made of polypropylene material that is also used in catheters and baby bottles.

Please bear in mind that it may take a little while before you start to get used to using the EaseCup during your period. The first few times, the EaseCup may not feel very comfortable yet or there may be some leakage. You probably experienced initial difficulties too when starting to use tampons while now you don't even think about it.

Most women must first practice a little before the insertion and removal of the EaseCup becomes easy. This is part of the course before starting to use the EaseCup. It's best to practice on a day when you don't have your period or when your discharge is less and when you don't have to work or have other distractions. You'll find it easier to relax on such days and it will be easier to insert the EaseCup and get used to using it.

After practicing for a bit, you will feel how the EaseCup adjusts to your body. The EaseCup is different when compared to tampons, menstrual sponges, monthly pads and menstrual caps that you are used to. The size of the EaseCup is suitable for (nearly) every woman because of its special design. Regardless of your age, the size of your vagina or the quantity of menstrual fluid.

Most women can use an EaseCup longer up to no fewer than 12 hours than traditional period products. The EaseCup does not cause dryness or irritation and has not been linked to toxic shock syndrome (TSS).

### How do you use the EaseCup?

1. The previous figure shows how you squeeze the EaseCup in order to insert it. Wash your hands and try to relax. We recommend sitting on the toilet in a relaxed manner with your knees apart to correctly insert the EaseCup.
2. Remove the EaseCup from its bag and hold it in such a way that the bottom side is hanging downwards. Next, squeeze the edge of the cup in the centre. Down from above, the cup's edge will now look like the figure right.
3. Carefully push the squeezed EaseCup completely into your vagina lengthwise. When you are sitting, your vagina should be nearly horizontal and slightly sloping downwards. Next, push the EaseCup using your finger as far as you can horizontally and backwards. Don't worry: you cannot push it too far and the EaseCup cannot disappear into your body. The EaseCup will automatically move to the correct place of the back of your vagina behind the pubic bone and under the mouth of the womb.

It is important that you do not push the EaseCup upwards like you do with a tampon. The EaseCup would then not be straight under

the mouth of the womb, which can lead to leakage. If you do push it up too far, remove the EaseCup again and reinsert it.

Once inserted, the EaseCup will adjust to your body. The EaseCup therefore offers a unique protection in line with your body that is not only comfortable, but does not leak!

Once you have inserted the EaseCup correctly, you will not know it is there. The EaseCup fits in exactly in the space at the back of your vagina under the mouth of your womb. If you do, however, feel the EaseCup after insertion in your vagina as you do with a tampon, the EaseCup is too much at the front in your vagina. Remove the EaseCup and try again.


We can leave the EaseCup inside for up to 12 hours. That is a lot longer than tampons or sanitary pads.

Every woman is of course different so choose a schedule that suits you with regard to inserting and discarding the EaseCup. Regardless of your personal schedule, do not forget to empty the EaseCup in the toilet every 12 hours. Next, you can dispose of the EaseCup with your household rubbish.

### Replacing the EaseCup

1. First wash your hands. Again sit on the toilet in a relaxed manner with your knees apart to correctly remove the EaseCup.

2. Next, insert a finger in your vagina. Feel where your pubic bone is with your finger at the top side in the direction of your stomach. The edge of the EaseCup can be found exactly behind the pubic bone. Hook your finger behind the edge of the EaseCup and, next, pull the EaseCup carefully towards you. Make sure the EaseCup is as horizontal as possible when you pull it out from your vagina towards the toilet.

Do not pull the EaseCup down when you remove it as you might do with a tampon. Try not to turn the EaseCup either. Simply pull it towards you slowly.

If you cannot hook your finger behind the edge to pull it, do not panic! Try to relax. You can move the EaseCup slightly by using your muscles. Use the muscles that you also use when pushing when going to the toilet or having a baby.

### Warnings:

- Are you using the EaseCup in combination with a cat? Take care: You may accidentally pull out the cat when you remove the EaseCup. If necessary, ask your GP for information before starting to use the EaseCup.
- If you have ever had toxic shock syndrome (TSS), only use the EaseCup after consulting your GP. As far as it is known, there is no risk of TSS when using the EaseCup.
- NEVER USE THE EaseCup FOR MORE THAN 12 HOURS CONSECUTIVELY WITHOUT EMPTYING AND REPLACING IT.
- Do not use the EaseCup immediately after going both or having a miscarriage. Consult your GP about the right moment to start using the EaseCup again. The EaseCup is not suitable as a replacement of maternity pads.
- You do not have to remove the EaseCup when you make love. The EaseCup does NOT, however, offer protection against pregnancy or sexually transmitted diseases (STD).
- If the use of the EaseCup should lead to pain or other complaints with regard to you or your partner, stop using it immediately and, if necessary, contact your GP.
- It is unlikely that the EaseCup will cause an allergic reaction, if your vagina does, however, feel itchy or irritated, stop using it.