## DONUT Pessary

Support for third-degree prolapse/procedentia  
(Stage III + Stage IV)

Rx Only

### AVAILABLE MODELS

|                |                                   |
|----------------|-----------------------------------|
| <b>REF</b>     | <b>MXPD0</b>                      |
| (Pessary Only) |                                   |
| <b>MXKPD0</b>  | (Pessary, TRIMO-SAN™, Jel Jector) |

### Single-Patient Use

CE 0086

Wash pessary with mild soap and thoroughly rinse prior to initial use.

### DESCRIPTION

The Milex® Donut Pessary is a medical device made from silicone that is inserted into the vagina to function as a supportive structure of the uterus, bladder and or rectum.

### WARNINGS

- Do not leave pessaries in place for long periods of time, as serious complications can occur which may require surgical intervention.
- Do not use these pessaries on a patient with a known silicone allergy.
- Chemicals in various vaginal preparations can interact with the pessary material, resulting in discoloration or deterioration of the pessary. TRIMO-SAN™ does not interact with the pessary material.
- Use only TRIMO-SAN. Other materials have not been tested for compatibility.

### INDICATIONS FOR USE

For effective support of third-degree prolapse or procedentia.

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### REPORT ANY OF THE FOLLOWING SYMPTOMS TO YOUR PHYSICIAN

- Any difficulty in urinating
- Any discomfort
- Any changes in the color or consistency of vaginal discharge
- Any increase in the amount of vaginal discharge or vaginal bleeding
- Any foul odor associated with vaginal discharge
- Vaginal itching
- If the pessary falls out

For medical emergencies and for all medically related advice, consult your healthcare professional.

### FOR THE PHYSICIAN / HEALTHCARE PROFESSIONAL

Review these instructions with the patient to establish use regimen.

### INSTRUCTIONS

- Wear dry gloves. When necessary, the entering end of the pessary can be coated with TRIMO-SAN.
- Use one finger to depress the perineum. Hold the DONUT almost parallel to the introitus. Guide the pessary into the vaginal vault using corkscrew motion (see Figure 1).
- The cervix should rest behind the DONUT (see Figure 2).

Figure 1

Figure 2

- Ask the patient to sit, stand and bear down slightly. If there is no leakage, and the patient is comfortable with the pessary in position, have her empty her bladder. A properly fitted pessary takes up slack in redundant tissue, holding the uterus higher in the vagina.
- If the patient can void without difficulty, pessary remains in position upon re-examination, and the patient is comfortable with the pessary in place, this is a good indication that the correct size may have been selected. Patient experience may vary.

**Note:** If the patient is unable to urinate with the pessary in position, remove it and fit her with the next smaller size. Repeat as necessary.

- Examine the patient while she is in the standing position to ensure the pessary has not shifted position. The patient should not feel the pessary once it is in position. The pessary should not be too loose as it may turn or be expelled and it should not be too tight as it may cause discomfort.

**7.** The healthcare professional should be able to sweep one finger between the pessary and the vaginal walls. If there is not enough space to do this, the next smaller size should be tried. If excessive space exists, the pessary will not be effective and may rotate or even be expelled.

- It is sometimes necessary to refit the patient with a different size or type of pessary after a period of time. Do not assume that a replacement will always be the same size as the previous one. Check the fitting to ensure continued patient comfort and relief of symptoms.

The useful life of a pessary is limited. Examine frequently for signs of deterioration (such as cracks or breaks in silicone outer surface). A pessary should be replaced if damaged.

- Ulcerations and erosions frequently occur in cases of complete prolapse due to irritation of the exteriorized cervix or vaginal wall.

Whenever possible, reducing the mass and treating the irritation are primary steps before using a pessary. Prolapse reduction may resolve cervical vaginal irritation. Verification of cervical cytology (pap) or biopsy as clinically indicated.

- During each visit, the vagina should be carefully inspected for evidence of pressure or allergic reaction. The patient should be questioned concerning douching, discharge, disturbance of bowel function or urination. It may be necessary to fit another size or an entirely different type of pessary.

- At the physician's discretion, the patient can be instructed in the proper removal, cleaning and reinsertion techniques for her own pessary. This process can be performed nightly or even weekly by the patient under ideal circumstances.

### TO REMOVE

Hook one finger inside the center of the pessary. Use the thumb and middle finger to compress the side of the DONUT (see Figure 3). Use the finger of other hand to press down on the perineum. Angle the pessary and pull gently through the introitus.

Figure 3

### RECOMMENDED FOLLOW-UP

- Have patient return within 24 hours for first examination.
- Have patient return for second examination within 3 days.
- Have patient return for examination every few months.

**Note:** The schedule of follow-up examinations may be altered to fit the needs of the individual patient at the discretion of the healthcare provider.

### EXPLANATION OF SYMBOLS

|            |                                    |
|------------|------------------------------------|
| <b>REF</b> | Reorder number                     |
| <b>LOT</b> | Batch code                         |
|            | Use-by date                        |
|            | Consult instructions for use       |
|            | Not made with natural rubber latex |

**Rx Only** CAUTION: U.S. Federal law restricts this device to sale by or on the order of a physician

|               |   |
|---------------|---|
| <b>EC REP</b> | Authorized Representative in the European Community |
|               | Manufacturer  |

Dawn® is a registered trademark of Procter & Gamble, Cincinnati, Ohio 45202  
TRIMO-SAN™ is a trademark of CooperSurgical, Inc.  
Milex® and CooperSurgical are registered trademarks of CooperSurgical, Inc.  
© 2017 CooperSurgical, Inc.

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