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EN- WonderMom offers women in their post-natal period stylish support with the 3-in-1 breast compresses. These reusable and toxin-free compresses are suitable for cold and hot use, and fit into your bra.

3-in-1 breast therapy pearl compresses

 To soothe pain and swelling in the breasts in case of: engorgement / general breastfeeding-related discomfort / menstrual pain / pain after a breast operation.

 Before feeding the baby, in case of: breast inflammation / blocked milk ducts / necessity to stimulate the milk flow and reflex.

 Fits every breastmilk pump. The milk flow will be stimulated, so that pumping will take less time.

User instructions:

 Place the compresses without covers (in a clean freezer bag) in the freezer at least 2 hours before use, or submerge the compresses in ice cold water until they have the desired temperature.

 Heat the coverless compresses. In warm (not boiling!) water. Put the compresses in warm water for 3-5 minutes with the heat source off, or put them in the microwave. Then keep heating for 5 seconds at a time, in order to reach the desired temperature. Never heat for longer than 30 seconds. Put the compresses in their covers when they have the right temperature.

 700 watt: 20 sec  1000 watt: 15 sec  1200 watt: 10 sec

Please note: test the temperature before use. The compresses may be used for a maximum of 20 minutes.

Washing: you can wash the compresses with water and soap. Wash the covers at max. 40 degrees.

Storage: store in a cool, dry place, out of direct sunlight.

Warning: For your own safety, always follow these instructions, and save the instructions. Always read the instructions before use. Do not put heavy objects on top of the compresses. Keep out of reach of children and animals. Only use with adult supervision. PLEASE NOTE: Do not exceed the recommended heating time. Always remove the grille from the microwave and use the rotating glass plate, never put the compresses on a metal surface. Always consult the instruction manual of your microwave. Overheating could damage the compresses and the surrounding area. Check for damage before use. Never apply to open wounds or sensitive skin. Never puncture the compresses and/or do not eat the contents. Do not fall asleep while using the compresses. Please contact a doctor in case of extreme, long-lasting or chronic pain. Only for external use. WonderMom is not accountable for any damage as a result of improper use or non-compliance with the instructions for the compresses.

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