Smart Body Composition Scale

 Manual / Handleiding

SACCOMA

Specifications

Size

278x278x27.5mm

TFT-Display size

2,8 inch

Power

Rechargeable lithium cell 1000mAh

Max. capacity

5 - 180kg

Divison

0.05kg

Charging time

2 hours

Connection

Wi-Fi and Bluetooth

Duration of use

160 days with regular use

Inhoud

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Foreword

This scale is now your personal health assistant. In order to guarantee an accurate measurement and a long service life of the scale, please read this manual before using and keep it for further reference.

This Smart Body Composition Scale could help you keep track of your personal health data. It records and analyzes the information from your scale, then transfer the data to your mobile device via wireless transmission technology. You can view your health data at any time of the day by downloading the free app.

 @SACOMA.NL

Thank you for your order!

Thank you for choosing our Smart Body Composition Scale! Your opinion matters. If you love our Composition Scale as much as we do, we would love to hear your feedback in a review.

Problems with your order?


The Smart Body Composition Scale is finished with care and is delivered to your home in perfect condition. If this is not the case, you can contact our customer service via: info@sacoma.nl
(Response within 12 hours on business days)

Discover more products online

Sacoma |

Display structure

Display description

C

error; step off
the scale


Batteries
almost empty

OVERLOAD

Remove the object
that is too heavy

1. User
2. Time
3. Bluetooth
4. WiFi
5. Battery
6. Unit
7. Weight
8. Measure functions
9. Data last 15
measurements

Icons display

Weight

BMI

Heart rate

Bone mass

Muscle
mass

Body water
rate

Body fat
rate

BMI tabel

Suitable for men and women.

Age	Under-weight	Healthy weight	Over-weight	Obesity	Severe obesity
18-100	<18.5	18.5-23.9	24-27.9	28-34.9	>35


Body Mass Index (BMI) = Weight (kg) : Height ² (m)

For example: BMI = 70 kg : (1.75 x 1.75m) = 22.86 (healthy weight)

Note: The BMI index is a commonly used standard to measure body weight and health status.

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Place & use


Step 1

Press the reset / power button to turn the scale on.


Step 2

Charge the scale for 2 hours with the supplied charging cable.


Step 3

Always place the scale on a hard, dry and flat surface. As soon as the TFT screen shows 0.0 kg, you can step on the scale.

Step 4

The lower part of the scale provides guidance. Therefore, make sure that your feet are always symmetrical on the scale.

Measurements

Step 1

gently press the scale to wake it up. When the scale shows "0.00", step on the scale with bare feet and make sure that your feet are in full contact with the electrodes.

note; please take of your shoes en socks, otherwise the scale can only measure your weight.

Step 2

Stand still while measuring. The weight value will flash three times and lock. Then the body fat is measured.

Note; if the display screen does not show any readings, please try to wet the soles of your feet and measure again.

Step 3

Wait until the scale vibrates and displays the measurement results and change curves of weight, BMI, body fat rate, body water rate, muscle mass, bone mass and heart rate in turn, then the measurement is complete.

note; do not step off the scale until the measurement is complete.

Step 4

the scale will automatically shut down and you can view multiple types of body data in the app.

App download

After you downloaded the App, do not forget to turn on your Bluetooth or Wi-Fi.

Download "Fitdays"

In the App Store,
Google Play or scan the
QR code.

App instructions

Create an account

1. Turn on Bluetooth or Wi-Fi and open the Fitdays app.
2. Click on "Register".
3. Go through all the fields and click on "Next".
4. Set your personal data and click on "Confirm".
5. Stand on the scale with your bare feet, the weighing will start automatically.
6. Wait until the weighing is complete. Your health data will be visible in the App.



Pairing

Koppelen via WiFi

1. Tap the scale icon to enter WiFi setting page.
2. Choose a stable 2,4 GHz WiFi and enter the correct WiFi password.
3. When “√” pops up in the app, you have successfully paired your scale with WiFi.

Note; Each user needs tot ap “account” > “Add user” to complete a personal profile, because the data you enter (gender, date of birth, height) is necessary to provide you with the most accurate information.

User data

1. Charts


Click on the “Charts” button to view and compare your data per day, week, month or year. This way you can see at a glance whether you have made progress or not.

2. User data

Then click on the “User data” button. Here you can click on any date you like and view and compare your health data.

3. Contrast

Click “+” in the top right corner and click on “Contrast”. Then click on the two desired weightings and click on  to start comparing the measurement data.

4. Share data

Click “+” in the top right corner and click on “Export”. Choose whether you want to export all data or only for the past month. Then click on  to share the information with your friends via Facebook or Instagram, among others.

Account

1. Add user


Click on the “Account” button and then click on “Add User”. Then go through all the steps and click on “Confirm”. You have now added a user.

2. Delete user

Click on the “Account” button. To delete a user, swipe the user to the left and select “Delete”.

3. Setting


From your account, click on “Settings”. Here you can set goals, change units, change your password and connect your Fitbit.

Sync with fitness apps




Apple Health

1. Open the 'Health' app on your iPhone and click on your account in the top right corner.
2. Click on "Apps" and then click on "Fitdays".
3. Turn on all categories to make the Fitdays app work optimally.


Google Fit

1. Go to settings and select "Google Fit".
2. Turn on the "Google Fit" button, then log in with your Google account.
3. You will be automatically redirected to the "Google Fit" page.


Fitbitt App

1. Go to settings and select "Fitbit".
2. Turn on the "Fitbit" button, then log in with your Fitbit account.
3. Select "Allow" to access the Fitbit app.
4. After pairing, the data is synced to the Fitbit app. You can view the data via the app as well as on the website.

Samsung Health app

1. Go to settings and select “Samsung Health”.
2. Turn on the “Samsung Health” button, then log in with your Samsung account.
3. You will be automatically redirected to the “Samsung Health” page.

Calibrate

Have you moved or turned the scales? Then we recommend that you calibrate the scale to ensure the accuracy of the measurements.

1. Place the scale on an even and dry surface.
2. Place one foot on the scale until the TFT screen responds.
3. The calibration process is complete when the TFT display shows "0.00kg".

Recognition function


After the first weighing, the scale has an automatic recognition function. When the weight deviates more than $\pm 2\text{KG}$ from the previous measurement, the app asks if the correct user has been selected. An envelope icon with a notification will appear on the app home page.

1. Claim weight

Click on the notification and click on "Claim". You can now link the weight measurement to the correct account.

2. Remove weight

Is the measurement incorrect? Then click on "Select", then check the incorrect measurement and click on "Delete".

FAQ

Why doesn't the scale work? Why does the data on screen disappear in a flash?

First check that the batteries are not empty. Are the batteries charged, place one feet on the scale for approx. 10 seconds. Contact us if the scale still does not switch on after completing these steps.

The scale cannot connect with App.

1. Check that the phone's software is iOS 8.0 or higher or Android6.0 or higher.
2. Download and open the newest version of the Fitdays App.
3. Enable Bluetooth on your phone. For Android 6.0 or higher, location needs to be enabled.

No body fat measured when weighing.

1. Step on with dry and bare feet.
2. Ensure Bluetooth or Wi-Fi is enabled and working.
3. Ensure personal data is entered.

Attention

Safety precautions

When using digital scale, certain basic safety precautions must be followed, including the following:

1. The Bluetooth scale should not be used to diagnose or treat any medical condition. You should consult a doctor when undertaking any diet or exercise program.
2. If you are pregnant or have a pacemaker or other internal device, you should not over rely on the data provided by this product. All data is for reference and is not a substitute for advice from a licensed health care professional.
3. Check the device before each use. Do not use the device if damaged. The continuous use of a damaged unit may cause injury or improper results.

Note

1. Scale is a high-precision measuring device. Never jump or stamp on the scale or disassemble. Please handle the scale carefully to avoid breakage.
2. Place the scale on a hard, flat surface during use. Using the scale on a soft surface like carpet or linoleum may affect the accuracy of the scale readings.
3. To ensure weighing accuracy, please gently step on the scale to wake it up and wait at least 10 seconds before weighing on the scale. You must also wake up the scale if it has been moved. Take your measurements at the same time each day to ensure the most accurate results.
4. Before using the scale, make sure you have entered all your personal data through the Fitdays app.
5. Use the scale while barefoot. The scale cannot measure body fat percentage if you're wearing shoes or socks.
6. Make sure your feet are dry before weighing yourself.
7. This scale is not waterproof; never immerse the scale in water. To clean the surface, use a damp cloth or glass cleaner. Do not use soap or other chemicals.

This device complies with part 15 of the FCC Rules. Operation ; subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation. Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation.

If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
 - Increase the separation between the equipment and receiver.
 - Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
 - Consult the dealer or an experienced radio/TV technician for help.
- The device has been evaluated to meet general RF exposure requirement.
- The device can be used in portable exposure condition without restriction.

FCC-ID:2ANDX-CK10A

100% satisfaction guarantee

We guarantee all our products are designed to the highest standard, we base our ethos on customer satisfaction. This means the customer experience and journey you have had with us should be exceptional. If you have experienced our customer service or product, please let us know by leaving us an honest review.

INFO@SACOMA.NL

COPYRIGHT © 2021 SACOMA. ALL RIGHTS RESERVED.

The material in this publication is protected under International and Federal Copyright Laws and Treaties, and as such, any unauthorized reprint or use of this material is strictly prohibited. No part of this book may be reproduced or transmitted in any form without the written permission of the author, except for the inclusion of brief quotations in a review. Reproduction or translation of any part of this work without the permission of the copyright holder is against the law.