Mondtape & Neusstrips Handleiding/Manual/Manuel

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Contents: With 33 pieces in each package, PurelyGoods provides enough mouth tape and nose strips for approximately one month of snore-free nights and a fresh mouth upon waking.

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klantenservice@purelygoods.nl

How to use?

- Clean the skin around your lips thoroughly with water or a mild cleanser and dry it well.
- Carefully remove the mouth tape from the plastic sheet without touching the adhesive side.
- Place the tape over your closed lips for a good seal and comfortable fit.
- Press the edges of the tape firmly for 10 seconds to ensure it stays in place throughout the night.
- Remove the tape the next morning by gently pressing against it with your tongue and slowly peeling it off.

Warnings:

- Do not use the mouth tape for more than 12 hours consecutively.
- Do not use on irritated, burned, or damaged skin around the mouth.
- Mouth tape is intended for adults, not for use in children.
- Discontinue use if a rash or irritation occurs.
- Each mouth tape is for single-use only.
- Not suitable for children under 5 years of age.

PurelyGoods® Mouth Tape is not intended for the treatment of sleep apnea.

How to use?

- Clean your nose thoroughly with alcohol or a mild soap and dry it completely to ensure optimal adhesion.
- Carefully remove the backing of the nose strip without touching the adhesive side.
- Place the strip across the bridge of your nose, just above the nostrils, for proper positioning and effective performance.
- Press the edges of the strip firmly for 30 seconds to ensure it stays in place throughout the night.
- Remove the strip the next morning by gently lifting one end and slowly peeling it off evenly from your nose.

Warnings

- Do not wear the nose strip for more than 12 hours.
- Do not use on nasal sores and/or burned or irritated skin.
- Nose strips are intended for adults, not for use in children.
- Discontinue use immediately if a rash or irritation occurs.
- Each nose strip is for single-use only.
- Not suitable for children under 5 years of age.

PurelyGoods® Nose Strips are not intended for the treatment of sleep apnea.

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