



## bEbird M9 S Smart Ear Cleaner User Manual

[Home](#) » [bebird](#) » bEbird M9 S Smart Ear Cleaner User Manual 

# bebird M9 S

## User Manual

### Contents

- [1 M9 S Smart Ear Cleaner](#)
- [2 Product Overview](#)
- [3 How to Replace Ear Spoon](#)
- [4 Accessories](#)
- [5 Hazardous Substance](#)
- [6 Troubleshooting](#)
- [7 Attention](#)
- [8 Product Parameters](#)
- [9 Product Parameters](#)
- [10 Documents / Resources](#)
- [11 Related Posts](#)

### M9 S Smart Ear Cleaner



Note: Do not use the charging equipment with output voltage exceeding 5V

How to connect the App

**Note:** Do not use the device unless you have successfully connected the device to your mobile phone through the Android/iOS App and are streaming live video footage

1. Download the APP

iDS: Search for “bebird” in App Store Android: Search for “bebird” in Google Play



OR scan the QR code to download the app



2. Open the cover Power on automatically




Open the cover  
Power on  
automatically




Green light  
flashes


3. Open the APP Connect device WiFi




4. Connect Product Wi-Fi "bebird-M9-XXXXXX"
5. Back to be bird APP Click on "GO " to use it



6. The visual screen



## Product Overview



## How to Replace Ear Spoon



If the ear spoon is damaged, please stop using it immediately. Please make sure the ear spoon is secured on the device before use.



Screw off the ear spoon and replace with a new one.

Please replace the ear tip with a new one every 1-2 months to avoid aging and breakage.

### Accessories

1

Ear spoon  
x 9



2

Charging  
Cable x 1



3

Standing  
base x 1



4

Cotton swab  
(for lens cleaning)  
x 3



### Hazardous Substance

Part name	Hazardous Material					
	(Pb)	(Hg)	(Cd)	(Cr(VI))	(PBB)	(PBDE)
Product Components	X	0	0			
Case	0	0	0	0	0	0
Plastic Libbe accessory		00		0	0	0

This form is formulated in accordance with SJIT 11364

O: represent the hazardous material ratio is below GB126572-2011

X: represent the hazardous material ratio is over GB126572-20

## Troubleshooting

Problem	Solution
No response after startup	Please use the incidental charging cable to charge
The LED light of the camera is not on	Interrupt the connection, restart the product to connect again
The app exits quickly	Delete the existing app and scan QR to reinstall it
The indicator light blinks quickly	Please use the incidental charging cable to charge
Abnormal display in app-	Interrupt the connection, restart the product to connect again
The Wifi connection was unstable and disconnected in a few seconds	Interrupt the connection, restart the product to connect again
The visual interface blur	Please wipe it with a double-ended cleaning brush or alcohol swab
Disorientation when in use	Please switch the left / right ear mode under the APP interface
If you still have problems with the above procedures, please contact our Customer Service Representatives directly for further assistance.	

## Attention

1. When cleaning the Lens, wipe It carefully with an alcohol cotton swab.
2. Do not use the device while you are moving. Use in a safe place to avoid being bumped by others.
3. This product is unsuitable for children under the age of 3
4. Children should not use this device unsupervised by an adult, to avoid accidental harm.
5. Do not drop or submerge the device in liquid, to avoid damage.
6. This product has a built-in lithium battery. not used frequently, it is recommended to charge the battery once a month
7. Do not expose the device to direct sunshine. Also keep the accessories out of direct sunshine to avoid softening.
8. During use, the device may warm up slightly (up to 35 degrees). This is normal and you can continue to use the device.
9. When using the device on a child, for the best experience, keep the child in a sitting position and use the device horizontally.

## **Product Parameters**

### **Attention**

1. When cleaning the Lens, wipe It carefully with an alcohol cotton swab.
2. Do not use the device while you are moving. Use in a safe place to avoid being bumped by others.
3. This product is unsuitable for children under the age of 3
4. Children should not use this device unsupervised by an adult, to avoid accidental harm.
5. Do not drop or submerge the device in liquid, to avoid damage.
6. This product has a built-in lithium battery. not used frequently, it is recommended to charge the battery once a month
7. Do not expose the device to direct sunshine. Also keep the accessories out of direct sunshine to avoid softening.
8. During use, the device may warm up slightly (up to 35 degrees). This is normal and you can continue to use the device.
9. When using the device on a child, for the best experience, keep the child in a sitting position and use the device horizontally.

## **Product Parameters**

### **FCC Warning**

This device complies with part 15 of the FCC rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

**NOTE:** This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following



measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

### Radiation Exposure Statement

The device has been evaluated to meet general RF exposure requirement. The device can be used in portable exposure condition without restriction.


<b>Brand</b>	<b>bebird</b>
Product Name	Smart Ear Cleaner
Model	M9 S
Network standard	IEEE802.11b/g/n
Antenna	Internal FPC antenna
Working frequency	2.4Ghz
Image transmission rate Image sensor	20fps
Working temperature	CMOS
Battery	#NAME?
Battery life	about 90 minutes
Charging time	about 60 minutes
Input Parameters	5V = 0.5A
Lens diameter	3.5mm
Optimum focal length	13.5mm
Pixel	10 Mega-pixels
Gyroscope	Gimbaled Intelligent Gyroscope



Warranty  
User Info  
User ,me  
Add  
Telephone number  
Seller'

Telephone number.

## Documents / Resources

	<p><a href="#">bEbird M9 S Smart Ear Cleaner</a> [pdf] User Manual 2AU63-M9PRO, 2AU63M9PRO, M9 S, Smart Ear Cleaner, M9 S Smart Ear Cleaner, Ear Cleaner</p>
---	--

[Manuals+](#)