COMPRESSOR NEBULIZER LD. Instruction Manual.

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
fig.3 rys.3 3. ábra fig.3 рис.3 Depending on the type of baffle, the particles of different sizes are distributed aerosols as follows:

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EN PARTS AND COMPONENTS (fig.1)

- COMPRESSOR.
- POWER SWITCH.
- ANGULAR HOLDER.
- CONNECTOR OF COMPRESSOR.
- INTAKE.
- SOCKET FOR THE FILTER.
- NEBULIZER LD-N10S.
- UPPER PART OF NEBULIZER.
- BAFFLE "B" (BLUE) LD-N002.
- NOZZLE.
- BOTTOM OF NEBULIZER.
- CONNECTOR OF NEBULIZER.
- ADULT MASK LD-N041.
- CHILD MASK LD-N040.
- INHALATION TUBE LD-N0S1.
- INHALATION MOUTHPIECE.
- ADULT NASAL NOZZLE LD-N0S8.
- CHILD NASAL NOZZLE LD-N0S9.
- BAFFLE «A» (YELLOW) LD-N001.
- BAFFLE "B" (BLUE) LD-N002.
- BAFFLE «C» (RED) LD-N003.
- FILTER INHALED.
- BAG.

NEBULIZER THERAPY – WHAT IS IT? Nebulizer is a device for formation and spraying of aerosol. The word "nebulizer" is derived from the Latin word "nebula" (fog, cloud) and was first used in 1874 for a device that turns a liquid substance into an aerosol for medical purposes. One of the first portable "aerosol apparatuses" was created by J. Sales-Girons in Paris in 1859. The first nebulizers were used as steam jet energy sources and were applied for inhalation the vapors of resins and antiseptics by tuberculosis patients. Presently, the term "inhaler" is often used instead of "nebulizer". The purpose of the nebulizer therapy is to quickly deliver to the respiratory passages a therapeutic dose of a preparation in aerosol form. Continuous supply of aerosol allows, within several minutes, creating high concentration of a medicine in the upper and lower respiratory passages and lungs, with low probability of any by-effects. Respectively, effective bronchodilation (bronchi expansion) is reached, and the need for hospitalization is eliminated or the hospital stay is reduced.

Little Doctor International (S) Pte. Ltd. offers you to use inhaler LD-220C, whose distinctive features are the possibility to use a wide range of medicines, low inhalation solution residual volume, and reliable and simple use. We thank you for your choice.

GENERAL INFORMATION Compressor nebulizer LD is designed for treating the diseases of respiratory passages and lungs by medicine solution aerosols.

This Instruction Manual is designed to assist the user with safe and effective operation of the Compressor Nebulizer LD.

Use this Device according to the rules described in this Manual. Operate the Device only as intended. Do not use the Device for any other purposes. Read and understand the whole Instruction Manual.

Functionally, the device consists of an air compressor and nebulizer (aerosol formation chamber). The air compressor, on/off power switch and air filter are united in one casing. From the air compressor, the compressed air is fed through a pipe to the nebulizer, where aerosol is formed. For cooling the compressor, air is force-feed into the device's casing.

SAFETY INFORMATION To assure the correct use of the product, basic safety measures should always be followed including the warnings and cautions listed in this instruction manual.

WARNING For regime of medication shall follow the instructions of your physician or licensed healthcare practitioner. Do not cover the compressor with a blanket, towel, or any other type of cover during using. This could result in the compressor overheating or malfunctioning.

Do not use the device where the device may be exposed to flammable gas or vapors. Do not use mineral water in the nebulizer for nebulizing purposes.

Always dispose of any remaining medication in the medication tank after each use. Use fresh medication each time you use the device.

Do not leave the device or its parts where it will be exposed to extreme temperatures or changes in humidity, such as leaving the device in a vehicle during warm or hot months, or where it will be exposed to direct sunlight.

CAUTION Limit the use of the device to 20 minutes at a time, and wait 40 minutes before using the device again.

Provide close supervision when this device is used by, on, or near infants, children or compromised individuals. Do not insert any object into the compressor.

Make sure that the nebulizer kit is correctly assembled, the air filter is properly installed, and the air tube is correctly connected to the compressor and the nebulizer kit. Air may leak from the air tube during use if not securely connected.

Inspect the compressor (main unit) and the nebulizer parts each time before using the device. Make sure no parts are damaged, the nozzle and air tube are not blocked and the compressor operates normally.

Do not use the device if the air tube is bent. Do not block the air filter cover. Do not alter the baffle, the nozzle in the medication tank or any part of the nebulizer kit.

Do not add more than 10ml of medication to the medication tank. Do not operate the device at temperatures greater than 40°C. Do not tilt the nebulizer kit so the angle of the kit is greater than 45°. Medication may flow into the mouth.

Do not shake the nebulizer kit while using the device. Do not stop the compressor, or any of the components to strong shocks, such as dropping on the floor. Do not use the compressor, or any of the components.

This device is approved for human use only. Do not disassemble or attempt to repair the device or components. Use the device only for its intended use as described in the instruction manual. Do not use attachments not recommended by the manufacturer.

Dispose of the device, components and optional accessories according to applicable local regulations. Unlawful disposal may cause environmental pollution.

Make sure that the air tube is securely attached to the compressor (main unit) and nebulizing parts, and does not come loose. Twist the air tube slightly when inserting it into the connectors to avoid the tube disconnecting during use.

RISK OF ELECTRICAL SHOCK Do not use the compressor (main unit) and the power cord while they are wet. Do not plug or unplug the power cord into the electrical outlet with wet hands.

Do not immerse the compressor (main unit) in water or other liquid. Do not spill water or other liquids on the compressor. These parts are not waterproof. If liquid spills on these parts, please unplug the power cord and wipe off the liquid with gauze or other soft absorbent material immediately.

Do not use or store the device in humid locations or outdoors. Use the device within the operating temperature and humidity. Do not overload power outlets. Plug the device into the appropriate voltage outlet.

CARE AND MAINTENANCE Regularly clean the device and all accessories. It is recommended that all the accessories should be wiped with a 3% solution of hydrogen peroxide with addition of 0.5% solution of a detergent (for example, a laundry powder). After that, the nebulizer should be washed by a rich jet of water. The mouthpieces and nose nozzles may be treated by boiling for 10 minutes or autoclaving at 150°C. After the treatment, wipe dry all parts of the device with a soft cloth. Regularly check the filter for dirt and replace it when needed. FILTER REPLACEMENT IS RECOMMENDED AT LEAST ONCE A YEAR. Repair the Device only in authorized organizations. On expiration of the warranted service life apply from time to time to authorized repair organizations to check the technical condition of the Device. Dispose of time to authorized repair organizations to check the technical condition of the Device. Dispose of the Device and its components according to the application local regulations. No special requirements to disposal of this Device are defined by the manufacturer.

USING THE DEVICE Preparing for inhalation. IMPORTANT: Before using the appliance for the first time it is necessary to make a full cleaning, as described in last paragraph «SAFETY INFORMATION».

- Place the nebuliser in front of you on the table. Make sure the device is turned off (power switch is in position «0»), and the power cord is not plugged into the mains.
- Remove the top of the nebulizer by turning it counterclockwise (Fig.2.a).
- Set the desired baffle.

Factory installed baffle inside the nebulizer is baffle «B» (blue), which is effective to affect the entire respiratory tract. For a more effective impact medicines on upper respiratory tract, set, instead of the blue baffle, red baffle «C». For a more effective impact on the lower respiratory tract - baffle «A» yellow color, which consists of two parts (Fig. 2.b, figures indicate the order of assembly). Graphics of the differential particle size distribution by mass for different nozzles are shown in Fig. 3. For effective delivery medicines to a specific site of respiratory tract, need to use the appropriate baffle.

4. Fill the bottom of the nebulizer inhalation solution. The dosage should not exceed the recommended by your doctor. The number of nebuliser solution is determined by the scratches on the case. The maximum reservoir volume of 10 ml.

5. Attach the nebuliser at the top, turning it clockwise until it stops. 6. Depending on the type of inhalation, using either a mouthpiece or nozzle or mask.

Hold the nebulizer vertically. IMPORTANT: Each patient is encouraged to use personal mouthpiece, a mask and / or nozzle for the nose. Plug the power cord to an electrical outlet.

CONNECT one end of the inhalation tube to the fitting of compressor, and others - to the fitting of nebulizer. Turn on the nebulizer, switching the power switch in position «I». NEBULIZER IS READY FOR INHALATION.

PERFORMING THE INHALATION The length of one treatment session should not exceed 20 minutes. Consult your attending physician about the length of the inhalation procedure.

You should always be calm and relaxed during the inhalation. Breathing should be slow and deep, so that the preparation could fill the lungs well and reach the deep portions of bronchi. Briefly hold your breath, and then exhale slowly. Do not attempt to breathe too rapidly. Make pauses if you feel that you need it.

BREATH-ACTUATED NEBULIZER (fig.4) The special design of the nebulizer in the form of chambers connected in a certain manner provides different ways of air streams during inhaling and exhaling. It allows obtaining the air stream with greatest aerosol concentration when inhaling and reducing aerosol loss when exhaling. The effectiveness of inhalation using the breath-actuated nebulizer is increased significantly.

COMPLETING THE INHALATION When the inhalation solution is used up and the inhalation time recommended by the doctor has expired, turn the device off by putting the tumblers in «0» position and unplugging it. After inhalation, breathe fresh air for some time for better treatment effect. After each application of the device, the residual preparation should be removed out of it. Clean and wash the device as described in last paragraph «SAFETY INFORMATION».

TECHNICAL SPECIFICATIONS	
Model	LD-220C
Type	Compressor Nebulizer
Rating	AC 230V, 50Hz
Power consumption	190VA
Extreme Pressure	210kPa–400kPa
Free Flow Range	≥7L/min
Nebulizing Pressure	60kPa–180kPa
Noise Level, approx.	~55 dB*
Max capacity of nebulizer kit	10ml
Nebulization rate, approximately	
Inhalation baffle «A» LD-N001	0.3 ml/min.*
Inhalation baffle «B» LD-N002	0.4 ml/min.*
Inhalation baffle «C» LD-N003	0.5 ml/min.*
Particle size (MMAD)	
Inhalation baffle «A» LD-N001	3.5 μm.*
Inhalation baffle «B» LD-N002	4.0 μm.*
Inhalation baffle «C» LD-N003	5.0 μm.*
Max. operating time	20 minutes
Operation mode	20 minutes on, 40 minutes off
Operation temperature, humidity/atmospheric pressure	10°C – 35°C 15%– 80% 86–106 kPa
Transport and storage temperature and humidity/atmospheric pressure	-20°C – -40°C 15%– 95% 50–106 kPa
Size (compressor)	147 x 95 x 137 mm
Net Weight (compressor)	~1300 g
Pollution Degrees	Degrees 2
Overvoltage Category	Category II
High Altitudes (m)	2000m
Year and month the manufacture is given in the bottom of the unit body in a serial number after symbol "A"	Year and month the manufacture is given in the bottom of the unit body in a serial number after symbol "A"

SYMBOL EXPLANATION: CE – Compliance with the Directive 93/42/EEC; E – Protection class II; I – Important: Read the instruction; R – BF type product; S – Storage/operation temperature, humidity and atmospheric pressure limitation; M – Manufacturer in the EU; SN – Serial number; K – Keep from moisture

Revision date of the present Manual is indicated on the last page as EXXX/YMM/NN, where YY is the year, MM is the month and NN is the number of revision. Technical characteristics may be changed without preliminary notification to improve the operation and quality of the product.

TROUBLESHOOTING TIPS No power on device when the power switch is on: Turn the power switch off. Plug the power plug into an electrical outlet. Turn the nebulizer on. No vibration or low nebulization rate when the power is on: Add the correct amount of prescribed medication to the medication cup.

Make sure the nebulizer kit is correctly assembled and the inhalation accessory is correctly attached. Hold the nebulizer kit correctly. Do not tilt the nebulizer kit so the angle of the kit is greater than 45 degrees. Make sure the air tube is correctly attached to the compressor and the nebulizer kit.

Make sure the air tube is not folded, kinked or bent. Inspect the air tube for any damage. Replace the air tube if damaged. The device is very hot: Do not cover the compressor with any type of cover during use. Turn the device off. Wait 40 minutes before using the device again.

WARRANTY The warranty liabilities are contained in the warranty card given at the sale of this device to a purchaser. The warranty does not apply to the consumables (masks, mouthpieces, inhalation tubes etc.).

CERTIFICATION AND STATE REGISTRATION This device manufacturing is certified according to international standard ISO 13485. IEC 10601-1, EN 10601-1, EN 10601-1-2, ISO 14971.

Manufactured under control and for Little Doctor International (S) Pte. Ltd. (7500A BEACH ROAD #11-313 THE PLAZA SINGAPORE 199591. Postal address: Yishun Central P.O. Box 9293 Singapore 917699). Manufacturer: Little Doctor Electronic (Nantong) Co. Ltd. (No. 8, Tongxing Road Economic & Technical Development Area, 226010 Nantong, Jiangsu, People's Republic of China).

Distributor in Europe: Little Doctor Europe Sp. z o.o., 57G Zawila Street, 30-390, Krakow, Poland, Sales Office phone: +48 12 2684746, 12 2684747, fax: +48 12 268 47 53. E-mail: biuro@littledoctor.pl www.LittleDoctor.pl Authorized Representative in the EU: Little Doctor Europe Sp. z o.o., 57G Zawila Street Krakow 30-390 Poland. For more information please visit www.littledoctor.sp

* - Measured by Little Doctor International (S) Pte. Ltd.

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Produs sub control și Little Doctor International (S) Pte.Ltd. (7500A BEACH ROAD #11-313 THE PLAZA SINGAPORE
* – Măsurată de Little Doctor International (S) Pte. Ltd.