
DE Inhalator Gebrauchsanweisung	4
EN Nebuliser Instructions for use.....	16
FR Inhalateur Mode d'emploi	27
ES Inhalador Instrucciones de uso.....	40
IT Aerosol Istruzioni per l'uso	52
TR Nebulizatör Kullanım kılavuzu	64
RU Ингалятор Инструкция по применению	75
PL Inhalator Instrukcja obsługi	89
NL Inhalator Gebruiksaanwijzing	101
DA Inhalator Betjeningsvejledning	113
SV Inhalator Bruksanvisning	124
NO Inhalator Bruksanvisning	135
FI Inhalaattori Käyttöohje	146

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Read these instructions for use carefully and keep them for later use. Make them accessible to other users and note the information they contain.

Contents

1. Included in delivery	16
2. Signs and symbols	16
3. Intended use	17
4. Warnings and safety notes	18
5. Device and atomiser set description	19
6. Initial use	20
7. Usage	20
8. Changing the filter	22
9. Cleaning and maintenance	22
10. Replacement parts and wearing parts	24
11. What if there are problems?	25
12. Disposal	25
13. Technical specifications	26
14. Warranty/service	26

1. INCLUDED IN DELIVERY

Check that the exterior of the cardboard delivery packaging is intact and make sure that all contents are present. Before use, ensure that there is no visible damage to the device or atomiser (=year pack) and that all packaging material has been removed. If you have any doubts, do not use the device and contact your retailer or the specified Customer Services address.

See device and atomiser description

- 1 x nebuliser
- 1 x atomiser
- 1 x compressed air hose
- 1 x mouthpiece
- 1 x silicone baby mask
- 1 x silicone children's mask
- 1 x angled fitting
- 10 x replacement filters
- 1 x storage bag
- 1 x set of these instructions for use

2. SIGNS AND SYMBOLS

The following symbols are used on the device, in these instructions for use, on the packaging and on the type plate for the device:

	Warning Warning indicating a risk of injury or damage to health
	Important Safety note indicating possible damage to the device/atomiser
	Product information Note on important information
	Application part, type BF
	Observe the instructions for use
	Protection class 2 device
	Manufacturer
I	On
O	Off
30 min. ON/ 30 min. OFF	30 minutes of operation, then 30 minutes break before operating again.
IP21	Protected against foreign objects ≥ 12.5 mm and against vertically falling drops of water
	CE labelling This product satisfies the requirements of the applicable European and national directives.
	Marking to identify the packaging material. A = Material code, B = Material number: 1-6 = Plastics, 20-22 = Paper and cardboard
	Separate the packaging elements and dispose of them in accordance with local regulations.
	Separate the product and packaging elements and dispose of them in accordance with local regulations.
	Disposal in accordance with the Waste Electrical and Electronic Equipment EC Directive - WEEE
	Permissible storage and transport temperature and humidity
	Medical device

	Unique Device Identifier (UDI) Identifier for unique product identification
	Item number
	Single patient – multiple use
	Serial number

- Systemic diffusion is reduced substantially.
- Only very low doses of the medication are required.
- Quick and effective therapeutic effect.
- Side effects are significantly reduced compared to systemic administration.
- Humidifying the airways.
- Loosening and making (bronchial) secretions more liquid.
- Releasing bronchospasms (spasmolysis).
- Relieving swollen and inflamed bronchial mucosa.
- Clearing secretions by coughing.
- Counteracting pathogens of respiratory infections.

3. INTENDED USE

Purpose

Nebulisers (including compressor, ultrasonic and mesh nebulisers) are medical devices for the nebulisation of liquids and liquid medication (aerosols). This device produces aerosols by combining compressed air and liquid medication. The aerosol treatment is suitable for treating the upper and lower airways. By nebulising and inhaling the medication prescribed/ recommended by your doctor, you can prevent diseases affecting the airways, or in the case that you contract such an illness, you can alleviate symptoms and speed up your recovery.

Target group

The inhaler is intended only for medical care in the domestic environment. The inhaler is not intended for health care for professional institutions.

The nebuliser is intended for use in the home healthcare environment, not in professional healthcare facilities. The nebuliser can be used on anyone over 6 month of age under supervision; it can be used for self-treatment by anyone over 12 years of age.

The application under supervision depends on the respective face form of the person to be treated. In this respect, the application is possible under supervision if necessary before or later only later. At the mask inhalation, make sure that the mask rests well and the eyes remain free.

Clinical benefits

Inhalation is the most effective way to administer medication for most respiratory system conditions.

The benefits are:

- The medication is transported directly to the target organs.
- The local bioavailability of the medication is significantly increased.

Indication

Examples of upper respiratory system diseases include:

- Nasal mucosal inflammation
- Allergic nasal mucosa inflammation
- Nasal sinus infection
- Pharyngitis
- Laryngitis

Examples of lower respiratory system diseases include:

- Bronchial asthma
- Bronchitis
- COPD (chronic obstructive pulmonary disease)
- Bronchiectasis
- Acute tracheobronchitis
- Cystic fibrosis
- Pneumonia

Contraindications

- The nebulizer is not intended for the treatment of life-threatening conditions.
- This device must not be used by children under the age of 12 and by people with reduced physical, sensory (e.g. reduced sensitivity to pain) or mental skills or a lack of experience or knowledge, unless they are supervised or have been instructed on how to use the device safely, and are fully aware of the consequent risks of use.
- Do not use the device on persons who are ventilated and/or unconscious.
- Check whether there are contraindications for use with the usual systems for aerosol treatment on the medication instruction leaflet.

- If the device does not work properly, or if you feel unwell or experience pain, stop using it immediately.

4. WARNINGS AND SAFETY NOTES

Warning

- The atomiser has no significant influence on the efficacy and safety of the medication administered and is not intended for the treatment of life-threatening conditions.
- If the device malfunctions, please refer to the “What if there are problems” section.
- The device is not a substitute for medical consultation and treatment. Always consult your doctor first in the event of pain or illness.
- If you have health concerns of any kind, consult your GP.
- Please observe the general hygiene measures when using the atomiser.
- You should always follow the instructions of your doctor regarding the type of medication to use, the dosage, and the frequency and duration of inhalation.
- Only use medication prescribed or recommended by your doctor or pharmacist.
- For treatment, only use parts indicated by your doctor according to the particular diagnosis.
- Keep the device away from your eyes when it is in use, as the mist of medication could be harmful.
- Never use the device near flammable gases, oxygen or nitrogen oxide.
- The device must be switched off and the plug pulled out before every cleaning procedure.
- Keep packaging material away from children (risk of suffocation).
- Children must not play with the device.
- Under no circumstances should you open or repair the device yourself, as faultless functionality can no longer be guaranteed thereafter. Failure to comply will invalidate the warranty.
- To avoid the risk of entanglement and strangulation, store cables and air lines out of the reach of small children.
- Do not use any additional parts that are not recommended by the manufacturer.
- The device must only be connected to the mains voltage that is specified on the type plate.
- Never submerge the device in water and do not use it in the bathroom. Under no circumstances may liquid enter the device.
- Protect the device from heavy impacts.
- Never touch the mains cable with wet hands, as you could get an electric shock.
- Do not pull the mains plug out of the socket using the mains cable.
- Do not crush or bend the mains cable, pull it over sharp-edged objects or leave it dangling down, and protect it from sources of heat.
- We recommend that the power cable is completely unrolled to avoid dangerous overheating.
- If the mains connection cable of this device is damaged, it must be disposed of. If it cannot be removed, the device must be disposed of.
- If the device is opened, there is a risk of electric shock. Disconnection from the power supply network is only guaranteed if the mains plug is unplugged.
- Modification of the device and atomiser is not permitted.
- If the device has been dropped, exposed to high levels of moisture or suffered any other damage, it must no longer be used. If in doubt, contact Customer Services or the retailer.
- The nebuliser may only be operated with compatible Beurer atomisers and with the appropriate Beurer atomiser set. The use of atomisers and atomiser sets made by other companies may result in less efficient treatment and could damage the device.
- Store the device and atomiser out of the reach of children and pets.

Important

- Power cuts, sudden interferences or other unfavourable conditions could lead to the device becoming inoperable. We therefore recommend that you obtain a replacement device or medication (the latter should be agreed with your doctor).
- Should you require an adapter or extension lead, this must meet the applicable safety requirements. The power limit and the maximum output specified on the adapter must not be exceeded.
- Never store the device (or the power cable) near to sources of heat.
- Do not use the device in a room in which a spray has previously been used. Air the room before carrying out the treatment.

- Do not allow objects to get into the cooling vents.
- Never use the device if it is making an abnormal sound.
- For hygiene reasons, it is essential that every user uses their own atomiser.
- Always disconnect the mains plug from the mains socket after use.
- Store the device in a location protected against climatic influences. The device must be stored in the environmental conditions specified.

FUSE

- The device is fitted with an overcurrent protection device. This may only be changed by authorised specialist personnel.

General notes

Important

- Only use the device
 - on people,
 - for the purpose for which it is designed (aerosol inhalation) and in the manner specified in these instructions for use.
- **Any form of improper use can be dangerous.**
- In the event of an acute emergency, providing first aid must take top priority.
- Apart from the medication, only use distilled water or a saline solution. Other liquids may cause a fault in the nebuliser or atomiser.
- This device is not intended for commercial or clinical use; it is designed exclusively for self-treatment in a private home.

Prior to initial use

Important

- Remove all packaging material before using the device.
- Protect the device against dust, dirt and humidity and never cover the device while it is in use.
- Do not operate the device in a very dusty area.
- Switch the device off immediately if it is faulty or not working properly.
- The manufacturer is not liable for damage resulting from improper or incorrect use.

Repairs

Note

- Under no circumstances should you open or repair the device yourself, as faultless functionality can no longer be guaranteed thereafter. Failure to comply with this instruction will void the warranty.
- The device is maintenance-free.
- For repairs, please contact Customer Services or an authorised retailer.

5. DEVICE AND ATOMISER SET DESCRIPTION

Nebuliser

The corresponding drawings are shown on page 3.

- | | |
|---------------------------------|---------------------------------|
| 1 Carrying handle | 4 Hose connection |
| 2 On/Off switch | 5 Filter cap with filter |
| 3 Mains connection cable | 6 Holder for atomiser |

Atomiser and atomiser set

- | | |
|--------------------------------|------------------------------------|
| 7 Compressed air hose | 12 Mouthpiece |
| 8 Atomiser | 13 Silicone children's mask |
| 9 Atomiser insert | 14 Silicone baby mask |
| 10 Valve | 15 Angled fitting |
| 11 Medication container | 16 Replacement filter |

Nasal douche

- | | |
|--|----------------------------|
| 17 Cover | 19 Saline container |
| 18 Collection container for nasal secretion | 20 Cover hole |

Note

The nebuliser can also be used by adults. Adults have the option of using a nasal douche and a separate atomiser set (includes mouthpiece, nosepiece, adult mask, children's mask, atomiser, compressed air hose, filter) for their nebuliser. Neither are included in the scope of delivery of the nebuliser [n/a]. Nevertheless, the appropriate use of the nasal douche, nosepiece, etc. is described below. You can find an overview of all replacement items in section "Replacement parts and wearing parts".

6. INITIAL USE

Setting up the device

Take the device out of the packaging. Place the device on a flat surface. Make sure that the vents are free from obstructions.

Before using the device for the first time

i Note

- Clean and disinfect the atomiser and replacement set before using them for the first time. See section “Cleaning and maintenance”.
- Connect the compressed air hose [7] to the medication container [11]. **A 1**
- Connect the other end of the compressed air hose [7] to the nebuliser hose connection [4] by turning it slightly. **A 2**

Mains connection

The device must only be connected to the mains voltage that is specified on the type plate.

- Insert the mains plug of the mains connector cable [3] into a suitable socket.
- To connect to the mains, fully insert the mains plug into the socket.

i Note

- Ensure that there is a mains socket close to the set-up area.
- Arrange the mains cable so that no-one will trip over it.
- To disconnect the nebuliser from the mains after inhalation, first switch off the device and then pull the mains plug out of the socket.

7. USAGE

! Important

- For hygiene reasons, it is essential to clean the atomiser [8] and the atomiser set parts after each treatment and to disinfect them after the last treatment of the day.
- The atomiser set parts may only be used by one person; use by several people is not recommended.
- If the therapy involves inhaling several different medications one after the other, please be aware that the atomiser [8] must be rinsed

under warm tap water following every usage. See section “Cleaning and maintenance”.

- Please observe the notes on changing the filter in these instructions for use.
- Check that the hose connectors are firmly attached to the nebuliser [4] and the atomiser [8] before each use of the device.
- Before use, check that the device is working correctly by briefly switching on the nebuliser (together with the connected atomiser, but without medication). If air comes out of the atomiser [8], the device is working.

7.1 Insert the atomiser insert

- Open the atomiser [8] by twisting the top section anticlockwise against the medication container [11]. Place the atomiser insert [9] in the medication container [11].
- Ensure that the cone for administering medication fits well on the cone for the air duct inside the atomiser [8].
- Insert the medication container [11] into the holder [6].

7.2 Fill the atomiser

- Fill with an isotonic saline solution or pour the medication directly into the medication container [11]. Avoid overfilling. The maximum recommended filling quantity is 8 ml.
- Use medication only as instructed by your doctor and ask about the appropriate inhalation period and quantity for you.
- If the prescribed quantity of medication is less than 2 ml, top this up to at least 4 ml with isotonic saline solution. Dilution is also necessary with viscous medications. Here too, please observe the instructions of your doctor.

7.3 Close the atomiser

- Close the atomiser [8] by twisting the top section clockwise against the medication container [11]. Ensure that it is connected correctly.

7.4 Insert the valve

- Press the valve [10] into the opening intended for this purpose on the upper side of the atomiser [8].

i Note

The valve system ensures better delivery of the medication to the lungs. Use the valve [10] only

with the mouthpiece **12** or masks **13** **14** and never with the nosepiece (n/a).

7.5 Connect the atomiser set to the atomiser

- Connect the atomiser **8** to the desired part of the atomiser set. The mouthpiece **12** or nosepiece [n/a] can be directly attached. The children's mask **13** and baby mask **14** must be attached using the angle piece **15** to the atomiser.

i Note

The most effective form of nebulisation is by using the mouthpiece. Nebulisation using a mask is only recommended if it is not possible to use a mouthpiece (e.g. for children who are not yet able to inhale using a mouthpiece).

When using a mask to inhale, take care to ensure the mask fits well and the eyes are unobstructed.

- Before the treatment, pull the atomiser **8** upwards out of the holder **6**.
- Start the nebuliser using the On/Off switch **2**.
- If there is mist pouring out, this indicates that the device is operating correctly.

7.6 Treatment

- Hold or place the baby or infant to be treated in the most upright position possible. Do not use the device when the baby or infant to be treated is lying down.
- When inhaling, older children should sit upright and relaxed at a table and not in an armchair to avoid compressing the airways and therefore impairing the effectiveness of the treatment.
- Breathe in the atomised medication deeply.

! Important

The device is not suitable for continuous operation; after 30 minutes of operation it must be switched off for 30 minutes.

i Note

During the treatment, hold the atomiser straight (vertically), otherwise the atomisation will not work and faultless functionality is not guaranteed.

! Important

Essential oils, cough syrups, gargling solutions and drops to be used as a rub or in a steam bath are wholly unsuitable for inhalation using

a nebuliser. These additives are often viscous and can impair the correct functioning of the device and therefore the effectiveness of the application in the long term. For individuals with a hypersensitive bronchial system, medications containing essential oils may under certain conditions cause an acute bronchospasm (a sudden cramp-like restriction of the bronchi with shortness of breath). Consult your doctor or pharmacist in relation to this matter.

7.7 Stop the inhalation

Once the mist is only coming out in an irregular flow or if the sound changes when inhaling, you can stop the treatment.

- Switch off the nebuliser after treatment using the On/Off switch **2** and disconnect it from the mains.
- Place the atomiser **8** back in its holder **6** after the treatment.

7.8 Nasal douche

The nasal douche is not included in the scope of delivery [n/a]. You can find an overview of all replacement items in section "Replacement parts and wearing parts".

You can use the nasal douche to clean your nasal cavities. The nasal douche creates a fine jet of liquid for treating respiratory tract problems or for humidifying the nasal mucous membranes. We recommend using the nasal douche for colds, dryness inside the nose, dust exposure, pollen allergies, chronic inflammation of the sinuses or upon the advice of your doctor.

! Warning

- Do not use the nasal douche in combination with medication or essential oils.

! Important

- Only use the nasal douche with isotonic saline solutions.
- Replace the nasal douche after a year.

The maximum volume of the nasal douche is 10 ml. If you are unsure about the correct dilution of the saline solution, please consult your pharmacist.

7.9 How to use the nasal douche

To use the nasal douche correctly, proceed as follows **B**:

1. Twist the cover **17** 90° anticlockwise and remove it in an upward direction.

- Lift the collection container **18** out.
- Fill the saline container **19** with saline solution.
- Place the collection container **18** in the saline container **19**.
- Place the cover **17** on the saline container **19**. Twist the cover **17** in a clockwise direction until it clicks firmly into place.
- Connect the compressed air hose **7** to the saline container **19** from below.
- Position the cover **17** carefully inside one of your nostrils. Breathe through your nose.
- Press the On/Off switch **2** on the back of the nebuliser. To start the application, place one finger over the cover hole **20** on the container.

Note

During the application, breathe slowly and evenly through your nose. Tilt your head slightly in the opposite direction to the nostril being treated so that the saline solution can flow deep into the nose.

- To pause or stop the application, take your finger off the cover hole **20** on the saline container **19**.
- Once your nose feels clear again, you can stop the application. The application should be stopped at the latest once the maximum filling quantity has been used.

7.10 Clean the components

See “Cleaning and maintenance”.

8. CHANGING THE FILTER

In normal operating conditions, the air filter must be replaced after approx. 500 operating hours or one year. Please check the air filter regularly (after 10-12 nebulisation procedures).

Replace the used filter if it is very dirty or clogged. If the filter has become damp, it must also be exchanged for a new filter.

Important

- Do not attempt to clean the used filter and reuse it.
- Only use the manufacturer’s original filter, otherwise your nebuliser may become damaged and sufficiently effective treatment cannot be guaranteed.
- Do not repair or maintain the air filter while it is being used by someone.

- Never operate the device without a filter.

To replace the filter, proceed as follows:

Important

- First switch the device off and disconnect it from the mains.
 - Allow the device to cool down.
- Pull off the filter cap **5** towards the front. **C**

Note

If the filter remains in the device after the cap has been removed, take the filter out of the device, e.g. with tweezers or similar.

- Re-insert the filter cap **5** with a new filter.
- Ensure that it is securely in place.

9. CLEANING AND MAINTENANCE

Atomiser and atomiser set

Warning

Adhere to the following hygiene instructions to avoid health risks.

- The atomiser **8** and atomiser set are designed for multiple use. Please note that different areas of application involve different requirements in terms of cleaning and hygienic preparation.

Notes

- Do not clean the atomiser or the atomiser set mechanically using a brush or similar device, as this could cause irreparable damage and it will mean that the best treatment results can no longer be guaranteed.
- Please consult your doctor about the additional requirements in terms of the hygienic preparation required (hand care, handling of medication/inhalation solutions) for high-risk groups (e.g. patients with cystic fibrosis).
- Ensure thorough drying after each cleaning or disinfection process. Residual moisture or wetness can represent an increased risk of bacterial growth.

Preparation

- Immediately after each treatment, all parts of the atomiser **8** and the atomiser set used

must be cleaned of residual medication and contamination.

- To do this, dismantle the atomiser [8], the masks [13] [14] and the nasal douche [n/a] into their individual parts.
- Remove the mouthpiece [12], the mask [13] [14] or the nosepiece [n/a] from the atomiser.
- Disassemble the nosepiece [n/a] if you have used it with the comfort attachment.
- Dismantle the atomiser by twisting the top section anticlockwise against the medication container [11].
- Remove the atomiser insert [9] from the medication container [11].
- Remove the valve [10] by pulling it out of the atomiser [8].
- Reassembly is carried out in reverse order.

Cleaning

Important

The device must be switched off, disconnected from the mains and allowed to cool down each time before cleaning.

The atomiser and the used atomiser set, such as the mouthpiece [12], mask [13] [14], atomiser [8], nasal douche etc. must be washed with hot but not boiling water after each use. Dry the parts carefully using a soft cloth. Put the parts together again when they are completely dry and place them in a dry, sealed container or disinfect them.

When cleaning, ensure that any residue is removed. Never use any substances for cleaning that could potentially be toxic if they came into contact with the skin or mucous membranes, or if they were swallowed or inhaled.

Use a soft, dry cloth and non-abrasive cleaning products to clean the device.

Do not use any abrasive cleaning products and never submerge the device in water.

Important

- Ensure that no water gets inside the device.
- Do not clean the device or replacement set in the dishwasher.
- Do not touch the device with wet hands when it is plugged in and do not allow water to spray on the device.
- Only operate the device if it is completely dry.
- Do not spray any liquid into the vents! If liquid penetrates the device this could cause

damage to the electrics or other nebuliser parts and lead to a malfunction.

Condensation, hose care

Condensation may form in the hose depending on the ambient conditions. It is essential to remove the moisture to prevent bacterial growth and ensure proper treatment. To do so, proceed as follows:

- Disconnect the compressed air hose [7] from the atomiser [8].
- The hose must remain in the nebuliser [4].
- Operate the nebuliser until the moisture is removed by the air passing through.
- In the event of heavy contamination, replace the hose.

Disinfection

Please carefully observe the points below when disinfecting your atomiser and atomiser set. We recommend disinfecting the individual parts on a daily basis after the last usage as a minimum measure.

(A small amount of clear vinegar and distilled water can be used.)

- First, clean the atomiser and atomiser set as described in the “Cleaning” section.
- Place the disassembled atomiser [8] and nasal douche [n/a], the mouthpiece [12], the silicone masks [13] [14] and the disassembled nosepiece [n/a] in boiling water for 5 minutes.
- For the remaining atomiser set parts, use a vinegar solution consisting of ¼ vinegar and ¾ distilled water. Make sure that the volume is sufficient to fully submerge parts, such as the PVC masks [n/a] and hose [7], in the solution. Atomiser set parts that can be boiled can also be disinfected in this way.
- Leave the parts in the vinegar solution for 30 minutes.
- Rinse the parts with water and dry them carefully with a soft cloth.

Important

Do not boil or autoclave the compressed air hose and the PVC masks [n/a].

- Put the parts together again when they are completely dry and place them in a dry, sealed container.

Note

Please ensure that the parts are completely dried after cleaning, otherwise the risk of bacterial growth is increased.

Drying

- Place the individual parts on a clean, dry and absorbent surface and leave them to dry completely (at least 4 hours).

Durability of materials

- As with any plastic parts, atomisers and replacement sets are affected by a certain amount of wear and tear when used and hygienically prepared on a frequent basis. Over time, this can lead to a change in the aerosol, which can have a negative effect on the efficiency of the treatment. We therefore recommend that you replace the atomiser and any replacement set after a year.
- The following point should be observed when selecting the cleaning product or disinfectant: only use a mild cleaning product or disinfectant in the concentration and volume prescribed by the manufacturer.

Storage

- Do not store in damp conditions (such as in a bathroom) and do not transport with any damp items.
- When storing and transporting, protect from prolonged direct sunlight.
- The atomiser set parts can be safely stored in the storage bag. Store the device in a dry place, ideally in the original packaging.

10. REPLACEMENT PARTS AND WEARING PARTS

Atomiser set parts and replacement parts are available from the corresponding service address (according to the service address list). Please state the corresponding order number.

Designation	Material	REF
IH 24 Kids/26 Kids atomiser set = year pack contains: Mouthpiece Silicone children's mask Silicone baby mask Atomiser Compressed air hose Filter	PP/silicone Silicone/PP Silicone/PP PP/silicone PVC PU	601.19

Designation	Material	REF
Standard IH 26 atomiser set = year pack contains: Mouthpiece Nosepiece Adult mask Children's mask Atomiser Compressed air hose Filter	PP/silicone PP/silicone PVC/aluminium PVC/aluminium PP/silicone PVC PU	601.28
Nasal douche	PP	601.37
Baby mask	PVC	601.31

Note

If the device is not used according to the instructions specified, perfect functionality cannot be guaranteed! We reserve the right to make technical changes to improve and develop the product. This device and its atomiser set comply with the European standards EN 60601-1 and EN 60601-1-2 (CISPR 11, IEC 61000-3-2, IEC 61000-3-3, IEC 61000-4-2, IEC 61000-4-3, IEC 61000-4-4, IEC 61000-4-5, IEC 61000-4-6, IEC 61000-4-7, IEC 61000-4-8, IEC 61000-4-11) as well as EN ISO 27427 and are subject to particular precautions with regard to electromagnetic compatibility. Please note that portable and mobile HF communication systems may interfere with this device. For more details, please contact our Customer Services at the address indicated. The device complies with the Regulation (EU) 2017/745 of the European Parliament and of the Council on medical devices, as well as the respective national provisions.

NOTES ON ELECTROMAGNETIC COMPATIBILITY

- The device is suitable for use in all environments listed in these instructions for use, including domestic environments.
- The use of the device may be limited in the presence of electromagnetic disturbances. This could result in issues such as the failure of the device.
- Avoid using this device directly next to other devices or stacked on top of other devices, as this could lead to faulty operation. If, however, it is necessary to use the device in the manner stated, this device as well as the other devices must be monitored to ensure they are working properly.

- The use of a atomiser set other than that specified or provided by the manufacturer of this device can lead to an increase in electromagnetic emissions or a decrease in the device's electromagnetic immunity; this can result in faulty operation.
- Keep portable RF communication devices (including peripheral equipment, such as antenna cables or external antennas) at least 30 cm away from all device parts, including all cables included in delivery.
- Failure to comply with the above can impair the performance of the device.

11. WHAT IF THERE ARE PROBLEMS?

Problem/question	Possible cause/remedy
The atomiser produces no or too little aerosol.	1. Too much or too little medication in the atomiser. Minimum: 2 ml, Maximum: 8 ml.
	2. Check nozzle for blockages. Clean nozzle if necessary (e.g. by rinsing out). Then start using the atomiser again. IMPORTANT: carefully pierce the fine holes only from the underside of the nozzle.
	3. Atomiser not held vertically.
The atomiser produces no or too little aerosol.	4. Unsuitable medication fluid added for nebulisation (e.g. too viscous). The medication fluid should be prescribed by a doctor.
The output is too low.	Kinked hose, clogged filter, too much inhalation solution.
What medications are suitable for inhaling?	Please consult your doctor in relation to this matter. As a rule, all medication that is suitable and approved for device inhalation can be inhaled.
There is inhalation solution residue in the atomiser.	This is normal and is due to technical reasons. Stop inhalation once the atomiser starts to make a notably different sound.

Problem/question	Possible cause/remedy
What should be taken into account when using the device with babies and children?	<ol style="list-style-type: none"> 1. On babies, the mask should cover the mouth and nose to ensure effective inhalation. 2. On children, the mask should also cover both the nose and mouth. It is not very useful to carry out nebulisation next to someone who is sleeping, as in this case not enough of the medication will reach the lungs. <p>Note: inhalation should only be carried out under the supervision of an adult and with their assistance, and the child should not be left alone.</p>
Why should the atomiser be replaced regularly?	<p>There are two reasons for this:</p> <ol style="list-style-type: none"> 1. To ensure a therapeutically effective particle spectrum, the nozzle hole must not exceed a specific diameter. Due to mechanical and thermal stresses, the plastic is subject to a certain degree of wear. The atomiser insert  is particularly delicate. This can also alter the composition of the aerosol droplets, which directly affects the effectiveness of the treatment. 2. Regularly changing the atomiser is also recommended for hygiene reasons.
Should each person have their own atomiser?	Yes, this is absolutely essential for hygiene reasons.

12. DISPOSAL

For environmental reasons, do not dispose of the device in the household waste.

Please dispose of the device in accordance with EC Directive – **WEEE** (Waste Electrical and Electronic Equipment).

If you have any questions, please contact the local authorities responsible for waste disposal.


13. TECHNICAL SPECIFICATIONS

Model	IH 24 Kids
Dimensions (W x H x D)	150 x 210 x 150 mm
Weight	1.28 kg
Operating pressure	Approx. 0.8 to 1.55 bar
Atomiser filling volume	Max. 8 ml Min. 2 ml
Medication flow rate	Approx. 0.35 ml/min
Sound pressure	Approx. 52 dBA (according to DIN EN 27427 Section 26)
Mains connection	230 V~; 50 Hz; 230 VA UK: 240 V~; 50 Hz; 240 VA Saudi Arabia: 220 V~; 60 Hz; 220 VA
Operating conditions	Temperature: +10 °C to +40 °C Relative humidity: 10% to 95% Ambient pressure: 700 hPa to 1060 hPa
Storage and transportation conditions	Temperature: -25 °C to +70 °C Relative humidity: 10% to 95% Ambient pressure: 700 hPa to 1060 hPa
Aerosol performance according to EN27427:2019 based on adult ventilatory pattern with sodium fluoride (NaF):	1) Aerosol delivery: 0.11 ml 2) Aerosol delivery rate: 0.07 ml/min 3) Fill volume dispensed in percent per min.: 3.5% 4) Residual volume: 1.68 ml 5) Particle size (MMAD): 3.07 µm 6) GSD (geometrical standard deviation): 2.03 7) RF (Respirable fraction < 5 µm): 75% 8) Large particle range (>5 µm): 25% 9) Medium particle range (2 to 5 µm): 47.8% 10) Small particle range (<2 µm): 27.2%

The serial number is located on the device or in the battery compartment.
Subject to technical changes.

Particle size diagram

Measurements were performed using a sodium fluoride solution with a “Next Generation Impactor” (NGI). The diagram may therefore not be applicable to suspensions or highly viscous medications. You can obtain more detailed information from the manufacturer of your medication.

14. WARRANTY/SERVICE

Further information on the warranty and warranty conditions can be found in the warranty leaflet supplied.

Notification of incidents

For users/patients in the European Union and identical regulation systems (EU Medical Device Regulation (MDR) 2017/745), the following applies: If during or through use of the product a major incident occurs, notify the manufacturer and/or their representative of this as well as the respective national authority of the member state in which the user/patient is located.

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