## **Newel Curved Rails**

**VY400/VY416**

## **Fixing and Maintenance Instructions**

## INTRODUCTION

Thank you for deciding to purchase your Newel Curved Rail. This range is manufactured from the finest materials and components available. When installed and used correctly they are designed to give many years of reliable trouble-free service.

**NB. This equipment must be installed by a competent person.**

## BEFORE USE

- Carefully remove all packaging. Avoid using any knives or other sharp instruments as this may damage the surface of the product.
- Check the product for any visible damage. If you see any damage or suspect a fault, do not use your product as it may compromise your safety, but contact your supplier for support.

## FIXING GRAB BARS/GRAB RAILS

Firstly check the substrate where you intend to fix your Newel Curved Rail, ensuring that it is sound. Also check the type of substrate you intend to fix to, e.g. brick/block, timber or concrete.

**Never drill through pre-stressed concrete to install a Newel Rail product.**

Always check the type of fixing you are going to use on this substrate for their compatibility.

**We suggest that for sound brick and block that a minimum screw length of 50 mm x size 10 is used with a suitable rawl plug and screws, in accordance with the guidelines issued by the fixings manufacturer.**

## PLASTER BOARD (STUD TYPE) LATHE AND PLASTER

When attempting to fix to plasterboard stud or lathe and plaster, we suggest that you check the position of the studs and span them with a suitable timber plate of a minimum thickness of 40 mm x 150 mm high, securing the plate to the studs. After ensuring the soundness of the studs (using a suitable counter sunk screw), fix your grab bar to the timber plate.

**Never install any grab bar or grab rail directly onto plaster board, MDF (medium density fibreboard), chipboard or other composite material. This will seriously compromise the safety of the user.**

**Always use all screw holes provided in the Newel Rail.**

## ELECTRICITY

Always check for concealed electrical cables when installing hinged arm supports. Also check your legal requirements under the NICE IC 16th Edition relating to Earth bonding of extraneous metal parts.

## **CLEANING**

The Newel Curved Rail should be cleaned using a non-abrasive cleaner or mild detergent with a soft cloth. Only use a mild detergent or disinfectant; never use an abrasive cleaner. Do not exceed 50°C, and never autoclave arm supports. Always ensure that you wipe dry all equipment and treat joints with a water dispersing agent, e.g. WD40.

## **REISSUE**

If you reissue or are about to reissue this product, please thoroughly check all components for their safety.

**If in any doubt, please do not issue or use, but immediately contact your supplier for service support.**

## **CARE & MAINTENANCE**

Please make a safety check of the product at regular intervals or if you have any concerns.

## **IMPORTANT INFORMATION**

The information given in this instruction booklet must not be taken as forming part of or establishing any contractual or other commitment by Aidapt Bathrooms Limited, Aidapt (Wales) Ltd or its agents or its subsidiaries and no warranty or representation concerning the information is given.

Please exercise common sense and do not take any unnecessary risks when using this product; as the user you must accept liability for safety when using the product.

Please do not hesitate to contact the person who issued this product to you or the manufacturer (detailed below) if you have any queries concerning the assembly/use of your product.

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**Altai France (Altai Europe Ltd)**