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IMPORTANT! READ THIS MANUAL BEFORE USE

In this smell training you will use so-called essential oils. You must handle these carefully. You should keep essential oils OUT OF THE REACH OF CHILDREN and DO NOT SWALLOW them. Although the bottles in this smell training only contain a low volume of essential oils, the greatest danger is that if they are swallowed, the oil can be expectorated into the oesophagus. The essential oil can end up in the lungs, causing a chemical inflammation in the lungs. Never leave a recently used bottle at a location that children have easy access to. DO NOT INDUCE VOMITING AFTER INGESTION due to the risk of damage to the oesophagus and/or lungs, but immediately call a doctor and take the packaging with you.

Small children like to taste from bottles. Although the bottles used for this smell training are equipped with a child-safe cap, it is still essential to store essential oils in places where children cannot reach them. Please return the vial and residual oil to the environmental service/waste disposal station of your municipality.

Smell training guide

Your attending physician expects a positive effect from this smell training. However, your doctor cannot tell you what final result you will personally achieve.

Scientific research indicates that the best result is achieved with 24 weeks of training. During the first twelve weeks you will train with the scents from SET 1 (blue). After 12 weeks you will continue the smell training with the scents from SET 2 (green). After these 12 weeks you continue the training with the scents from SET 3 (orange). In total you train 36 weeks.

Performing the smell training

For the training you use the scents in this kit. This allows you to perform the smell training as follows:

1. Take one of the bottles.
2. Hold it about three centimetres from the underside of your nose. It is important that the oil does not touch your nose.
3. Relax and calmly take a deep breath through your nose. If you inhale too quickly, the smell will not be breathed in properly.
4. Repeat this four times; you will then have inhaled the fragrance five times.
5. Rest for a few minutes before proceeding with the next fragrance.
6. Repeat the same process for all scents. You can choose the order of scents yourself.
7. Preferably repeat the above steps twice a day, for example in the morning and in the evening. Close the bottles after use to prevent the scent from evaporating.

Storage and shelf life of the set after opening

Store at room temperature. This set can be kept for three months after opening.

PLEASE NOTE! An (odourless) diluent, namely isopropyl palmitate, has been added to the oil in the Rose and Jasmine vial. This is because Rose and Jasmine would make the smell training sets unnecessarily expensive in their undiluted state. However, this diluent will solidify/crystallize if it gets too cold (e.g. if you keep the bottle in the refrigerator or if the smell training set is delivered during cold weather). As soon as the oil returns to room temperature, it will liquefy again, but this has no effect on the smell quality.