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
Neck stretcher and back brace

Our mission

We strive for a world where people get the most out of life. We do this by creating sustainable products where experience and convenience are key, along with excellent service. Thanks to our products you are one step closer to living your ideal life.

Why the combination of these two products?

Your upper back and neck are closely connected. Often, neck complaints find its origin in the upper back, the logical result of an often prolonged bad posture. The back brace supports the position of your shoulders, which will be pulled back, simulating a more natural position for your upper back. The neck stretcher ensures that your neck muscles slowly return to natural length, which makes the neck tension disappear. Be informed: the current tensions you may be experiencing right now may be the result of sometimes years of bad posture. Hence, it may take several days to weeks before you will see optimal results with the use of these products.

How to use the neck stretcher

Lay the neck stretcher down on the floor, preferably on a fitness mat because this is more comfortable. Make sure that the narrow side of the stretcher faces towards your upper back. Lie down - with your neck on the stretcher - carefully positioning it until it feels comfortable for you. The first few sessions will most likely feel a bit awkward, because your neck is not used to this position. It will soon feel familiar and relaxing though!

Guidelines:

- Use for about 8 to 12 minutes a time
- 1 to 2 times a day

Using the back brace

Putting on the back brace is a very similar process to how you would put on a backpack. The difference is that with this back brace, it is important that you tighten the brace enough. As a result, you will immediately feel the supportive effect as your shoulders are pulled back a bit.

Guidelines:

- Start with once a day 30-60 minutes during everyday activities, like office work / sitting behind a desk
- Steadily build up to 2-3 of these sessions per day

WARNING: Read carefully before use

- Do not exceed recommended usage guidelines.
- Consult a doctor before use if you have any injuries that may be affected by the use of these products.

Disclaimer: No medical claims are warranted or implied by the use of this product. All rights reserved.

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