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ENGLISH

Dear Customer, thank you for choosing a MEDLUXY product. This MEDLUXY device complies with the provisions of European Regulation on Medical Devices, whereby it is classified as a medical product of class I.

Please carefully consider this instruction and read all the provisions prior to use. Follow the given care instructions.

Intended use: Clamping, grasping, holding small objects, tissue, blood vessels, wires and so on. Field: general medicine, eg medical practices.

How to hold: Artery clamps usually have two looped, scissor-like handles, a slender but blunt nose, and a lock to ensure the instrument holds its tight grasp on the blood vessel. The lock allows surgeons and medics to leave multiple clamps on a critically bleeding patient. The forceps' metal teeth, located near the finger loops, slide into place as the instrument's handles are squeezed together. A doctor has only to shift his or her fingers slightly to unlock the teeth and allow the instrument to release its grip. The forceps' tips can be either straight or curved.

Target patient group: All

Improper use/contraindications to the use of the product: Any other use or use outside the above intended use shall be considered improper. The manufacturer is not responsible for damage resulting from such use. In this case, all risks are assumed solely by the user.

Warnings: Keep out of reach of children. Use under supervision of an adult. **Secondary treatment limitations:** Frequent secondary treatment has little effect on these instrument. The service life of the product is usually defined based on the degree of wear and use-related damage. Ultrasonic cleaning is not permitted! High-speed autoclaving in a hot air sterilization will damage the instrument!

Additional information: If used and stored properly, these MEDLUXY products will serve you reliably for many years.

Instructions:
Preparation for cleaning: After using the scissors, we recommend immediately rinsing them under running water or in a slightly alkaline solution to avoid possible drying of residues after use (e.g., blood). In order to properly carry out steam sterilization (see "Cleaning, Sterilization" section), it is always necessary to carry out a thorough pre-cleaning!
Manual cleaning: According to the recommendation of the National Institute for Public Health and Environment (RIVM), the mechanical treatment is preferable. Manual cleaning is not recommended.

Automatic cleaning: Equipment: CM 310 (Maquet), cleaning agent (Neodisher®FA forte 0.4%/Neodisher®Z 0.2%)
 G7828 (Miele), cleaning agent (Mucapur®XL 0.4%/Mucapur®Z 0.15%)
 WD 390 (Belimed), cleaning agent (Mucapur®AF 0.5%/Mucapur®Z 0.1%)

1. Thoroughly wash the instruments under running water before mechanical treatment to avoid the remnants of cleaning/disinfecting agents entering the machine. 2. Place the instruments in a special container. 3. Place the instrument container in the cleaning and disinfection device so that the water jets directly fall on the instruments. 4. Follow the instructions for the detergent dosage. 5. Start the Vario TD program followed by thermal disinfection. Thermal disinfection is provided in accordance with the established value of A0 (the

disinfection time, in seconds, at a certain temperature) and in accordance with national regulations. 6. When the program is completed, remove the instruments from the cleaning and disinfection device and dry them (National Institute for Public Health and Environment (RIVM) recommends to dry with compressed air). Special attention should be paid to the places of the instrument container that are difficult to reach during drying. 7. Visual inspection for integrity and cleanliness with a suitable magnification object (based on experience, 8x magnification is sufficient for visual inspection). If there are residues of contamination on the instruments after mechanical treatment, repeat cleaning and disinfection until they are completely cleaned. 8. Attention! Only mechanical cleaning (without confirmed disinfection) requires subsequent thermal disinfection in a steam sterilizer placed in suitable trays or mesh containers. **Disinfection:** For example, Helipur can be used to disinfect the instruments. The appropriate concentrations can be found in the manufacturer's instructions.

Sterilization: The instruments can be successfully sterilized with the following parameters: Fractional pre-vacuum (three times), sterilization temperature: 134°C, exposure time: at least 3 min. (full cycle), drying time: 10 min.

Maintenance: Damaged parts shall be rejected or replaced. Surgical lubricant can be applied to moving parts, such as hinges.

Functional control and check: Check the ease of movement of the moving parts of the scissors. Scissors assembly: Perform a visual inspection for damage and wear.

Transportation and storage: Keep the instrument away from dust, humidity and contamination during storage and transportation. Temperature: -10°C to +50°C
 Relative humidity: from 30% to 75%

Disposal: Scissors (stainless steel and plastic) must be handed over to a recycling center after prior preparation. Please observe the local regulations.

Warranty: The warranty total 12-months from the date of purchase if used properly and in accordance with the User's Manual. If any issues arise or repair is necessary, please contact your dealer.

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