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
Stand der Information: Februar 2012

For eFlow[®]rapid Type 178G1005

Read these instructions for use carefully and completely before using the device. Keep these instructions in a safe place.

Risks caused by failure to comply with the instructions:

If these instructions are not followed, injuries or damage to the product cannot be ruled out.

Table of contents

1	IMPORTANT INFORMATION	37
	Safety instructions	37
	Intended use.....	37
	Patient groups.....	37
	Medications.....	38
	Liability	39
	Operating conditions.....	39
	Ambient conditions	39
	Operating position.....	39
	Terms and conditions of warranty	39
	Life cycle	40
	Materials used.....	40
	Service and manufacturer.....	40
2	COMPONENTS	41
3	PREPARATION FOR INHALATION	42
	Connection to the power supply	42
	Mains power operation (100 V~ to 240 V~).....	42
	Operation with batteries.....	43
	Operation with rechargeable batteries.....	43
	Assembling the nebuliser handset	44
	Connecting the nebuliser handset.....	46
4	INHALATION.....	47
	Filling and sealing the medication reservoir.....	47
	Performing inhalation.....	48
	Inhalation with mask	48
	Pause function	49
	Display indicators during inhalation.....	49
	End of the inhalation session.....	50

5	HYGIENIC REPROCESSING	51
	General.....	51
	Dismantling the nebuliser handset	52
	Hygienic reprocessing at home	53
	Cleaning	53
	Disinfection.....	54
	Drying, storage, transport.....	55
	Hygienic reprocessing in the hospital and doctor's practice	55
	Material resistance	55
	Cleaning and disinfection	56
	Sterilisation.....	56
	Storage.....	56
6	STORAGE AND TRANSPORTATION CONDITIONS	57
7	TROUBLESHOOTING	58
8	DISPOSAL	62
9	SPARE PARTS AND ACCESSORIES	62
10	TECHNICAL DATA	63
	General.....	63
	Classification according to DIN EN 60601-1.....	63
	Information concerning electromagnetic compatibility	64
11	EXPLANATION OF SYMBOLS	65

1 IMPORTANT INFORMATION

Your eFlow[®]*rapid* inhalation system is an innovative device for the treatment of respiratory diseases. Throughout its development, particular care was taken to consider patients' needs for safe, fast, and above all efficient treatment.

You should have a medical examination before treating illnesses.

Safety instructions

- The device must not be left unattended while it is in operation, or used in damp areas or where there is a risk of explosion.
- The device must not be used or operated if the power adapter housing or the power cable show visible signs of damage, or if a fault is suspected following a fall or similar impact. This may cause injury to individuals. The device must be inspected for damage before every inhalation session.
- Please remove the power plug from the socket to ensure that the power is completely cut off.
- Never remove the power plug from the socket with wet hands. There may be a danger of electric shock.
- Do not allow pets (e.g., rodents) near the cables.

Intended use

The eFlow[®]*rapid* is a reusable electronic inhalation system for the treatment of diseases of the respiratory system and lungs. It is suitable for temporary oral inhalation of liquid medication prescribed or recommended by a doctor for use in the home, hospital or doctor's practice.

The eBase Controller can also be used to operate the drug-specific Altera[®] nebuliser handset. The Altera[®] nebuliser handset is solely designed for inhalation of the medication Cayston[®] (aztreonam lysine). The instructions for cleaning the eFlow[®]*rapid* nebuliser handset apply correspondingly to the Altera[®] nebuliser handset.

Patient groups

The eFlow[®]*rapid* assures highly effective, rapid inhalation treatment for patients of all ages, from infants to adults.

Caution:

In home applications, the nebuliser handset is intended for use by one person only and must not be used by anyone other than that patient! Use by several patients is only permitted in clinics and doctors' practices, and only in conjunction with appropriate sterilisation procedures.

Children:

- Supervised treatment: 2 to 8 years
- Treatment after instruction:
8 years and older

Adults:

To be used as directed (or after competent instruction).

Caution:

- Babies, children and anyone who requires assistance must be supervised constantly by an adult during inhalation therapy. This is the only way to ensure safe and effective treatment. These individuals often underestimate the hazards associated with the equipment (e.g. risk of strangulation with the cable), and this can result in risk of injury.
- The product contains small parts. Small parts can block the respiratory tract and pose a choking hazard. Therefore, be sure to keep them out of the reach of children at all times.
- The device must not be used while oxygen is being administered in a closed environment such as an oxygen tent.
- The nebuliser handsets are only suitable for use by patients who can breathe unassisted and are conscious. This is the only way to ensure effective treatment and avoid a choking hazard.

Medications

The eFlow[®]rapid nebuliser handset must only be used to nebulise permitted inhalation solutions and suspensions from the following groups of agents:

- Antibiotics
- Cromone (DNCG)
- Anticholinergics and b2-sympathomimetics (medications for dilating the bronchia)
- Corticosteroids (anti-inflammatory medications)
- Mucolytics (break down or hydrolyse mucus)
- Saline solutions for inhalation

Note:

i In all cases, follow the instructions for use of the respective medication.

Depending on the composition of the medication, mixing certain medicines can cause adverse chemical or physical reactions to the substances.

When a nebuliser handset (including aerosol head) is used for several medications one after the other, the residual volume of the previous medication must be emptied completely and all nebuliser handset parts must be washed thoroughly with warm tap water before refilling.

Caution:

Using other liquids, such as essential oils, can cause serious health problems. Medications that are approved for inhalation with a medication-specific nebuliser handset must not be used with an eFlow[®]rapid nebuliser handset (risk of incorrect dosage).

Liability

The eFlow[®]*rapid* is an electrically operated device.

Whenever the device is used, a thorough knowledge of the instructions for use and compliance with the same is assumed.

The device must only be operated in accordance with its intended use. The device must only be serviced by PARI Pharma GmbH or a service location expressly authorised by PARI Pharma GmbH. Only original PARI accessories must be used to operate the device.

PARI Pharma GmbH is not responsible for damage or malfunctions caused by the operator using the machine improperly or contrary to the instructions.

Operating conditions

Ambient conditions

- Ambient temperature: +10 °C to +40 °C
- Relative humidity of ambient air: 15% to 93% (non-condensing)
- Air pressure: 700 hPa to 1060 hPa

Medications intended for nebulising must be at room temperature.

Operating position

The nebuliser handset must be held horizontally during inhalation.

Caution:

If the nebuliser handset is tilted, the dosage may not be accurate.

Terms and conditions of warranty

During the warranty period, we or our customer service will repair any defects resulting from faults in material or workmanship free of charge. The owner shall not be entitled to cancel the sale or demand a partial or full refund of the purchase price. This warranty does not cover damage caused by improper use of the device.

The warranty shall become void if repairs are carried out on the device by unauthorised persons.

No compensation will be paid for consequential or direct damage.

In the event of a claim, bring the complete device to your dealer or call our Service Center. Wearing parts such as the nebuliser and the aerosol head are not covered by the warranty.

The warranty period is 2 years. The warranty period shall commence on the date of purchase.

Life cycle

The nebuliser handset is designed for multiple uses. The individual components of the eFlow[®]rapid nebuliser handset are subject to changing stresses during treatment and during preparation for hygienic reprocessing. The frequency and duration of use are decisive factors. The distinction must also be made between use at home (without a change of patients) and use in a

clinic or a doctor's practice (change of patients is possible).

When the life cycle has expired, it is recommended that the components be replaced to ensure proper functioning.

The device components and batteries/rechargeable batteries must be disposed of in accordance with local regulations (see Chapter 8, DISPOSAL).

Component	Expected life cycle	Usage profile
Control unit	3 to 5 years	
Nebuliser handset (without aerosol head)	12 months	Used at home with two inhalation treatments for a total of 30 minutes' inhalation and one disinfection cycle per day
Aerosol head	6 months	The life cycle of the aerosol head may be as short as 3 months if it is used more often and for longer periods than the values indicated in the profile above. The components must then be replaced correspondingly sooner.
	3 months	

Materials used


- Nebuliser handset:
Polypropylene, thermoplastic elastomers, polyoxymethylene
- Mouthpiece:
Polypropylene, silicone rubber

The eFlow[®]rapid contains no parts made from natural rubber (latex).

Service and manufacturer

Manufacturer:	PARI Pharma GmbH Moosstr. 3 82319 Starnberg Germany
Contact for technical problems or questions regarding the device:	Local service partner Contact: www.eflowrapid.info

2 COMPONENTS

- 1 Control unit
- 2 Nebuliser handset incl. aerosol head (2 parts):
- 3 Aerosol head**
- 4 Medication cap (pre-fitted), consisting of:
 - 4a Cap seal
 - 4b Cap
- 5 Medication reservoir and
- 5a Aerosol chamber (pre-fitted)
- 6 Inspiratory valve
- 7 Mouthpiece with
- 7a Expiratory valve (pre-fitted)

- 8 Connection cord (connection between control unit and nebuliser handset)

- 9 International power adapter

- 9a Interchangeable adapter (3 items)

- Transport case and nebuliser bag
- easyCare cleaning aid for aerosol heads

Check that all components are contained in your package. If they are not, contact your local service partner.

**TouchSpray® Technology made under licence from the Technology Partnership PLC

3 PREPARATION FOR INHALATION

Caution:

Follow the hygiene instructions and ensure that the nebuliser handset is also cleaned and disinfected before it is used for the first time (see Chapter 5, HYGIENIC REPROCESSING).


Connection to the power supply

The device can be operated with rechargeable batteries, non-rechargeable batteries or connected to the mains power supply via the power adapter supplied.

Mains power operation (100 V~ to 240 V~)

Caution:

- Only use the power adapter supplied.
- If the power adapter is damaged or faulty, discontinue use immediately.
- Plug the connector of the power adapter into the control unit:


- Plug the power plug on the power adapter into the socket.

When the adapter is plugged into the socket, the device is ready for use.

Changing the country-specific adapter:

- Disconnect the power adapter from the socket.
- Push the interchangeable adapter upwards from the bottom until it detaches from the power adapter.
- Attach the country-specific adapter (interchangeable adapter, 9a) by inserting and sliding it into place:


- Make sure that the adapter engages in the power adapter housing.


Operation with batteries

Recommended: 1.5 V Mignon AA LR6/Alkaline or equivalent

- Open the battery compartment by lifting the clip:

- Insert the batteries as indicated by the polarity symbols:

- Close the battery compartment by snapping the battery compartment cover onto the control unit.

When the device is running on batteries, the LED on the control unit flashes red at one-second intervals to indicate that power is low. The device switches itself off automatically shortly afterwards.

Batteries vary considerably in terms of storage life and output depending on battery quality. The life cycle of about 90 minutes can only be achieved with the recommended batteries (4x). It is recommended to carry a spare pack of batteries for trips

or leisure activities when connection to mains power is unavailable.

Caution: Risk of accident if operated in a motor vehicle (e.g. car, motor home, boat)

For your own safety, do NOT use the inhalation device while driving.

Park the vehicle and switch the engine off before starting the inhalation treatment.

Operation with rechargeable batteries

The control unit can also be operated with standard rechargeable batteries.

Recommended:

1.2 V Mignon AA rechargeable, at least 2100 mAh or equivalent


Prepare the control unit as for operation with batteries.

Batteries cannot be recharged in the control unit!

Assembling the nebuliser handset


i **Note:**
Check the parts of your nebuliser handset before every application. Replace all broken, misshapen and heavily discoloured parts. Please also follow the instructions for assembly given below. Damaged parts and incorrectly assembled nebuliser handsets can impair the function of the nebuliser handset and thus also affect the treatment.

- To open the medication reservoir (5) and aerosol chamber (5a), gently raise the locking clip on the reservoir.


i **Note:**
The medication reservoir and aerosol chamber are already pre-assembled and must not be separated at the hinged joint (danger of breakage).

- Grasp the engaging hook on the side of the aerosol head (3) - do not touch the membrane in the middle - and insert the two metal connectors in the slots on the medication reservoir:


i **Note:**
The aerosol head must snap into place. Please ensure that the aerosol head is fitted the right way round.

- Place the inspiratory valve (6) flat against the rim of the aerosol chamber (5a):

i **Note:**
Ensure that the valve is inserted in the right direction! The inspiratory valve must lie flat against the aerosol chamber and both valve wings must be flush.


- Press the aerosol chamber and the medication reservoir together until the clip on the reservoir engages in the locking lug of the aerosol chamber:

i **Note:**
If the button closure does not close, open the nebuliser handset by lifting the clip and ensure that the inspiratory valve is seated correctly.


Fit the medication cap (4), consisting of the cap seal with soft sealing collar and cap:

- Hold the cap (4b) with the small opening upwards and push the cap seal (4a) onto the four locking lugs from below until it snaps audibly into place.


i **Note:**
Check that the soft collar is fully extended and without creases.

- Attach the expiratory valve (pre-assembled, 7a) in the slot provided for that purpose in the mouthpiece (7):


- Insert the mouthpiece (7) with the expiratory valve pointing upwards in the connector to the aerosol chamber:


Connecting the nebuliser handset

- Mount the nebuliser adapter of the connection cord (8) on the hinge axis of the aerosol chamber (5a) and move upwards to engage the metal contacts:


- Insert the round plug on the other end of the connection cord into the socket on the front of the control unit (1):

4 INHALATION

Filling and sealing the medication reservoir


- Add the quantity of medication prescribed by your doctor to the medication reservoir:


Caution:

Ensure that the medication is filled no higher than the upper scale graduation mark (6 ml). Do not overfill under any circumstances!

- Close the medication reservoir: Place the cap (4) on the medication reservoir so that the slots in the side of the cap are positioned above the notches in the reservoir. Press gently and turn the cap clockwise as far as it will go:

- **Note:** The closing mechanism is functioning correctly if the cap seal (4a) rises with the turning motion to form a seal.


- Check that all parts are connected tightly and that the medication reservoir is sealed. Otherwise nebulisation may be impaired.

Performing inhalation

- Hold the nebuliser handset in your hand.
- Sit in an upright position and relax. This makes inhaling easier and improves deposition of the medication in the respiratory tract.
- Hold the mouthpiece between your teeth and enclose with your lips. The lips should not touch the blue expiratory valve:

- Press the ON/OFF button on the control unit to start producing the aerosol.

A green LED beside the ON/OFF button lights up and an audible signal (1 beep) is emitted to indicate proper functioning.

- Inhale and exhale as slowly and deeply as possible through the mouthpiece or mask.
The mouthpiece can also remain in the mouth while exhaling. Do not breathe through your nose. You may use a nose clip if your doctor agrees.

i **Note:**
It is normal for aerosol to escape through the expiratory valve in the mouthpiece during exhalation; this is not a malfunction.

Inhalation with mask

Inhalation with the mouthpiece is the most effective form of inhalation, since it results in minimum medication loss on the way to the lung. The PARI SMARTMASK[®] or PARI SMARTMASK[®] Kids (see Chapter 9, SPARE PARTS AND ACCESSORIES) is recommended only for patients who are unable to inhale with a mouthpiece.

- Put the mask on so that it forms an airtight seal over the mouth and nose.
 - Adults: PARI SMARTMASK[®]
 - Children aged 2 and older: PARI SMARTMASK[®] Kids
- Perform the inhalation as described in Section "Performing inhalation" (Chapter INHALATION).

i **Note:**
Follow the instructions for use of the mask!

Pause function

- To activate the Pause function, press the ON/OFF button for about 1 second. Pause mode is activated if the LED flashes green.
- To continue inhalation, press the ON/OFF button for about 1 second once again.

i *Note:*
The Pause function can be used after an operating time of 5 seconds.

If a malfunction occurs, refer to Chapter 7, TROUBLESHOOTING.

The device can be switched off manually at any time during operation by pressing the ON/OFF button again. To confirm that it has been switched off, the device beeps once and the LED lights red.

As long as there is enough liquid in the medication reservoir, the nebulisation can be restarted by pressing the ON/OFF button again.


Display indicators during inhalation

The following indicators appear on the display during inhalation:


Start

During inhalation


Nebulising with battery power


Nebulising with mains power

Activated Pause function

Interrupting nebulisation during battery operation


Interrupting nebulisation during mains operation

End of the inhalation session

The device is switched off automatically when the medication dosage that is available for nebulising has been used up, or after a maximum of 20 minutes.

The exact treatment time depends on the type and quantity of the medication used. This should not be longer than 4 minutes for a 2.5 ml volume of a saline solution. If this time is exceeded significantly, please read Chapter 7, TROUBLESHOOTING and Section "Life cycle" (Chapter IMPORTANT INFORMATION).


When inhalation has been completed successfully, the following indicator appears on the display:

Caution:

Note that about 1 ml of medication cannot be nebulised and remains unused in the medication reservoir; this residue must be discarded. This is intentional and does not represent a malfunction.

The device switches itself off as soon as the membrane of the aerosol head is no longer coated with medication.

Caution:

If the device switches off prematurely, and significantly more than 1 ml of medication is left in the medication reservoir, please read Chapter 7, TROUBLESHOOTING.

Note:

If the device does not switch off automatically, even though the medication dosage that is available for nebulising has been used up, switch the device off manually (see also Chapter 7, TROUBLESHOOTING).

- When the inhalation session has ended, disconnect the power adapter plug from the socket!

5 HYGIENIC REPROCESSING

General

Warning:

- To prevent any health risk, e.g. infection from a contaminated nebuliser handset, it is imperative to comply with the following hygiene instructions.
- Before cleaning, switch the device off and remove the power adapter from the socket.

The eFlow[®] *rapid* is designed for multiple uses. Be aware that different requirements regarding hygienic reprocessing are in effect for different fields of application:

- At home (single patient use): The nebuliser handset (including the aerosol head) must be cleaned immediately after each use and disinfected at the end of the day.
- Hospital/medical practice (possible change in patient): The nebuliser handset (incl. aerosol head) must be cleaned, disinfected and sterilised before every use.

If several different medications are to be inhaled successively as part of the inhalation treatment, the following should also be noted:

All medication residue must be rinsed out of the nebuliser handset with warm, running tap water after every application. This helps to prevent undesirable cross-reactions between medication residues.

After the last application in each session, the nebuliser handset must be cleaned, disinfected at least once a day, and also sterilised if it is to be used for another patient. This is essential to prevent the growth of bacteria.

Additional requirements concerning the necessary hygienic preparations (hand care, handling of medications or inhalation solutions) in high-risk groups (e.g. patients with cystic fibrosis) can be obtained from the respective self-help groups.

Please ensure adequate drying after each cleaning, disinfection and/or sterilisation. Condensation or residual wetness may present a risk due to the growth of bacteria.

Proof of suitability of the nebuliser handset for effective cleaning, disinfection and sterilisation has been obtained through an independent testing laboratory using the recommended procedures. Using the indicated alternatives remains the responsibility of the user.

Please check the parts of the nebuliser handset regularly and replace any defective (broken, misshapen, discoloured) parts. Replace the nebuliser handset at least every year and the aerosol head at least every six months (see Section "Life cycle", Chapter IMPORTANT INFORMATION).

Dismantling the nebuliser handset

Dismantle the nebuliser handset into its individual parts in the following order (see Figure on page 41):

- Detach the nebuliser handset (2) from the connection cord adapter (8):
- Remove the mouthpiece, the mask or the Filter/Valve Set from the nebuliser handset.
- Carefully pull the expiratory valve (7a) out of the slot in the mouthpiece.

Note:

The expiratory valve is held in place by a thin fixing clip so that it cannot get lost. This clip must not be detached from the mouthpiece.


- Take off the medication cap (4) and shake out the medication residue.
- Press the cap seal (4a) downwards and out of the cap (4b).
- Lift the locking clip on the aerosol chamber (5a) so that the nebuliser handset can be opened up.

Note:

The medication reservoir and aerosol chamber are already pre-assembled and must not be separated at the hinged joint (danger of breakage).

- Press lightly on the engaging hooks on the sides of the aerosol head (3) to release it from the retainer on the medication reservoir (5):

- Remove the inspiratory valve (6) from the aerosol chamber.

If the device is being processed for hygienic re-use at home, read the next Section "Hygienic reprocessing at home". If the device is being processed for hygienic re-use in the hospital or a doctor's practice, read Section "Hygienic reprocessing in the hospital and doctor's practice".

Hygienic reprocessing at home

Cleaning

- Place all nebuliser handset parts as well as the aerosol head in warm (approx. 40 °C) tap water with a little dishwashing soap (dosage according to the instructions of the manufacturer) for 5 minutes.

Cleaning the aerosol head


- After immersing in water, rinse the front and back of the aerosol head for approx. 1 minute under warm running tap water (approx. 40 °C).

- If parts are heavily soiled, clean the aerosol head again by shaking in the water in the bowl and then rinsing with running tap water.

Caution:

- Do not put the aerosol head in a microwave oven.
- Do not attempt to clean the aerosol head in a dishwasher.
- Mechanical cleaning of the aerosol head with brushes or abrasives may cause irreparable damage!


Cleaning other parts of the nebuliser handset (not the aerosol head)

- After immersing in water, clean these nebuliser handset parts thoroughly, if necessary using a soft, **clean** toothbrush, and then rinse off with warm (approx. 40 °C), running tap water. You can remove excess water more quickly by shaking all parts.

Cleaning the control unit and connection cord

- Switch the control unit off and disconnect the power cord and the connection cord from the device.

- Clean the control unit housing and the connection cord with a damp cloth.

Caution:

Do not hold the control unit under running water and do not use liquid cleaning agents!

If liquid gets into the control unit, the electronics may be damaged and the control unit may no longer function properly.

If liquid does get into the control unit, contact the local service partner immediately.


Disinfection

At the end of the day, disinfect the dismantled nebuliser handset (including the aerosol head) immediately after cleaning.

Using a standard commercial disinfectant for baby bottles (not a microwave oven)

Disinfection must take at least 15 minutes. For instructions on disinfecting and the amount of water required, it is essential to read the instructions for use of the disin-

factor concerned. Always check carefully that it is clean and fully functional.


Possible alternative:**Boiling in distilled water**

Place the individual parts of your nebuliser handset (incl. aerosol head) in boiling water for at least 5 minutes. Use only distilled water. Ensure that the water level in the pot remains high enough to prevent the nebuliser handset parts from coming into contact with the hot base of the pot.

Drying, storage, transport

A damp environment encourages the growth of bacteria. Therefore, remove the nebuliser handset parts and accessories from the pan or disinfectant as soon as disinfection has finished. The risk of infection is reduced when the parts are dried completely.

- Place all nebuliser handset parts on a dry, clean and absorbent surface and allow to dry completely (at least 4 hours):

**Note:**

Do not dry in damp rooms (e.g. the bathroom).

- Between uses, and in particular during long breaks in treatment, place the nebuliser handset in the accompanying nebuliser bag and keep it in a dry, dust-free environment (e.g. not in the bathroom).
- Reassemble the nebuliser handset just prior to the next treatment. This is to ensure that the silicone seals do not become deformed over time.
- For transport, pack the nebuliser handset together with the control unit, power adapter and connection cord in the transport case provided for this purpose.

Hygienic reprocessing in the hospital and doctor's practice

When using the eFlow[®]rapid for several patients, the nebuliser handset must be cleaned, disinfected and sterilised before each change in patient.

Ensure that only adequately validated device-specific and product-specific procedures are used for cleaning, disinfection and sterilisation, and that the validated parameters are observed in each cycle.

The cleaning and disinfection procedure must be recognised as effective (e.g. in the list of tested and recognised disinfectants and procedures issued by the Robert Koch Institute/DGHM) (Deutsche Gesellschaft für Hygiene und Mikrobiologie - German society for hygiene and microbiology) and already validated in all cases.

If another procedure is used, its effectiveness must be confirmed during validation. CE-approved chemical disinfectants may also be used, provided they satisfy the requirements for material resistance.

In addition, follow the hygiene regulations of the hospital or the medical practice.

Material resistance

The nebuliser handset unit can withstand temperatures up to 121 °C.

Note the following when selecting cleaners/disinfectants:

- Essentially, the aldehydic disinfectant group is suitable for disinfecting this nebuliser handset (hospital/doctor's surgery).
- The materials resistance of this nebuliser handset has not been tested against other groups of cleaning agents and disinfectants.
- When selecting the chemicals to be used, make sure that they are suitable for cleaning and disinfecting medical devices made from the materials specified in

Section "Materials used" (Chapter IMPORTANT INFORMATION).

Cleaning and disinfection

Clean and disinfect immediately after use. Ideally, a mechanical procedure (instrument dishwasher) should be used.

Prepare for the process of cleaning with disinfection as described in Chapter "Dismantling the nebuliser handset".


Recommended procedure: Thermal disinfection

- Place the nebuliser handset components in the instrument dishwasher.
- Select the 93 °C programme (treatment time 10 minutes).

The effectiveness of this procedure has been validated using a G7736 disinfectant made by Miele and using neodisher[®]MA (Chemische Fabrik Dr. Weigert, Hamburg) as a cleaning agent in conjunction with deionised water as a neutraliser.

Recommended procedure: Manual cleaning/disinfection

- Place the nebuliser handset components in a 0.5% solution of Bodedex[®]forte for 5 minutes.
- Then place the components in a 4% solution of Korsolex[®]extra disinfectant solution for 15 minutes.
- Then thoroughly rinse the nebuliser handset components under warm running water and leave them to dry completely on a clean, dry and absorbent surface (for at least 4 hours):

Sterilisation

After cleaning/disinfection, pack the dismantled nebuliser handset in sterilisation packaging (disposable sterilisation packaging, e.g. foil/paper sterilisation bags). The sterilisation packaging must conform to DIN EN 11607 and be suitable for steam sterilisation. Then sterilise the nebuliser handset using the following sterilisation procedure:

Recommended procedure: Steam sterilisation

Validated in accordance with DIN EN ISO 17665-1. Sterilisation temperature: 121 °C (to be maintained for: 30 min.).

The resistance of the aerosol head has been tested for 50 cycles in an autoclave.

Caution:

Using temperatures higher than 121 °C for sterilisation will quickly damage the aerosol head/nebuliser handset.

Storage

Keep the sterilised nebuliser handset in a dry, dust-free place protected from contamination.

6 STORAGE AND TRANSPORTATION CONDITIONS

Ambient conditions for the safe transport and storage of the eFlow[®] *rapid*:


- Temperature: -25 °C to +70 °C
- Relative humidity of ambient air: 0% to 93% (non-condensing)
- Air pressure: 500 hPa to 1060 hPa

Extreme fluctuations in the ambient temperature can cause condensation to form and impair the functioning of the device.

Protect from continuous direct sunlight during storage and transportation.

Do not store the nebuliser handset, control unit and power adapter in damp rooms (e.g., bathrooms) or transport together with damp objects.


If the device is not expected to be used for a prolonged period, the batteries must be removed:

7 TROUBLESHOOTING

Faults in installation or operation (e.g. nebuliser handset not connected, no medication added, power supply interrupted) are displayed approx. 5 to 6 times in flashing mode. The device then switches off. If several causes of malfunction exist, they will be

displayed in the following order each time the device is switched on again. Please follow the step-by-step instructions for correcting such malfunctions (see following list):

Display indicator (flashes alternately)	Possible cause	Solution
	Battery flat.	Insert new batteries or operate with the power adapter connected.
	Cable connection to the control unit or on the nebuliser handset interrupted.	Check cable connections.
	No medication added.	Add medication.
	Inhalation interrupted.	Resume inhalation by restarting the device (press ON/OFF button).
	Battery charge level 50%.	Have new batteries on hand or operate with the power adapter connected.
	The maximum operating time of 20 min. for one application has been exceeded.	Resume inhalation by restarting the device (press ON/OFF button).

The indicator light (LED) always flashes whenever an error mode is indicated in the display, as described on page 58:

Fault	Possible cause/Remedy
The device cannot be activated (no beep, no red or green light signal).	<ul style="list-style-type: none"> - To switch on, press the ON/OFF button for approx. 2 seconds. - Is power plug plugged into mains and the device (for mains operation)? - Are batteries inserted correctly? - Check the charge level of the batteries! - Check the connection between the control unit and the nebuliser handset.
Indicator light (LED) flashes green/red and goes out after a few seconds.	Check the connection between the control unit and the nebuliser handset.
No aerosol seems to be created when the device is activated, or the device switches itself off after a few seconds.	<p>Has the medication reservoir been filled with medication?</p> <p>To check, open the cap on the medication reservoir again, fill with medication if necessary, and close securely again (see section "Connecting the nebuliser handset", Chapter PREPARATION FOR INHALATION).</p>
The indicator light (LED) lights up red during operation.	The batteries are almost flat. Replace them or use the power adapter.
<p>The device switches off during operation, even though there is still sufficient medication in the medication reservoir.</p> <p>E.g.</p> <ul style="list-style-type: none"> - shortly after switching on or - after the automatic shutoff time has elapsed (20 minutes). 	<ul style="list-style-type: none"> - Check the charge level of the batteries! - The nebuliser handset was not being held horizontally. Hold the nebuliser handset horizontally. - Press the ON/OFF button again to resume treatment, so that you can inhale the prescribed quantity of medication.
The device does not switch itself off, even though there is no more than about 1 ml medication left in the medication reservoir.	To switch the device off, press the ON/OFF button.

Fault	Possible cause/Remedy
Overly long inhalation times	<p>If the inhalation time with the same fill volume and the same substance becomes considerably longer, the following points should be checked after the device is switched off.</p> <ul style="list-style-type: none"> - Is the cap of the medication reservoir correctly fitted and securely in place? - Has the aerosol head been prepared for hygienic reprocessing as described in Chapter 5? - Has the aerosol head been physically damaged? - Has the life cycle been exceeded (see section "Life cycle", Chapter IMPORTANT INFORMATION)? <p>After checking these points, a functional test for aerosol heads should be carried out as described in the following.</p> <p>Functional test for aerosol heads in the event of overly long nebulisation times</p> <ul style="list-style-type: none"> • Measure the nebulisation time with 2.5 ml isotonic sodium chloride solution (0.9%) without inhaling. <p>If it takes longer than 4 minutes to nebulise, clean the aerosol head with the easycare cleaning aid (item no. 078G6100). Additional cleaning with the easycare cleaning aid can help to shorten the nebulisation time as the nebulising performance of the aerosol head becomes less than optimal.</p> <p>Please follow the instructions for use of the easycare cleaning aid.</p> <p>If the nebulisation time is longer than 5 minutes when using 2.5 ml isotonic saline solution (0.9%) despite cleaning with the easycare cleaning aid, the aerosol head must be replaced.</p> <p>Note: <i>Nebulisation times for certain medications may differ from the values indicated above for isotonic sodium chloride solution (0.9%).</i></p>

Fault	Possible cause/Remedy
Aerosol continuously escapes in large volumes from the opening slots in the medication reservoir.	Check that the nebuliser handset has been assembled correctly (see Section "Assembling the nebuliser handset", Chapter PREPARATION FOR INHALATION). Check carefully that both wings of the inspiratory valve are lying flat and have not been pressed through.
Little or no aerosol escapes below an ambient temperature of 10 °C (see Section "Ambient conditions", Chapter IMPORTANT INFORMATION).	<ul style="list-style-type: none"> - The ambient temperature should be at least 10 °C. - Wait until the device has stabilised at the lowest ambient temperature.
The medication cap cannot be removed from the medication reservoir after the inhalation treatment.	<p>The negative pressure in the aerosol chamber is too great.</p> <p>Open the nebuliser handset and hold so that the aerosol head points upwards. Grasp the aerosol head by the locking lugs on the side and raise it gently. Air will flow in, indicated by a quiet hissing. The cap can now be removed easily.</p>

If the fault is not eliminated after following these steps, notify the local service partner immediately.

8 DISPOSAL

This product falls within the scope of the European Council Directive on Waste Electrical and Electronic Equipment (WEEE¹⁾) and is included in Product Category 8: Medical Devices for the purposes thereof. Accordingly, live parts

1) Directive 2002/96/EC of the EUROPEAN PARLIAMENT AND THE EUROPEAN COUNCIL dated 27 January 2003 on waste electrical and electronic equipment.

(such as the control unit, the cable, the power adapter and the batteries) must not be disposed of with normal domestic waste. The device components and (rechargeable) batteries must be disposed of in accordance with local regulations.

Recycling helps to reduce the consumption of raw materials and protect the environment.

9 SPARE PARTS AND ACCESSORIES

Description	Part no.
Nebuliser handset (incl. aerosol head)	678G8222
Aerosol head	678B2620
Mouthpiece	078B3600
PARI universal mouthpiece (alternative)	022E3050
AC power supply international (100-240 V~, 50/60 Hz)	078B7106
Connection cord	178G6009
Carrying case	078E8001
PARI Filter/Valve Set	041G0500
PARI SMARTMASK [®] (adults)	041G0730
PARI SMARTMASK [®] Kids (children aged 2 and older)	078G5000
easycare cleaning aid for aerosol heads	078G6100

10 TECHNICAL DATA

General

Electrical connection:	Power adapter (REF 078B7106, Type FW7555M/12): Input: 100 V - 240 V~, 50/60 Hz Output: 12 V ---
Operation with rechargeable or disposable batteries:	4 x 1.2 V (rechargeable batteries) 4 x 1.5 V (disposable batteries)
Nebuliser handset weight:	approx. 55 g
Weight eFlow [®] rapid (incl. batteries):	approx. 300 g
Nebuliser handset dimensions (W x H x D):	5.0 cm x 6.3 cm x 14.5 cm
Control unit housing dimensions:	H 4.0 cm, Ø 11.6 cm
Min. fill volume:	2.0 ml
Max. fill volume:	6.0 ml
Aerosol output ^{a)}	0.5 ml ^{b)}
Aerosol output rate ^{a)}	0.33 ml/min ^{b)}
MMAD (mass median aerodynamic diameter) ^{a)}	4.1 µm ^{b)}

a) At 2 ml fill volume, compliant with Annex CC of EN 13544-1

b) May vary depending on medication and the aerosol head used.

Classification according to DIN EN 60601-1

Type of electric shock protection:	Protection class II
Degree of protection from electric shock of part used:	Type BF
Degree of protection against ingress of water as per IEC 529 (IP rating):	No protection
Degree of protection when used in the presence of flammable mixtures of anaesthetics with air or with oxygen or nitrous oxide:	No protection
Operating mode:	Continuous operation

Information concerning electromagnetic compatibility

Electrical medical equipment is subject to special precautionary measures with regard to electromagnetic compatibility (EMC). It must only be installed and operated in accordance with the EMC notes contained at the end of these instructions for use.

Portable and mobile RF-communication devices can affect medical electrical equipment.

Using accessories, a different power adapter or different power cords than those specified, with the exception of the power adapter and power cords converters and power cords that the manufacturer of the medical electrical device sells as spare parts for internal components, can result in higher emission levels or lower the resistance to interference of the device.

The device must not be placed directly beside or on top of other devices for operation. If the medical electrical device must be placed beside or on top of other devices to operate it, it should be monitored constantly to ensure that it is operating properly in the arrangement used.

Detailed information about the electromagnetic compatibility of the eFlow[®]*rapid* is included in an attachment to these instructions for use.

11 EXPLANATION OF SYMBOLS

Follow the instructions for use (background = blue, symbol = white)

Protect device from moisture

Alternating current

Direct current

Protection class of the part used: Type BF

Protection class II appliance

CE-marking for medical devices with identification number of the notified body

The medical device was distributed commercially after 13 August 2005.
This product must not be disposed of with normal domestic waste. The symbol of a crossed-through refuse bin indicates that separate collection is required.

Manufacturer

We reserve the right to make technical changes.

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Electromagnetic compatibility - Guidance and manufacturer's declaration DIN EN 60601-1-2:2007 (IEC 60601-1-2:2007)

Medical electrical devices are subject to special precautionary measures in particular regarding the EMV with the installation and the operation.
Portable and mobile HF-communication devices e.g. mobile phone can affect medical electrical devices.

A use of other accessories and lines than the indicated, can lead to an increased sending or a reduced noise immunity of the equipment. The equipment has to be operated exclusively with original accessories.

The device should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the device should be observed to verify normal operation in the configuration in which it will be used.

Electromagnetic environment

The eFlow[®]rapid is intended for use in the electromagnetic environment specified below. The customer or the user of the eFlow[®]rapid should assure that it is used in such environment.

Electromagnetic emissions

Emission test	Compliance	Electromagnetic environment – Guidance
RF emissions CISPR 11	Group 1	The device uses RF energy only for its internal function. Therefore, its RF-emission is very low and not likely to cause any interference nearby electronic equipment.
	Class B	The eFlow [®] rapid is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations / flicker emissions IEC 61000-3-3	Passed	

Table 1 of DIN EN 60601-1-2:2007 (IEC 60601-1-2:2007)

Electromagnetic immunity

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient /burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input / output lines	±2 kV for power supply lines ±1 kV for input / output lines	Mains power quality should be similar to that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV common mode ±2 kV differential mode	±1 kV common mode ±2 kV differential mode	Mains power quality should be similar to that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% $U_T^a)$ (>95% dip of U_T) for 1/2 cycle 40% U_T (60% dip of U_T) for 5 cycles 70% U_T (30% dip of U_T) for 25 cycles <5% U_T (>95% dip of U_T) for 5 s	<5% U_T (>95% dip of U_T) for 1/2 cycle 40% U_T (60% dip of U_T) for 5 cycles 70% U_T (30% dip of U_T) for 25 cycles <5% U_T (>95% dip of U_T) for 5 s	Mains power quality should be similar to that of a typical commercial or hospital environment. When the user of the eFlow [®] rapid continued function also calls in the event of disruption of supply, it is recommended the eFlow [®] rapid from an uninterruptible power supply or a battery.
Power frequency (50 Hz/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Conducted RF IEC 61000-4-6	3 V_{eff} 150 kHz to 80 MHz	10 V_{eff} 150 kHz to 80 MHz	Portable and mobile RF communication equipment should be used no closer to any part of the eFlow [®] rapid, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.

a) U_T is the mains voltage prior to application of the test level.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – Guidance
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2,5 GHz	10 V/m	<p>Recommended separation distance:</p> $d = 0.35 * \text{SQRT} (P)$ $d = 0.35 * \text{SQRT} (P) \text{ 80 MHz to 800 MHz}$ $d = 0.70 * \text{SQRT} (P) \text{ 800 MHz to 2.5 GHz}$ <p>Where „P“ is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and „d“ is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey^{a)}, should be less than the compliance level in each frequency range^{b)}.</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 

Table 2 and 4 of DIN EN 60601-1-2:2007 (IEC 60601-1-2:2007)

- a) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy.
- To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the eFlow[®]rapid is used exceeds the applicable RF compliance level above, the eFlow[®]rapid should be observed to verify normal operation.
- If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the eFlow[®]rapid.
- b) Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 10 V/m.

Note 1: At 80 Hz and 800 MHz, the higher frequency range applies.

Note 2: This guidance may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Separation distances

The eFlow[®]rapid is intended for use in the electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the eFlow[®]rapid can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communication equipment (transmitters) and the eFlow[®]rapid as recommended below, according to the maximum output power of the communications equipment.

Recommended separation distances for portable and mobile RF communication equipment and the eFlow[®]rapid

Rated maximum output power of transmitter [W]	Separation distance according to frequency of transmitter [m]		
	150 kHz to 80 MHz $d = 0.35 * \text{SQRT} (P)$	80 MHz to 800 MHz $d = 0.35 * \text{SQRT} (P)$	800 MHz to 2.5 GHz $d = 0.70 * \text{SQRT} (P)$
0.01	0.04	0.04	0.07
0.1	0.11	0.11	0.22
1	0.35	0.35	0.70
10	1.11	1.11	2.21
100	3.50	3.50	7.00

Table 6 of DIN EN 60601-1-2:2007 (IEC 60601-1-2:2007)

For transmitters rated at a maximum output power not listed above, the recommended separation distance „d“ in meters [m] can be determined using the equation applicable to the frequency of the transmitter, where „P“ is the maximum output power rating of the transmitter in watts [W] according to the transmitter manufacturer.

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.

Note 2: These guidances may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

CE 0123

PARI Pharma
Advancing Aerosol Therapies

PARI Pharma GmbH
Moosstrasse 3 • 82319 Starnberg • Germany
E-Mail: info@pari.de • www.eflowrapid.info