## Warnings

- Put small part of tape on your skin for testing before applying the whole part
- Slightly "fold in" your lips before applying the tape
- Moisten your lips slightly before applying the tape
- If you have dry lips, it could indicate dehydration and you may need to drink more water
- Remove the tape carefully when taking it off in the morning
- Use your tongue to first moisten your lips and then push out the tape from within
- If you have irritating redness or swollen, please stop using and go for a doctor
- Don't use if you have sensitive skin and get sore lips
- Don't use if you have a cold
- take off from above to under

This product is not intended for :

- children under the age of 6
- those with breathing disorders
- heart conditions or nasal congestion
- risk of vomiting
- with the use of alcohol or sedatives
- those unable to remove the tape on their own
- not intended to diagnose, treat, cure, or prevent any disease and/or health condition.

## Instructions

1. clean the lips and skin around
2. dry the lips and skin around
3. remove the sticker from the sheet
4. gently press your lips together
5. stick the tape on the mouth and press the edges well for 30 seconds

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