

# Nokia Steel HR™

Activity - Sleep - Heart Rate  
Installation and Operating Instructions



# Table of Contents

---

Minimum Requirements.....	4
Internet access .....	4
Smartphone or device .....	4
OS version.....	4
Withings Health Mate® app version .....	4
Description .....	5
Overview .....	5
Box Contents .....	5
Watch screens description .....	5
Setting Up the Withings Steel HR™ .....	7
Installing the Withings Health Mate® App.....	7
Installing my Withings Steel HR™ .....	7
Setting the Time .....	20
Customizing the settings of my Withings Steel HR™ .....	21
Changing the screen order.....	21
Enabling notifications.....	22
Tracking my activity.....	24
Overview .....	24
Checking my progress .....	24
Changing my activity goal .....	25
Tracking my swimming sessions .....	26
Tracking my sleep .....	27
Overview .....	27
Retrieving my sleep data.....	27
Setting an Alarm .....	28
Tracking my heart rate .....	31
Overview .....	31
Wearing my Withings Steel HR™ .....	31
Retrieving my heart rate data .....	31
Launching the continuous heart rate mode .....	32

Getting more active.....	33
Competing with my friends.....	33
Competing with myself .....	34
Syncing your Data.....	36
Dissociating your Withings Steel HR™.....	37
Cleaning and Maintenance .....	39
Cleaning my Withings Steel HR™.....	39
Charging the battery .....	39
Replacing the Wristband .....	39
Calibrating the Hands.....	41
Updating the Firmware Version.....	43
User Guide Copyright.....	44
Legal Notice .....	44
Personal Data.....	44
Safety Instructions .....	44
Use and Storage .....	44
Heart Rate Measurement .....	45
Water Resistance .....	45
Service and Maintenance .....	45
Specifications .....	46
Document Release Overview .....	48
Warranty.....	49
Regulatory Statements.....	50

## Important Notice



**By using your Nokia Steel HR™ you expressly agree to the Withings Services Terms and Conditions available on our [website](#).**

## Minimum Requirements

---

### Internet access

Access to the Internet is required to:

- Download the Health Mate® app,
- Set up your Nokia Steel HR™,
- Upgrade the firmware.

For more details, refer to [“Specifications” on page 48](#).

### Smartphone or device

An iOS device or an Android device (with operational Bluetooth Low Energy and Wi-Fi or 3G/4G) is required to:

- Set up your Nokia Steel HR™,
- Retrieve data from your Nokia Steel HR™,
- Interact with your Nokia Steel HR™,
- Interact with people.

For more details, refer to [“Specifications” on page 48](#).

### OS version

You need iOS 8 (or higher) or Android 6.0 (or higher) on your device to install and run the Health Mate® App.

For any device, BLE (Bluetooth Low Energy) compatibility is required to install your Nokia Steel HR™.

### Health Mate® app version

You always need the latest version of the Health Mate® app to have all the latest features available.

# Description

---

## Overview

The Nokia Steel HR™ does it all and still manages to look cool and stylish while doing it. Not only is it able to track your activity and your sleep, it is also water-resistant up to 5ATM, allowing you to track your swimming sessions. However what sets it apart is that it's able

to track your heart rate during your workout sessions, or even continuously if you choose to. Last but not least, it allows you to set a vibration-based alarm that wakes you up effectively without waking up your entire household.

You can easily check your activity by taking a look at the activity hand which goes from 0% up to 100% when you reach your daily step goal. By default, the daily step goal is 10,000 steps. You can, however, change that goal anytime you want through the Health Mate® app. Whenever you go above your daily step goal, the activity hand goes back to the 0% mark and starts moving toward the 100% mark all over again.

An OLED screen also allows you to see additional information directly from your watch: heart rate, calories burned, distance, and alarm time. Any additional information you might need is available from the Health Mate® app, along with the full history of the data collected by your Nokia® products.

## Box Contents



(1) Nokia Steel HR™



(2) Charging cable



(3) Quick start guide

## Watch screens description

The OLED display of your Nokia Steel HR™ allows you to display a number of screens showing additional information. Following are the screens available:

- Date and time: Displays the date and time. This screen cannot be disabled.
- Heart Rate: Displays your latest heart rate measurement.

- Steps: Displays your current number of steps for the day.
- Distance: Displays your current distance for the day.
- Active Calories: Displays the number of calories you've actively burned for the day.
- Alarm: Displays the time of your next alarm.

# Setting Up the Nokia Steel HR™

---

## Installing the Health Mate® App

If the Health Mate® app is not already installed on your device, perform the following steps:

1. Type [go.withings.com](http://go.withings.com) in your device's web browser.
2. Tap **Download on the App Store**.
3. Tap **FREE**.  
The installation process starts.

## Installing my Nokia Steel HR™

### I'm new to Nokia®

If you don't already have a Nokia® account, you have to create one in order to set up and use your Nokia Steel HR™.

If you already have an account, refer to “[I already have a Nokia® account](#)” on page 14.


To install your Nokia Steel HR™, perform the following steps:

1. Open the Health Mate® app.

2. Tap **Join Withings**.



3. Tap **I'm setting up my device**.



4. Tap **Activité Watches**:



5. Tap **Steel HR**.



6. Tap **Install now**.




7. Press the button on the side of your Nokia Steel HR™ and tap **Next**.



8. Tap **Install now**.



9. Enter the code displayed on your watch and tap **Next**.



10. Tap **Next**.



11. Enter your email address, password, and then confirm your password. Tap **Next**.

12. Enter your first name, surname (last name) and birthday. Tap **Next**.

13. Select your gender and enter your current height and weight. Tap **Create**.

14. When prompted to, slide your finger on the wheel displayed in the app until the minute hand reaches the 0 mark. Tap **Next**.



15. Slide your finger on the wheel displayed in the app until the hour hand reaches the 0 mark. Tap **Next**.



16. Slide your finger on the wheel displayed in the app until the activity hand reaches the 0 mark. Tap **Next**.




Note that the activity hand only moves in one direction.

17. Tap **OK**.

18. Tap **Next** to watch the tutorials.



19. Tap **Enable notifications** if you want to receive notifications on your watch and tap **Next**.



Go to [“Enabling notifications”](#) on page 23 for more information on notifications.

20. Tap **Close**.



21. Tap and slide the daily step goal up or down to select the goal you want.  
The range for your goal is between 3,000 and 15,000 steps a day.



22. Tap **OK**.

### **I already have a Nokia® account**

If you already have a Nokia® account, you can get started on the installation process.  
To do so, perform the following steps:

1. Open the Health Mate® app.

2. Tap **My devices**.
3. Tap **+**.



4. Tap **Activité Watches**:



5. Tap **Steel HR**.



6. Tap **Install now**.




7. Press the button on the side of your Nokia Steel HR™ and tap **Next**.



8. Tap **Install now**.



9. Enter the code displayed on your watch and tap **Next**.



10. Tap **Next**.



11. Select the user who will be using the Nokia Steel HR™.



12. When prompted to, slide your finger on the wheel displayed in the app until the minute hand reaches the 0 mark. Tap **Next**.



13. Slide your finger on the wheel displayed in the app until the hour hand reaches the 0 mark. Tap **Next**.



14. Slide your finger on the wheel displayed in the app until the activity hand reaches the 0 mark. Tap **Next**.




Note that the activity hand only moves in one direction.

15. Tap **OK**.

16. Tap **Next** to watch the tutorials.



17. Tap **Enable notifications** if you want to receive notifications on your watch and tap **Next**.



Go to [“Enabling notifications”](#) on page 23 for more information on notifications.

18. Tap **Close**.



19. Tap and slide the daily step goal up or down to select the goal you want.

The range for your goal is between 3,000 and 15,000 steps a day.



20. Tap **OK**.

## Setting the Time

The Nokia Steel HR™ is set up using your mobile device, which allows it to be automatically set to the current time, based on the time zone you are in. If you ever have to change time zones while wearing your Nokia Steel HR™, all you have to do is open the Health Mate® app to sync your watch with your device and the correct time will be displayed.

Note that you won't have to worry about daylight savings either. In this case as well, syncing your Nokia Steel HR™ with your device will be enough to update to the correct time.

## Wearing my Nokia Steel HR™

The Nokia Steel HR™ should be worn about a centimeter below the wristbone. To ensure that you get accurate heart rate readings, the wristband should be tight around your wrist in a way that remains comfortable.

# Customizing the settings of my Nokia Steel HR™


## Changing the screen order

You can change the order of the screens of your Nokia Steel HR™, or even remove some of them if you wish. To do so, perform the following steps:

1. Open the Health Mate® app.
2. Go to **My devices**.
3. Tap **Steel HR**.




4. Tap **Screens customization**.




5. Select the screens you want to display:

- Tap the blue circles on the left to hide a screen.
- Tap the empty circles on the left to display a screen.



6. Tap and hold the icon on the right to drag a screen up or down.




7. Tap the back button to confirm.

## Enabling notifications


If you want to get notified on your Nokia Steel HR™ whenever you have an incoming call or when you receive a text message, you can enable the notifications option. To do so, perform the following steps:

1. Open the Health Mate® app.
2. Go to **My devices**.


3. Tap **Steel HR**.



4. Tap **Notifications**.



5. Tap **Enable notifications** and select the notifications you want to enable.



6. Tap the back button to confirm.

# Tracking my activity

## Overview

Your Nokia Steel HR™ allows you to track your activity: number of steps, run tracking and swim tracking. Whether you are walking, running or swimming, you don't have to do anything to start tracking it. All you have to do is wear your Nokia Steel HR™, it will do the rest.



(1) Visual representation of the number of steps: each bar corresponds to a 30-minute period

(2) Number of steps during the day

(3) Percentage of the daily step goal achieved

(4) Estimated distance walked during the day

(5) Estimated calories burned during the day (active calories and total number of calories)

The integrated accelerometer allows your Nokia Steel HR™ to automatically count your number of steps and to recognize when you are running.

Your running activity, your distance and your calories are calculated thanks to an algorithm based on your number of steps and your personal data (age, height and weight).

## Checking my progress

In order to check your progress in terms of number of steps, all you have to do is have a look at the activity hand. Depending on its position, you can tell what percentage of your daily step goal you have achieved. Once you reach the 100% mark, your watch vibrates and the cycle starts all over again, from 0% to 100%.


To get more information (such as the exact number of steps, the distance, or the calories you have burned), simply press the button located on the right-hand side. The different screens of your Nokia Steel HR™ will be displayed one after the other.

To get even more details, simply log in to the Health Mate® app. The app offers many graphs that allow you to easily see your history and your progress from the very first time you use it.

## Changing my activity goal


This feature allows you to set a custom daily step goal. To change your daily step goal, perform the following steps:

1. Open the Health Mate® app.
2. Go to your **Timeline**.
3. Tap **Goal**.



If this screen is not displayed, swipe left or right in the top widget (colored part at the top of the screen) until you find it.

4. Tap **Goal** once more.



5. Tap and slide the daily step goal up or down to select the goal you want.



6. Tap **OK**.
7. Sync your Nokia Steel HR™ by opening the Health Mate® app and refreshing your Timeline by sliding your finger from top to bottom.

## Tracking my swimming sessions

Your Nokia Steel HR™ allows you to track your swimming sessions. All you have to do is wear your Nokia Steel HR™ while swimming. It will automatically detect that you're swimming and track your swimming sessions without any action required on your part.

Once you're done swimming, the Health Mate® app will let you know how long your session was and how many calories you've burned.

# Tracking my sleep

## Overview

Your Nokia Steel HR™ lets you track your sleep without requiring any action on your part. All you have to do is wear your Nokia Steel HR™ when you go to sleep: it will do the rest on its own.

If you do not want to track your sleep, simply don't wear your Nokia Steel HR™ to sleep.

## Retrieving my sleep data

You can retrieve the sleep data collected by your Nokia Steel HR™ in your **Timeline** in the Health Mate® app.



- |  |                                  |                                     |
|--|----------------------------------|-------------------------------------|
| (1) Breakdown of your night:                 | (2) Time spent laying            | (5) Resting heart rate              |
| - Time spent awake (orange)                  | (3) Total time spent sleeping    | (6) Time it took you to fall asleep |
| - Duration of light sleep cycle (light blue) | (4) Percentage of the sleep goal | (7) Number of wake-ups              |
| - Duration of deep sleep cycle (dark blue)   |                                  | (8) Heart rate graph                |


The Health Mate® app keeps your entire history, allowing you to easily track how your sleep data evolves over time.

# Setting an Alarm

## Enabling the Alarm


In order to enable and set an alarm, perform the following steps:

1. Open the Health Mate® app.
2. Go to your **Timeline**.
3. Tap **Set your alarm**.




If this screen is not displayed, swipe left or right in the top widget (colored part at the top of the screen) until you find it.

4. Put your Nokia Steel HR™ close enough to your mobile device.
5. Tap the **Alarm** toggle to enable the alarm.




6. Tap and slide the wake up time up or down to select the alarm time.



You can only select multiples of five minutes.

7. If necessary, select a smart wake-up time when setting an alarm.  
Your watch will then find the best moment to wake you up during that time period to make for a smoother wake-up experience.
8. Tap **OK**.
9. Tap **Repeat on** if you want to select on which days the alarm should be triggered.




10. Tap the back button to confirm.

## Disabling the Alarm


1. In the Health Mate® app, go to your **Timeline**.

2. Tap your alarm time.



If this screen is not displayed, swipe left or right in the top widget (colored part at the top of the screen) until you find it.

3. Put your Nokia Steel HR™ close enough to your mobile device.
4. Tap the **Alarm** toggle to disable the alarm.



5. Tap the back arrow at the top left of the screen to go back to your **Timeline**.  
Your alarm is disabled.

## Stopping the Alarm

You can stop the alarm by pressing the button of your Nokia Steel HR™. It will stop right away instead of vibrating eight times before stopping on its own.

# Tracking my heart rate

---

## Overview


Your Nokia Steel HR™ can help you track your heart rate. By default, it takes a heart rate measurement every ten minutes, allowing you to see how your heart rate fluctuates during the day.

You also have the option of using the continuous heart rate mode to get your continuous heart rate, which is especially interesting during a workout session.

In this case, a measurement is taken every second. Note that the continuous heart rate mode also enables a timer so that you know exactly how long you've been working out.

## Wearing my Nokia Steel HR™

To make sure your heart rate can be effectively tracked by your Nokia Steel HR™, make sure the wristband of your Nokia Steel HR™ is tight enough while still remaining comfortable. It should also be worn about a centimeter (or half an inch) away from your wrist bone.



## Retrieving my heart rate data

You can see your latest heart rate measurement directly from your Nokia Steel HR™. To do so, simply press and release the button to go through the screens until the heart rate screen is displayed.

Note that going to the heart rate screen will trigger a new heart rate measurement.

If you want to see the full history of your heart rate measurement, you can log in to the Health Mate® app.

## Launching the continuous heart rate mode

You can easily launch the continuous heart rate mode before a workout session. To do so, perform the following steps:

1. Press the button of your Nokia Steel HR™ for about a second until it vibrates.  
The continuous heart rate mode is on
2. Press the button until it vibrates again when you want to stop the continuous heart rate mode.

Note that only the heart rate screen appears on the display of your Nokia Steel HR™ when the continuous heart rate mode is on.

## Competing with my friends

### Overview


Your Nokia Steel HR™ and the Health Mate® app help you stay motivated by bringing a social dimension to the equation. You can add your friends and family to your leaderboard and compete with them to be the one who walks the most every week. It doesn't matter if they have a Nokia® product or not: you can add anyone as long as they have the Health Mate® app on their smartphone, thanks to the in-app tracking that can act as an activity tracker.

This feature is only available to the main user of the account.

### Adding someone to my leaderboard

You can invite people to join your leaderboard. To do so, perform the following steps:


1. In the Health Mate® app, go to **Leaderboard**.
2. Tap ....
3. Tap **Challenge a friend**.
4. Select how you want to send your invitation.  
Several options are available, such as message, email or social media websites.



## Interacting with someone in my leaderboard

Once you have friends in your leaderboard, you can interact with them by tapping their name and selecting one of the three actions available:

- **Cheer** to encourage them,
- **Taunt** to show off,
- **Message** to chat with them.




Adding a social component to your activity tracking is a great way to stay motivated and push your limits.

## Competing with myself

### Badges

The Health Mate® app offers a system of badges, virtual rewards that you unlock after certain achievements. There are two kinds of badges:


- daily steps badges whenever you reach a certain number of steps on a given day,
- cumulative distance badges when you reach an overall distance.



Badges are a fun way to stay motivated and to be rewarded for your efforts.

## Insights

Once you have been using your Nokia Steel HR™ for a few days, the Health Mate® app will provide you with Insights. This feature looks at your activity and lets you know when you are more active than usual. Whenever you are less active than usual, the Insights encourage you to keep going so that you can reach your health goal.



**It is not possible to disable the Insights.**

# Syncing your Data

---

Thanks to background synchronization, you should never have to worry about syncing your data, as long as your Bluetooth is always enabled. Background sync is automatically triggered when:

- You have taken more than 1,000 steps since the last sync.
- You have reached your daily step goal.
- More than 6 hours have passed since the last sync took place.
- You change time zones or when daylight savings start.

You can also sync your data manually. All you have to do to sync your data is open the Health Mate® app. When you do, your Nokia Steel HR™ automatically syncs your data with the app. Make sure your Nokia Steel HR™ is close enough to your mobile device and the Bluetooth is activated on your device or you will not be able to synchronize your data.



**Your Nokia Steel HR™ has a limited memory. Make sure you open the app regularly so that you do not lose any of your data.**

# Dissociating your Nokia Steel HR™


Dissociating your Nokia Steel HR™ allows you to remove it from your account. Once it is no longer linked to your account, it can be installed on another account or on the profile of another user of your account.




**Dissociating your Nokia Steel HR™ will not remove any of your data from the Health Mate® app.**

To dissociate your Nokia Steel HR™, perform the following steps:


1. Press the button of your Nokia Steel HR™ for 5 to 6 seconds.
2. Press the button three times until “Factory reset” is displayed.
3. Press the button for one second until “Reset in 5 clicks” is displayed.
4. Press the button five times to effectively reset your watch.
5. Go to the Bluetooth settings of your mobile device and “forget” your Nokia Steel HR™.
6. Open the Health Mate® app.
7. Go to **My devices**.
8. Tap **Steel HR**.



9. Tap **Dissociate this product.**



10. Tap **Remove from my account.**



## Cleaning my Nokia Steel HR™

1. Use a lint-free cloth moistened with water to clean the glass and casing of your Nokia Steel HR™.
2. Run the silicone band under water and rub it with hypoallergenic soap to clean it.
3. Dry the silicone wristband with a soft cloth.

## Charging the battery

Your Nokia Steel HR™ can operate for up to 25 days on a single charge, depending on your use of the heart rate mode. Using the continuous heart rate mode all the time instead of just during activities will deplete the battery much faster.

To charge your Nokia Steel HR™, perform the following steps:

1. Connect the USB end of the charging cable to a power source .
2. Position the other end of the charging cable on the back of your Nokia Steel HR™ in such a way that the magnets align.
3. Leave your Nokia Steel HR™ to charge for up to two hours.  
You can check the percentage of the battery on the display of the watch while it's charging.

Note that if you don't charge your Nokia Steel HR™ right away, the power reserve mode will allow you to keep using it for 20 additional days. However, you will only be able to read the time and track your activity in this mode. You will need to charge your Nokia Steel HR™ to use the heart rate mode.

## Replacing the Wristband

You can easily replace the wristband of your Nokia Steel HR™ thanks to the quick release spring bars with which it is attached. The 36-mm model requires wristbands that are 18-mm wide, while the 40-mm model requires wristbands that are 20-mm wide.

To change the wristband, perform the following steps:

1. Turn your Nokia Steel HR™ over.
2. Locate the knob of each spring bar.

3. Pull the knob of the first spring bar to the left side and remove the wristband.



4. Insert one side of the spring bar of the new wristband.
5. Put the wristband in place while pulling the knob to the left side.




6. Release the knob.  
The wristband is now properly attached to your Nokia Steel HR™.
7. Repeat the previous steps for the other side.


## Calibrating the Hands

If you ever feel that the time displayed on your Nokia Steel HR™ is incorrect, you may need to recalibrate the hands. To do so, perform the following steps:

1. In the Health Mate® app, go to **My devices**.
2. Tap your Nokia Steel HR™.



3. Tap **Set the time**.



4. Tap **Start calibration** when the app is done updating your Nokia Steel HR™.



5. When prompted to, slide your finger on the wheel displayed in the app until the minute hand reaches the 0 mark. Tap **Next**.



6. Slide your finger on the wheel displayed in the app until the hour hand reaches the 0 mark. Tap **Next**.



- Slide your finger on the wheel displayed in the app until the activity hand reaches the 0 mark. Tap **Next**.



Note that the activity hand only moves in one direction.

- Tap **Done**.  
Your Nokia Steel HR™ is now properly calibrated and the hands move to their correct position to display the current time.

## Updating the Firmware Version

We regularly offer new firmware versions so that you can benefit from the latest features we add to our products. When a new firmware version is available, the Health Mate® app usually prompts you to install it through your **Timeline**.




**Make sure you sync your Nokia Steel HR™ before updating the firmware version. Failing to do so will result in a loss of data.**

However, you can check if updates are available and force their installation whenever you want. To do so, perform the following steps:

- In the Health Mate® app, go to **My devices**.
- Tap your Nokia Steel HR™.

3. Tap **Check for update**.



4. Put your Nokia Steel HR™ close enough to your mobile device.  
If your Nokia Steel HR™ is not up to date, the Health Mate® app will prompt you to install a new firmware version.

5. Tap **Upgrade**.

6. Tap **Close** when the upgrade is complete.

**Note:** You can check which firmware version your Nokia Steel HR™ uses in **My devices > Steel HR**.

# User Guide Copyright

---

This User Guide is protected by Intellectual Property laws and copyright. Any reproduction, modification, representation and/or publication, without Nokia® prior approval is strictly prohibited. You may print this User Guide for your personal use exclusively. For any question please contact Withings at : [contact.withings.com](https://contact.withings.com).

## Legal Notice

---



**By using your Nokia Steel HR™ you expressly agree to the Withings Services Terms and Conditions available on our website.**

**This product is not a medical device and should not be used to diagnose or treat any medical condition.**

## Personal Data

- Make sure you have read our privacy policy, which you can find on our [website](#).
- Make sure your Nokia® account password is secure enough to restrict access to your account. It should be at least eight-character long, have mixed case, and use a combination of alphanumeric and special characters.
- We recommend that you protect access to the Health Mate® app with an additional passcode and/or Touch ID. You can do so from the app by going to **Settings** in the app and activating **Touch ID and Passcode**.
- We recommend that you export your data on a regular basis to local storage on your home computer. You can do so from the Health Mate® [web interface](#) by clicking **Measurement table**, selecting the tab corresponding to the data you want to export and clicking **Download**.

## Safety Instructions

---

### Use and Storage

- Use the product at temperatures between -10°C and 45°C (14°F and 113°F).
- Product may be stored in a clean, dry location between -20°C and 85°C (4°F and 185°F) when not in use.
- Avoid any extreme temperature changes.
- Avoid any kind of shocks, be they physical, thermal or otherwise. They may damage internal and external components of your watch.
- Keep the glass case of your Nokia Steel HR™ away from any sharp or hard objects to prevent any scratches.

- Keep your Nokia Steel HR™ and its accessories out of reach from children, as it contains small parts and may present a choking hazard.

## Heart Rate Measurement

- Skin perfusion, or how much blood is flowing through the skin, varies among individuals and according to the environment. In certain circumstances, the Withings Steel HR may not be able to provide an accurate heart rate measurement.
- The technology used by the Withings Steel HR is based on the light absorbed by your skin and blood. Any permanent or temporary changes to your skin, such as tattoos, sweat, or even hair on the wrist may affect the measurements.
- Environmental factors may affect the ability of the Withings Steel HR to take an accurate measurement. For example, water may obstruct or modify the light dispersion, and cold weather may impact skin perfusion.
- The Nokia Steel HR™ needs to stay securely in place to provide a heart rate measurement. Activities that cause force to the wrist, such as boxing or volleyball, may not allow the Withings Steel HR to provide an accurate heart rate measurement.

## Water Resistance

- Your Nokia Steel HR™ can withstand a pressure similar to that of a 50 m (165 ft) depth. You can, therefore, use it when dishwashing, splashing or swimming. However, it can't be used for diving or snorkeling.
- Do not press the button of your watch when it is under water
- Avoid long exposure to salt water. Make sure you rinse your watch with fresh water after spending time swimming in salted water.
- Avoid exposure to any chemical products such as solvents, detergents, perfumes or cosmetics. They may damage the wristband or the casing of the watch.
- Avoid extreme atmospheres (e.g., saunas or steam rooms).

## Service and Maintenance

- Do not attempt to repair or modify your Nokia Steel HR™ yourself. Service and maintenance not described in this User Guide should be performed by professional technicians.
- Exposure to high or low temperatures may shorten battery life.
- The aging of the materials may affect the water resistance of the watch. Water resistance cannot be permanently guaranteed.

# Specifications

---

## Nokia Steel HR™ (36 mm)

- Stainless steel case
- Mineral glass
- Watch diameter: 36.3 mm diameter (1.4 in.)
- Watch thickness: 12.5 mm (0.5 in.)
- Weight: 39 g (1.38 oz)
- Wristbands: Fits wrists up to 200 mm. (7.9 in.)
- Water resistant up to 50 m (165 ft)

## Nokia Steel HR™ (40 mm)

- Stainless steel case
- Mineral glass
- Watch diameter: 39.5 mm (1.6 in) diameter
- Watch thickness: 12.5 mm (0.5 in.)
- Weight: 49 g (1.72 oz)
- Wristbands: Fits wrists up to 230 mm. (9.1 in.)
- Water resistant up to 50 m (165 ft)

## Display

- Analog dial displaying hours and minutes
- Subdial displaying your progress toward your daily step goal
- OLED screen displaying additional information

## Metrics and Features

- Steps (distance and calories)
- Running (distance and calories)
- Swimming (session duration and calories)
- Heart rate
- Sleep
- Silent alarm
- Timer

## Battery

- Up to 4 days in continuous heart rate mode
- Up to 25 days in normal use
- Power reserve mode of an extra 20 days
- Charging time of two hours

## Sensors

- Day and night motion sensor
- High precision MEMS 3-axis accelerometer
- Low power consumption

## Connectivity

- Bluetooth Low Energy (Bluetooth Smart Ready)

## Compatible Devices

- iPhone 4s and higher
- iPod Touch 5th Gen and higher
- iPad 3rd Gen and higher
- BLE (Bluetooth low Energy) compatible Android device

The Nokia Steel HR™ cannot be set up from a computer.

## Compatible OS

- iOS 8.0 and higher
- Android 6.0 and higher

# Document Release Overview

---

Release date		Modifications
October 2016	v1.0	First release
January 2017	v2.0	Update of: “Tracking my activity” on page 25
May 2017	v3.0	Nokia graphic chart:



**The screenshots in this manual are used for explanatory purposes. Your actual screens may differ from the screenshots in this manual.**

# Warranty

---

## Nokia® One (1) Year Limited Warranty - Nokia Steel HR™

Nokia® Inc. of 16192 Coastal Highway, Lewes, DE 19958 (“Nokia®”) warrants the Nokia® branded hardware product (“Nokia Steel HR™”) against defects in materials and workmanship when used normally in accordance with Nokia®’ published guidelines for a period of ONE (1) YEAR from the date of original retail purchase by the end-user purchaser (“Warranty Period”). Nokia®’ published guidelines include but are not limited to information contained in technical specifications, safety instructions or quick start guide. Withings does not warrant that the operation of the Nokia Steel HR™ will be uninterrupted or error-free. Nokia® is not responsible for damage arising from failure to follow instructions relating to the Nokia Steel HR™’s use.

### **The warranty does not extend to:**

- Normal wear and tear such as alteration of the color and/or the material of the wristband, deteriorated glass case, or peeling of the watch case plating.
- Any damage to the watch and its parts resulting from abnormal or abusive use, negligence or lack of care, or accidental damage (such as dents).
- Indirect damages of any kind.

# Regulatory Statements

---

FFCC ID: XNAHWA03

This device complies with Part 15 of the FCC Rules and with Industry Canada licence-exempt RSS standard(s).

Operation is subject to the following two conditions:

- (1) This device may not cause harmful interference, and
- (2) This device must accept any interference received, including interference that may cause undesired operation.

Changes or modifications made to this equipment not expressly approved by Nokia® may void the FCC authorization to operate this equipment.

Radiofrequency radiation exposure Information:

The radiated output power of the device is far below the FCC radio frequency exposure limits. Nevertheless, the device shall be used in such a manner that the potential for human contact during normal operation is minimized.

This device complies with the R&TTE Directive.

A copy of the EU Declaration of Conformity is available online at [www.withings.com/compliance](http://www.withings.com/compliance)