MANUAL DE INSTRUCCIONES

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
MODELS: Pirámide, Ópera, Obelisco and Esfinge

ENGLISH

This wheelchair is designed for occasional or frequent use, and can be used indoors and outdoors by a single occupant of up to 110 kg. The wheelchair can be manually moved in forward and reverse. All Mobiclinic wheelchairs are manufactured to the highest quality and are CE compliant.

Folding and unfolding the wheelchair

To fold, simply lift up under the center edge of the seat upholstery. To unfold, tilt the chair slightly to one side to rise the wheels on the opposite side off the floor. Then press down on one or both seat rails. Before folding ensure that the footrest are either locked into position and the plastic footrests flipped up or remove them completely.


Applying the wheel locks

Push forwards on the lock tips, (or pull back on the pull to lock type), until the locks snap into the locked position. Do not attempt to enter or exit the wheelchair without having the locks securely engaged. Always activate both brakes simultaneously.

Footrest length adjustment

The position of the footplate on either standard footrests or elevating legrests is adjustable to fit the user's leg length. The footplate should be adjusted to support the weight of the user's foot and lower leg in such a position that permits weight bearing by the thighs. A footrest adjustment that is too long will result in a line of pressure under the thigh, at the front edge of the seat upholstery. An adjustment that is too short will raise the user's knees and cause excessive weight to be born by the buttocks. Either situation increases the risk of pressure sores. The knees and hips of the user should be at approximately the same level. The lowest edge of the footplate must be at least 7 cm above the floor to provide safe clearance when negotiating ramps or inclines. This adjustment is made by loosening the adjustment bolt with a wrench, and desired position, then re-tightening the bolt securely.

Using the stepper tube and kerbs

Use the stepper tube to raise the front casters (when mounting a kerb for example). To use, push down on the stepper tube with a foot. Do not raise the front casters by pushing down on the push handles as this could result in damage to the wheelchair.

To mount a kerb

Approach the kerb head on. Then the attendant uses the stepper tube to raise the front casters, and lowers the front casters on the raised kerb. Finally the attendant should push the wheelchair forward, lifting it up slightly to mount the kerb if required.

To go down a kerb

Line up the front casters with the edge of the kerb. The attendant uses the stepper tube to raise the front casters and tip the occupant slightly back. Keeping the castors raised, slowly lower the wheelchair down the kerb.

Warning: When going down a kerb, the front casters must be raised to prevent the occupant from falling out. The above are only recommendations and may not be suitable for all events.


IMPORTANT POINTS TO REMEMBER

The user or caregiver should perform basic safety checks on the wheelchair at frequent intervals.


1. Check hand grips and the rubber tips on the tipping levers to ensure that they are tight and secure.
2. Check the locks for proper adjustment to confirm that they lock the large wheels securely when engaged.
3. Look over all nuts, bolts, and attaching hardware for proper tightness.
4. If the wheelchair is equipped with pneumatic tires, check for proper tire pressure.
5. Check for proper footrest length adjustment.
6. If any unsafe conditions are found please discontinue use of the wheelchair.

Remember to engage the wheelchair lock before transferring to or from the chair.

Do not lean forward in the wheelchair unless both feet are flat on the floor. It's particularly important for users with heavy leg casts using elevating legrests.

Ensure that brakes are locked when getting on/off the wheelchair.

Do not step on the foot plate when getting off the wheelchair, or else the chair will tip forward.

Do not place the wheelchair near any sources of fire.

CLEANING AND MAINTENANCE

To keep your product clean, you can wipe down the surfaces with a damp cloth, never wet! After that and as a part of your product cleaning, be sure to dry everything thoroughly.

Do not use solvents, bleaches, synthetic detergents, wax glazes, sprays or other abrasive products which may damage it.

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