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ENGLISH

Dear Customer, thank you for choosing a MEDLUXY product. This MEDLUXY device complies with the provisions of European Regulation on Medical Devices, whereby it is classified as a medical product of class I.

Please carefully consider this instruction and read all the provisions prior to use. Follow the given care instructions. This product is for professional use only. Specially formulated for use in Dentistry.

Intended use: Used to restore a tooth by placing and carving a dental material to the tooth structure. Restorative dental hand instruments include – Amalgam Carrier, Composite placement instruments, carvers, pluggers and Amalgam Well.

How to use Amalgam Carrier: Push the tip of the carrier into the amalgam compound in order to fill the hollow cylinder of the carrier tip. Press the carrier plunger to deposit the amalgam compound on to the area being treated.

Target patient group: All

Improper use/contraindications to the use of the product: Any other use or use out-side the above intended use shall be considered improper. The manufacturer is not responsible for damage resulting from such use. In this case, all risks are assumed solely by the user.

Warnings: Keep out of reach of children. Use under supervision of an adult. **Secondary treatment limitations:** Frequent secondary treatment has little effect on these instrument. The service life of the product is usually defined based on the degree of wear and use-related damage. Ultrasonic cleaning is not permitted! High-speed autoclaving an hot air sterilization will damage the instrument!

Additional information: If used and stored properly, these MEDLUXY products will serve you reliably for many years.

Instructions: **Preparation for cleaning:** After using, we recommend immediately rinsing them under running water or in a slightly alkaline solution to avoid possible drying of residues after use (e.g., blood). In order to properly carry out steam sterilization (see "Cleaning, Sterilization" section), it is always necessary to carry out a thorough pre-cleaning!

Manual cleaning: According to the recommendation of the National Institute for Public Health and Environment (RIVM), the mechanical treatment is preferable. Manual cleaning is not recommended.

Automatic cleaning: 1. Thoroughly wash the instruments under running water before mechanical treatment to avoid the remnants of cleaning/disinfecting agents entering the machine. 2. Place the instruments in a special container. 3. Place the instrument container in the cleaning and disinfection device so that the water jets directly fall on the instruments. 4. Follow the instructions for the detergent dosage. 5. Start the Vario TD program followed by thermal disinfection. Thermal disinfection is provided in accordance with the established value of A0 (the disinfection time, in seconds, at a certain temperature) and in accordance with national regulations. 6. When the program is completed, remove the instruments from the cleaning and disinfection device and dry them (National Institute for Public Health and Environment (RIVM) recommends to dry with

compressed air). Special attention should be paid to the places of the instrument container that are difficult to reach during drying. 7. Visual inspection for integrity and cleanliness with a suitable magnification object (based on experience, 8x magnification is sufficient for visual inspection). If there are residues of contamination on the instruments after mechanical treatment, repeat cleaning and disinfection until they are completely cleaned. 8. Attention! Only mechanical cleaning (without confirmed disinfection) requires subsequent thermal disinfection in a steam sterilizer placed in suitable trays or mesh containers. **Disinfection:** For example, Helipur can be used to disinfect the instruments. The appropriate concentrations can be found in the manufacturer's instructions. **Sterilization:** The instruments can be successfully sterilized with the following parameters: Fractional pre-vacuum (three times), sterilization temperature: 134°C, exposure time: at least 3 min. (full cycle), drying time: 10 min.

Functional control and check: Perform a visual inspection for damage and wear.

Transportation and storage: Keep the instrument away from dust, humidity and contamination during storage and transportation. Temperature: -10°C to +50°C Relative humidity: from 30% to 75%

Disposal: Medical Instruments (stainless steel and plastic) must be handed over to a recycling centre after prior preparation. Please observe the local regulations.

Warranty: The warranty total 12-months from the date of purchase if used properly and in accordance with the User's Manual. If any issues arise or repair is necessary, please contact your dealer.

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