
- DE Bedienungsanleitung**
- EN Operating Instructions**
- FR Mode d'emploi**
- NL Handleiding**
- IT Istruzioni per l'uso**
- ES Instrucciones de uso**
- PT Manual de utilização**

MANUAL DOWNLOAD:

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MICROSCOPE GUIDE:

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
MICROSCOPE FAQ:

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EXPERIMENTS:

www.bresser.de/downloads

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RISK to your child!

Aids with sharp edges and tips are sometimes used with this device. Please store the device and all of its accessories and aids out of the reach of children. There is a risk of **INJURY**.

Children should only use this device under supervision. Keep packaging materials (plastic bags, rubber bands, etc.) away from children. There is a risk of **SUFFOCATION**.

The chemicals and liquids provided should be kept out of reach of children. Do not drink the chemicals! Hands should be washed thoroughly under running water after use. In case of accidental contact with the eyes or mouth rinse with water. Seek medical treatment for ailments arising from contact with the chemical substances and take the chemicals with you to the doctor.

Fire/Burning RISK!

Never subject the device - especially the lenses - to direct sunlight. Light ray concentration can cause fires and/or burns.

TIPS on cleaning

Clean the lens (objective and eye-piece) only with the cloth supplied or some other soft lint-free cloth (e.g. micro-fibre). Do not use excessive pressure - this may scratch the lens.

Dampen the cleaning cloth with a spectacle cleaning fluid and use it on very dirty lenses. Keep instructions and packaging as they contain important information.

DISPOSAL

Dispose of the packaging material/s as legally required. Consult the local authority on the matter if necessary.

EC Declaration of Conformity

Bresser GmbH has issued a "Declaration of Conformity" in accordance with applicable guidelines and corresponding standards. The full text of the EU declaration of conformity is available at the following internet address:

www.bresser.de/download/8859460/CE/8859460_CE.pdf

UKCA Declaration of Conformity

Bresser UK Ltd. has issued a "Declaration of Conformity" in accordance with applicable guidelines and corresponding standards. The full text of the UKCA declaration of conformity is available at the following internet address:

www.bresser.de/download/8859480/UKCA/8859480_UKCA.pdf

Bresser UK Ltd. • Suite 3G, Eden House, Enterprise Way • Edenbridge, Kent TN8 6HF • GB

Your experiment set is made up of the following components:

- 1 20 prepared slides with QR Code
- 2 6 Blank slides
- 3 Slide carrying case
- 4 7 Vials of specimens
 - Shrimp eggs (Artemia)
 - Sea salt
 - Yeast
 - Dry shrimp
 - Textil (Stoff/Gewebe)
- 5 2 Empty vials
- 6 Scalpel
- 7 Tweezers
- 8 Dissecting needle
- 9 Magnifying glass
- 10 Petri dish
- 11 Pipette
- 12 Slide covers and adhesive labels
- 13 Prepared slide "Butterfly Wings"
- 14 12 stones and QR Code

What you need to know

Caution!

When using your experiment set, wear an apron (or old clothing) and rubber gloves! Always carry out your experiments under the guidance and with the help of an adult!

Your experiment set includes numerous accessories that help you carry out different experi-

ments. You can use all the parts independently of the microscope that you use for your research.

Experiments

Use the following web link to find interesting experiments you can try out. How do you hatch brine shrimp? How to make salt crystals and much more.

<http://www.bresser.de/downloads>

You can do them with the individual parts of your experiment set.

Prepared Slides

The prepared slides contain a QR code. Let your parents show you how to read a QR code on your smartphone to learn more about the specimen.

You will need a transmitted light microscope to view the prepared slides (e.g. 8855000, 8855002).

Brine shrimp eggs (Artemia)

These are special eggs that are capable of survival even when they are dried out. As a result, they are particularly ideal for raising brine shrimp.

Sea salt

You can observe the sea salt using your micro-

scope. Salt is a crystal that looks very interesting through a microscope. Apart from that, you will also need the sea salt to make a salt solution, in which you will hatch your brine shrimp.

Yeast

The yeast in your experiment set is included as feed for the brine shrimp. If you feed them regularly, you will be able to watch how they grow.

Prepared slides "Textile" (material/fabric), "Dry Shrimp" and "Butterfly Wings"

Different prepared slides that you can investigate using your microscope.

Empty vials

You can keep a few samples in the empty vials (e.g. leaves, and, etc.)

Slide covers

You use the slide covers to cover specimens that you have placed on a slide in order to observe with a microscope.

Adhesive labels

You can use these labels to identify your permanent prepared specimens. You can also use them to label the empty vials.

Slides

You place the specimen on a slide. After you have covered it with some water or gum media, you can

place a cover slip on top. Then you place the slide under the clips on the microscope.

Scalpel

The scalpel is a sharp knife that you can use to slice specimens/samples.

Tweezers

These are a type of pincers that you can use to pick up smaller objects.

Dissecting needle

A dissecting needle can be used for many things. You can use it to scratch the surface of a specimen, to fix objects in place or to stir liquids.

Pipette

You can use this instrument to soak up liquids in small amounts and squirt them out again. To soak up liquids, push the upper part (head) of the pipette together and put the opening on the lower part in the liquid. Let the head go and the pipette will fill with the liquid.

When you want to release some of the liquid, simply push the head together again. The more carefully you press, the less liquid will come out of the opening.

Magnifying glass

With the magnifying glass, you can observe objects with 2 times the magnification.

Stones

12 stones and QR Code. Let your parents show you how to read a QR code with your smartphone to learn more about the stones.

To view the stones you need a reflected light microscope.

How do I make my own specimens?

Take the object that you want to observe and place it on a glass slide (17). Then, add a few drops of distilled water on the object using a pipette. Now, place a cover slip vertically at the edge of the drop of water, so that the water runs along the edge of the cover slip. Then, slowly lower the cover slip over the water drops.

After use

To make your experiment set - and the enjoyment of it - last longer, you should make note of the following for after each use.

1. Clean the slide covers, slides and the microscope instruments carefully with water and some soap.

Caution!

Slide covers, slides and the instruments have sharp edges and are in some cases pointed. Make sure to always let adults help you when cleaning, so that you don't hurt yourself!

Warning! The slide covers are particularly thin and could break. It is best to place them on a piece of paper on a hard surface for cleaning.

2. Close all the vials tightly after use. Some containers have liquids that evaporate very easily, and could dry up.

Warranty & Service

The regular warranty period is 2 years and begins on the day of purchase. To benefit from an extended voluntary warranty period as stated on the gift box, registration on our website is required.

You can consult the full warranty terms as well as information on extending the warranty period and details of our services at www.bresser.de/warranty_terms.

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