C CARE

POA300

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C CARE

We thank you for buying our C-CARE Pulse Oximeter POA300 and thank you for the confidence you have placed in our brand and products.

We are convinced that with this purchase you have acquired a high quality product that you will enjoy using for many years.

Please read these operating instructions carefully before using the product, even if you are familiar with the type of product. It gives you useful tips for getting the most out of your product.

If, after reading this manual, you still have questions about this product and how to use it, please send an e-mail to the address below; we will answer your question immediately.

CRESTA INTERNATIONAL B.V.
ROLBRUGWEG 4
NL-1332 AS ALMERE
NEDERLAND
info@cresta-international.com
www.cresta-care.com
Tel: + 31 20 6932932

INTRODUCTION

General

This chapter contains a general description of the Oximeter POA300 including a brief description of the product and the product features

Applied use

The Oximeter is a portable, non-invasive device intended for random monitoring of functional oxygen saturation of your arterial hemoglobin (SpO₂) and your pulse rate (PR). This portable device is intended for use with adult patients in clinical settings and in the home environment. It is not suitable for continuous monitoring.

Short product description

The Oximeter is based on digital technology and is intended for non-invasive sample measurements of functional oxygen saturation of arterial hemoglobin (SpO₂). The advanced DSP algorithm * can minimize the influence of motion artifacts and improve the measurement accuracy of low perfusion *.

The Oximeter can be used to measure human SpO₂ and heart rate through the finger. The product is suitable for use in family, hospital -including clinical use in internist / surgery, anesthesia, pediatrics etc.- oxygen bars, in social medical organizations, physical care in sports etc.

Note: This product is a medical device for home use in the private environment, do not take any therapeutic measures based on a self-performed measurement. Never change the dosage of any medication prescribed by your doctor and always follow the instructions of the treating doctor.

We recommend that you submit any medical questions, complaints or symptoms to the attending physician or other healthcare professional in a timely manner in order to receive professional examinations, diagnostics and possibly medical care.

CAUTION: THIS IS A MEDICAL PRODUCT THAT MUST BE USED ACCURATELY ACCORDING TO THE INSTRUCTIONS. READ THESE INSTRUCTIONS FOR USE CAREFULLY FIRST.

SAFE USE

Instructions for the safe use and operation of the POA300 Oximeter.

- Do not try to open this Pulse Oximeter. Only qualified personnel can carry out any repairs if required.
- Periodic sensor site changes may be necessary for extended use or depending on the patients condition. Change the sensor location and check the skin health, circulatory status and correct alignment at least every two hours. If you have no results on one finger, then try another finger.
- SPO2O measurements can be adversely affected by the presence of strong ambient light. Protect its sensor area with, for example, a towel from direct sunlight if necessary.

The following causes can interfere with the POA300 oximeter's test accuracy:

- Presence of high frequency electrosurgical equipment.
- Placement of a sensor in a site with an arterial blood pressure cuff catheter, or intravascular line.
- The patient has high blood pressure, severe narrowed arteries, severe anemia or hypothermia.
- The patient is in cardiac arrest or in shock.
- Nail polish or fake nails can cause inaccurate SPO2 measurements.
- The device should be kept in the operating environment for at least 10 minutes to allow the product to adapt to the ambient temperature before measuring
- The device is non-sterile and not intended to be sterilized

WARNINGS

- Although this Medical Equipment complies with the requirements of the EN 60601-1-2 standard regarding electromagnetic compatibility, electrical equipment may cause interference. If malfunction is suspected, move the equipment away from the sensitive device or contact us.
- Portable and mobile RF communications equipment can affect the normal operation of this instrument.
- HAZARD OF EXPLOSION - Do not use the pulse oximeter in a flammable environment where concentrations of flammable anesthetics or other flammable materials are present.
- Do not remove batteries by throwing them into an open fire; this could cause the batteries to explode.
- Do not attempt to recharge normal dry cell batteries. These could leak or possibly cause a fire or the batteries may explode.
- Do not use the pulse oximeter in an MRI or CT scan environment.
- Do not make any changes to the product without the permission of the manufacturer.
- If changes are made, extensive research and testing is necessary to ensure that the equipment can be used safely.
- Stay away from active high-frequency surgical equipment and the shielded radiofrequency area of a medical X-ray system, where the intensity of electromagnetic disturbances is high.
- Use of this equipment adjacent to or stacked with other equipment should be avoided as it may cause improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify normal operation.
- Use of accessories, converters and cables other than those specified or supplied by the manufacturer of these equipment can lead to increased or decreased electromagnetic emissions, electromagnetic immunity of this equipment and result in improper operation.
- Portable RF communications equipment - including peripherals such as antenna cables and external antennas - should be used not closer than 30 cm /12 inches to any part of the POA300 pulse oximeter. This is valid also for cables specified by the manufacturer. Otherwise, performance degradation may occur.
- High pressure sterilization is not possible with this device

C CARE

If applicable: a list of all cables and maximum cable lengths, converters and others accessories, which are replaceable by the responsible organization and are also likely to affect the compliance of the Medical Equipment or Medical System with the requirements of Clause 7 - Emissions - and Clause 8 - Immunity, Accessories can be specified generically (e.g. shielded cable, load impedance) or specific (e.g. by manufacturer and equipment or type reference).

If applicable: the performance of the Medical Device or Medical System is defined as the "essential performance", or the description of what the user can expect when using it.

This "essential performance" can be lost or deteriorated as a result of Electromagnetic Interference.

Caution: Keep the operating environment free from dust, vibrations, corrosive or flammable materials and extreme temperatures and humidity.

Note: remove the batteries from the device if you are not going to use it for an extended period of time.

Note: Only use the device with the battery cover closed

Caution: Dispose of used batteries only according to local regulations. For your correct local regulations, please contact your municipality.

Note: keep this product away from children and pets; this to prevent swallowing of small parts.

C CARE

DEFINITIONS AND SYMBOLS

SYMBOL	DESCRIPTION	SYMBOL	DESCRIPTION
	Type BF equipment		Batch code
	Information about the manufacturer		Production date
	Temperature limitation		Serial no.
	Do not dispose of with normal household waste		Information about the EU authorised representative
	Read the manual	NOTE	Information to prevent damage
IP22	Anti dust & anti water classification		Important to know
Warning		Information to protect patients and medical staff for injury	
		The pulse oximeter does not sound an alarm for readings outside the normal limits.	

PRODUCT FEATURES

- Lightweight, easy to use
- Manually adjust the interface
- OLED color display, simultaneous display for value testing and plethysmography.
- Low perfusion: 0.3%. (The advanced DSP algorithm can improve the measurement accuracy, under the condition of low perfusion.)
- Image and sound reminder function
- Real-time sampling
- Low battery indicator.
- Automatic shut-off
- Standard operation on 2 AAA 1.5 V Alkaline batteries, good for 20 hours of continuous use

Note: This device is not suitable for measurements on children under 3 years old, the accuracy of the measurement cannot be guaranteed

Note: The pulse oximeter is intended only as an aid in patient assessment. It should be used in conjunction with other methods of assessing clinical signs and symptoms.

Note: The patient is the intended user and can perform maintenance on the equipment themselves

Note: A function tester cannot be used to assess the accuracy of a pulse oximeter or sensor. Clinical tests are used to determine the SpO₂ accuracy. The measured arterial SpO₂ value (SpO₂) from the sensor is compared with the value of arterial hemoglobin oxygen (SaO₂) determined from blood samples with a laboratory CO-oximeter. The accuracy of the sensors compared to the CO-oximeter samples measured over the SpO₂ range of 70 - 100%.

DSP algorithm: algorithm for digital signal processor.

- Low perfusion: In physiology, perfusion is the process by which a body delivers blood to a capillary bed in its biological tissue. Under the condition of low perfusion, the measurement of non-invasive pulse-blood oxygen saturation may show low accuracy.
- Plethysmograph: is an instrument for measuring changes in volume in an organ or the whole body (usually due to fluctuations in the amount of blood or air it contains).
- PI (Perfusion Index) is the ratio of the pulsating blood flow to the non-pulsating static blood flow in the peripheral tissue, such as fingertip, toe, or earlobe. The Perfusion Index is an indication of the pulse strength on the sensor.

INSTALLATION, CONFIGURATION AND OPERATION

Description of the front

1. On / Off button
2. Display of SPO2 / PR data and plethysmography
3. Battery compartment

Display

After you turn on the POA300, the display will look like:

Parameter setting

When the device is under the measurement interface, press the Direction button for 1 second to enter the Menu page (Figure 3.3.1 and Figure 3.3.2).

There are two submenus to choose from:

Menu 1: Reminder installation

Press the direction button for 1 second and go to Set Reminder. The user can adjust the setting by moving the "" symbol to the back of the sound reminder, beep, recovery or brightness.

Sound reminder

Press the direction button for 1 second, move the symbol "" to the back of the sound reminder and long press the direction button to turn it on / off.

Note: If the measured value exceeds the maximum or minimum value of SPO2 or PR, sound will follow when sound reminder is turned on.

- **Beep**
Press and hold the directional button for 1 second, move it "" symbol towards the back of Beep, long press the directional button to turn it on / off.

Note: When the beep is on, the sound emitted during the test will indicate the sound of the heartbeat

- **Recovery**
When the symbol "" is displayed after "Restore", a long press on the direction button can be changed to "OK", which will restore the factory settings.
- **Demo**
Press and hold the directional button for 1 second, move the "" symbol to the back of the demo and long press the directional button to turn it on / off.
- **Brightness**
When the "" symbol is displayed at "Brightness", long press the Direction button to change the brightness value from 1 to 5.

Menu 2: Limit setting

When the * symbol is displayed in the Reminder Settings, long press the Direction button until you enter the Reminder Limit setting menu (Figure 3.3.2). The user can press the direction button to select the items. Press the direction button for 1 second to change the data you need.

On the Settings menu page Reminder Limit (Figure 3.3.2), when the * symbol is displayed after the "+/-". Press and hold the direction button for 1 second to change the "+" to "-" or the "-" to "+".

When "+" is displayed on the right side, press the Direction button for 1 second, move the "" after the Spo2 Hi or PR Hi setting to increase the value to a higher value (until it reaches the highest value) .

When "-" is displayed on the right side, press the Direction button for 1 second, move the "" after the Spo2 Lo or PR Lo value setting, can decrease the value to a lower value until it has reached the lowest value.

Remind Setup	*
Sound Reminder	on
Beep	off
Demo	on
Restore	OK
Brightness	4
Exit	

FIG. 3.3.1.

Limit Setup	*
SpO2 Hi	100
SpO2 Lo	94
PR Hi	130
PR Lo	50
+/-	+
Exit	

FIG. 3.3.2

C CARE


Attention:

The sound reminder has a delay of 1 second after the incorrect result is detected. The user can set the limit value to 98 or 99 to check if it is normal to set a sound reminder.

SERVICE

Installing the batteries

Place 2 AAA batteries in the battery holders in the correct + and - polar direction and close the cover


WARNING: Do not attempt to recharge regular alkaline batteries, they may leak and catch fire or even cause an explosion.

Turn on the oximeter.

Pulse oximeters react strongly to movements, Keep your hands as still as possible during the measurement. Good blood flow is necessary for an accurate measurement.

If you have cold hands or if your circulation is not optimal, rub your hands gently before measuring to stimulate the blood flow.

Insert one of your fingers into the rubber opening of the pulse oximeter. Place the finger as far back as possible with the nail surface up, then release the clamp.

Read the corresponding data from the display screen

C CARE


WHAT IS THE MEANING OF THE DISPLAYED DATA?

The oxygen saturation SpO₂ of the blood indicates how much of the red blood pigment - hemoglobin - is loaded with oxygen. The normal value in humans is between 90 and 96% SpO₂. A too low value can indicate the presence of certain diseases such as a heart defect, circulatory problems, asthma or certain lung conditions.

A too high value can be caused, for example, by rapid and deep breathing, which, however, entails the risk of an insufficient carbon dioxide level in the blood. The result obtained with this device is in no way suitable for establishing or confirming diagnoses - always consult your doctor about this.

OLED display description

The display interface of OLED can rotate four directions with six different display modes after pressing the power button for less than 0.5 second. It is shown as below:

Attention:

1. When the battery is low, the indicator displays the symbol in the OLED display to remind users to replace the battery.
2. Plethysmography can be considered correct if the wave fluctuates regularly.

CLEANING AND DISINFECTION

Cleaning

Turn off the power and remove the batteries before cleaning. Keep the surface of the appliance clean, free of dust and dirt. Clean the surface - including the OLED display screen - of the device with a dry and soft cloth. Use 75% medical alcohol to clean the surface. Use a dry cloth with little alcohol to prevent alcohol from entering into the device.

Disinfection

Disinfect the device after each patient's use if multiple patients are using the device in the hospital.

Use 75% medical alcohol to clean the surface which has been in contact with the patient.

Note: Do not use strong solvents such as acetone

Note: Do not use abrasives or steel wool

Caution: Do not allow liquid to enter into the product and do not immerse the device in any liquid.

Note: Do not pour any liquid over the device while cleaning

Note: Do not leave any cleaning agent on the surface

C CARE

TROUBLESHOOTING AND MAINTENANCE

PROBLEM	CAUSE	SOLUTION
Display Oxyhemoglobin or heart beat incorrect.	Finger not correctly positioned Blood flow too low to measure	Reapply finger Try several times. If you are not sure whether it is the product, go to your doctor in time for a diagnosis
Display Oxyhemoglobin or heart beat unstable	Finger not inserted deep enough Your body may vibrate	Reapply the finger Try to move less. Relax yourself
Display Oxyhemoglobin or heart beat abnormal and causes a sound reminder.	Finger not inserted correctly Patients SpO2 or heart rate is abnormal	Reinsert the finger Visit a doctor for further examination
Pulse oximeter does not turn on.	No or batteries with insufficient voltage installed Batteries inserted incorrectly Product is damaged	Replace batteries Reinstall the batteries Please contact customer service
The display turns off unexpectedly	The product turns itself off after 16 seconds of inactivity Battery power too low	Normal use Possibly replace the batteries

C CARE

SPECIFICATIONS

PRODUCT	PULSE OXIMETER
MODEL	POA300
ANTI -ELECTRIC SHOCK TYPE	INTERNALLY POWERED EQUIPMENT
ANTI-ELECTRIC SHOCJ EQUIPMENT DEGREE	TYPE BF
EMC TYPE	TYPE B CLASS 1
ENCLOSURE DEGREE OF INGRESS PROTECTION	IP22
INTERNAL POWER	2 X AAA 1,5 V ALKALINE BATTERIES
POWER CONSUMPTION	BELOW 45 mA
DISPLAY	0,96" OLED
SPO2 DISPLAY	35 - 100%
HEART RATE DISPLAY	30 - 250 BPM
PI DISPLAY	0 - 20%
RESOLUTION	SPO2: 1%
	HEART BEAT: 1 BPM
	P.I.: 0,1%
MEASUREMENT ACCURACY	SPO2: $\pm 3\%$ (70 - 100%) UNDER 70% NOT SPECIFIED
	HEART BEAT: ± 2 BPM
	PI: 0,1% (0 - 1%) 1,0% (1 - 20%)
OPERATING ENVIRONMENT	TEMPERATURE: 5°C - 40°C (41°F - 104°F) HUMIDITY: 15 - 85% NON CONDENSING AIRPRESSURE: 86 - 106 Kpa
TRANSPORT AND STORAGE ENVIRONMENT	TEMPERATURE: -25°C - 55°C (-13°F - 131°F) HUMIDITY: 15 - 85% NON CONDENSING
MEASUREMENTS	63 X 36 X 34 MM
WEIGHT	55 \pm 2 GRAM (INCL. BATTERIES)
ACCESSORIES	2 X AAA BATTERIES HANG STRING USER MANUAL

MANUFACTURER'S DECLARATION OF THE EMC

1. All necessary instructions for maintaining basic safety and essential performance with respect to electromagnetic disturbances throughout its life.
2. Guidance and Manufacturer's Declaration - Electromagnetic Emissions and Immunity.

Table 1 - Guidance and manufacturer's declaration - electromagnetic emissions	
Emissions test	Compliance
RF emissions CISPR 11	Group 1
RF emissions CISPR 11	Class B
Harmonic emissions IEC 61000-3-2	Not applicable
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not applicable

Table 2- Guidance and manufacturer's declaration - electromagnetic immunity		
Immunity Test	IEC 60601-1-2 Test level	Compliance level
Electrostatic discharge (ESD) IEC 61000-4-2	38 kV contact 32 kV, 34 kV, 38 kV, 315 kV air	38 kV contact 32 kV, 34 kV, 38 kV, 315 kV air
Electrical fast transient/burst IEC 61000-4-4	Power supply lines: 32 kV Input/output lines 31 kV	Not applicable
Surge IEC 61000-4-5	Line(s) to line(s): 31 kV. Line(s) to earth: 32 kV. 100 kHz repetition frequency	Not applicable
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% 1 cycle and 70% 25/30 cycles Single phase: at 0 0% 300 cycles	Not applicable
Power frequency magnetic field IEC 61000-4-8	30 A/m 50Hz/60Hz	30 A/m 50Hz/60Hz
Conducted RF IEC61000-4-6	150KHz to 80MHz: 3Vrms 6Vrms (in ISM and amateur radio bands) 80% Am at 1kHz	Not applicable
Radiated RF IEC61000-4-3	10 V/m 80 MHz - 2,7 GHz 80 % AM at 1 kHz	10 V/m 80 MHz - 2,7 GHz 80 % AM at 1 kHz
Noe: U_T is the A.C. mains voltage prior to application of the test level.		

Table 3: Guidance and manufacturer's declaration – electromagnetic immunity

	Test Frequency (MHz)	Band (MHz)	Service	Modulation	Modulation(W)	Distance(m)	Immunity test Level (V/m)
Radiated RF IEC61000-4-3 (Test specifications for enclosure port immunity to RF wireless communication equipment)	385	380-390	TETRA 400	Pulse Modulation 18 Hz	1.8	0.3	27
	450	430-470	GMRS 460, FRS 460	FM 35kHz deviation 1 kHz sine	2	0.3	28
	710	704-787	LTE Band 13, 17	Pulse modulation 217 Hz	0,2	0.3	9
	745						
	780						
	810	800-960	GSM 800/900, TETRA 800, IDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18 Hz	2	0.3	28
	870						
	930						
	1720	1700-1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band1,3,4, 25; UMTS	Pulse modulation 217 Hz	2	0.3	28
	1845						
	1970						
	2450	2400-2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation 217 Hz	2	0.3	28
	5240	5100-5800	WLAN 802.11 a/n	Pulse modulation 217 Hz	0,2	0.3	9
	5240						
5785							

C CARE

DISPOSAL

Follow applicable regulations when disposing of the pulse oximeter and batteries. This pulse oximeter should not be disposed of with THE normal household waste.

All users are required to return all electrical - or electronic equipment, regardless of whether it contains toxic substances, to a municipal or commercial collection point for disposal in an environmentally sound manner.

Remove the batteries before disposing of the pulse oximeter. Do not dispose of old batteries with the household waste, but hand them in at a collection point for batteries at your local recycling point or at a store. For the exact local regulations and collecting points please contact your municipality.


The waste bin symbol below will remind you of this.

MANUFACTURER:

Shenzhen Aeon Technology Co., Ltd.
RM6H02, Block 27-29, Tianxia IC Industrial
Park, Majialong, No.133 of Yiyuan road, Nantou Street,
Nanshan District, Shenzhen, CHINA.

Shanghai International Holding Corp. GmbH (Europe)
Eiffestrasse 80, 20537 Hamburg, Germany

Printing errors and modifications reserved

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Europe - EU Declaration of Conformity

According to the Medical Device Directive. Classification: (MDD, Annex IX): II a, Rule 10.

UMDNS CODE: 17148 Pulse oximeter

GMDNS CODE: 37770

Manufactured for:

Name: CRESTA INTERNATIONAL B.V.
Address: Roibrugweg 4 - 1332 AS Almere
The Netherlands

Equipment:

EAN number: 8711974340130
Report reference: A300
Marketing name: C CARE POA300
Class: Class IIa, rule 10

Hereby, Cresta International BV herewith declares that the stated medical device product meets the transportation into national law, the provisions of Council Directive 93/42/EEC of 14.6.1993 concerning medical devices; including, at 21.3.2010, the amendments by Council Directive 2007/47/EEC.

All supporting documentations are retained under premises of the manufacturer and the notify body.

Directives:

Medical Device Directive: Council Directive 93/42/EEC Annex II, excluding Section 4

ISO13485: 2016

IEC60601-1: 2005 + Corr.1: 2006 + Corr.2: 2007 + A1: 2012

EN60601-1-2: 2015

IEC60601-1-11: 2015

ISO 80601-2-61: 2017

EN ISO14971:2012

ISO 10993-5: 2018

ISO 10993-10: 2010

EN 62304:2006/A1:2015

EN 62366-1: 2015

ISO 15223-1: 2016

EN1041: 2008

MEDDEV 2.7.1: 2016 REV 4

RoHS Directive (2011/65/EU)

REACH Directive (1907/2006/EC)

Standard:

All applied harmonized Standards (published in the Official Journal of the European Communities)

Certificate:

Registration: G1 073767 0013 REV. 00

MANUFACTURER:

Shenzhen Aeon Technology Co., LTD, RM6H02, Block 27-29, Tianxia IC Industrial Park, Majialong, No. 133 of Yiyuan Road, Nantou Street, Nanshan District, Shenzhen, China.

NOTIFIED BODY:

TÜV SÜD Product Service GmbH - Rikferstrasse 65 - D-80339 München - Germany - CE 123

EC REP:

Authorized representative:

Shanghai International Holding Corp. (Europe);

Eiffestrasse 80, D-20537 Hamburg, Germany

Signed for and on behalf of:	CRESTA INTERNATIONAL B.V.	 
Place:	Almere	
Date:	October 1st, 2020	
Name:	A. Benay	
Title:	Director	