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- Do not use test strips after the use-by date. If expired test strips are used, test results may be inaccurate. The use-by date can be found on the container beside the hourglass icon 🕒.
- The test strips can be stored for up to three months after opening the container. Write this use-by date (date of opening + 3 months 🕒) on the label. The storage time is reduced if it extends beyond the use-by date (see date beside the hourglass icon 🕒).
- Do not use the test strips if one of the use-by dates (🕒/🕒) has expired.
- The test strips can be touched anywhere with clean, dry hands.
- Use the test strip for measurement immediately on removal from the container.
- Do not bend, cut or change test strips in any way.
- Do not use test strips that have come into contact with liquids for measurements.

How they work

Test strips enable a quantitative measurement of the glucose level in whole blood. When the slot for blood taking comes into contact with a drop of blood, it is automatically filled by simple capillary action. The blood is drawn into by the absorption slot in the test strip and the monitor measures the glucose content of the blood.

The test is based on the measurement of an electric current that is generated by the chemical reaction of the glucose with the enzyme glucose dehydrogenase (Aspergillus oryzae) on the strip. During the reaction, a mediator transports electrons through the electrode surface and so generates a current.

The monitor analyses the current. The current flow is proportional to the glucose content in the blood sample.

The results are shown on the display of the blood glucose monitor. Only a very small volume of blood is required (0.6 microlitre) and the actual measurement takes about five seconds. The test strips can measure glucose values from 20 to 630 mg/dL (1.1 to 35.0 mmol/L).

Checking test strips with the control solution

When you open a new container of test strips, you will need to check that your monitoring system works correctly with them. Use Beurer control solution LEVEL 3 + LEVEL 4 (REF 464.16). If the test result is in the control range, the system is operating correctly and your procedure is correct. The control range for the function test with control solution is printed on the test strip container. See the Instructions for Use of your monitoring system and the package insert of the control solution for details of testing with the control solution.

Measuring blood glucose

Read the Instructions for Use for your monitor on taking the blood sample, measuring the blood glucose level, and evaluating the test results.

Chemical components of the test-strip sensor:

- FAD-glucose dehydrogenase 6 %
- Potassium ferricyanide 56 %
- Non-reactive components 38 %

Restrictions

• A hematocrit value (proportion of red blood cells) between 30 % and 55 % has no significant influence on the measurement results.

• A very high or very low haematocrit value (proportion of red blood cells) may cause inaccurate measurements. If the haematocrit value is very high (over 55 %), the blood glucose readings may be too low; if the haematocrit value is very low (under 30 %), it may be too high. If you do not know your hematocrit value, ask your doctor.

- Do not use test strips to monitor the blood glucose of newborns.
- Use only fresh capillary whole blood. Do not use serum or plasma.
- Use capillary blood without squeezing the puncture site. If a limb is crushed, the blood is diluted with tissue fluid and this may lead to an incorrect measurement result.
- Do not use the test strips at altitudes above 7010 m.

Disposing of test strips

⚠️ WARNUNG
<p>When disposing of test strips, strictly observe the generally applicable precautions for working with blood. Dispose of all blood samples and materials that you or your patients come into contact with correctly to prevent injury and infection of other persons.</p>

Comparison of measured values with laboratory values

Präzision

Three lots of the blood glucose test strips have been tested to assess the precision of the GL44 and GL50 evo blood glucose monitors. This includes a repeat assessment using venous blood and a laboratory precision assessment using the control material. The blood glucose content of the venous blood samples ranges from 42.7 to 418.0 mg/dL (2.4 to 23.2 mmol/L) and control material from three concentrations is used.

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