FOR A GOOD REASON

GRUNDIG

**USER MANUAL
BENUTZERHANDBUCH
MANUEL DE L'UTILISATEUR
GEBRUIKSAANWIJZING
MANUAL DE INSTRUCCIONES
MANUALE DELL'UTENTE**

CE
0123

Model:MDI231
Item no.: 871125224865

EN:3-in-1 Ear and Forehead infrared thermometer
DE: 3-in-1 Ohr- und Stirn-Infrarot-Thermometer
FR: Thermomètre infrarouge auriculaire/frontal 3-en-1
NL:3-in-1 infrarood oor- en voorhoofdthermometer
IT:Termometro a infrarossi da fronte e orecchie 3 in 1
ES: Termómetro infrarrojo frontal y auricular 3 en 1

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GRUNDIG

EN

3-in-1 Ear and Forehead infrared thermometer

Model:MDI231

Item no.: 871125224865

CE
0123

USER MANUAL

Catalog

1. Summary of infrared thermometer	4
2. Safety precautions	5
3. Intended use	9
4. Temperature measurement mode and range description	9
5. Feature	11
6. Overall description	13
7. Operation instruction	14
8. Care and cleaning	22
9. Maintenance	23
10. Calibration	23
11. Storage	24
12. Accessories	24
13. Trouble-shooting	25
14. Specifications	26
15. Standard list	28
16. Disposal	29
17. Normalized symbols	30
18. EMC declaration	31

1. Summary of infrared thermometer

Thank you for choosing our Infrared Thermometer. The MDI231 Infrared Thermometer is used to measure an objects temperature based on the relationship between temperature and measurable infrared radiation. Simply aim the unit's probe toward the surface to be measured to obtain a quick and accurate temperature. To ensure proper use, please be sure to read this user manual carefully, paying close attention to the safety precautions.

- In order to use this product correctly, please read the user manual before use.
- In order to properly use this product, please carefully read the full text of this manual before using, in particular the "safety precautions" section.
- Please keep the instructions on the side for easy checking.

Basic principles:

All objects above absolute zero temperature emit certain percentage of infrared radiation energy based on its temperature. The amount of the radiation energy and the distribution of the wavelength have very close relationship. When human body temperature in 36-37°C, it emits wavelength 9-13um of infrared radiation. Based on this principle, according

to the relationship between surface forehead temperature and eardrum of ear temperature, we are able to measure the human body actual temperature.

2. Safety precautions

Warning !

- Use of this thermometer is not intended as a substitute for consultation with your physician. It is dangerous for user to perform a self-evaluation and self-treatment based on the measuring result. Be sure to follow doctor's instruction.
- Keep the thermometer out of reach of children. For accidental swallow of battery or other component, please consult the doctor at once.
- The patient who suffer the external otitis and media otitis or other ear disease that will be prohibited to using the product to measure ear temperature. It may worsen the lesion.
- Don't throw the battery into fire.

Notice !

- The device is precision instrument, don't drop, tramp or impose any vibration or impact on the thermometer.
- Do not touch the lens of the probe with your fingers, and do not disassemble the device by yourself.
- Before measuring forehead temperature, make sure the hair removed,

sweat dried.

- After you do some exercise, eating and bathing, you should stay still indoor about 30 minutes before measurement.
- To make the measurement data reliable and stable, when ambient temperature varies a lot, the thermometer should be placed indoors for about 30 minutes before using.
- When we measure somebody continuously, the temperature should be measured every minute, if you need to measure yourself continuously for a short time, there are some slight errors when you read the temperature, which is a normal phenomenon. At this time, we should choose the average. We recommend that you measure yourself continuously maximum of three in a unit of time, and because the temperature of the human will conduct to the thermometer, it may affect the accuracy of measurement.
- There is no absolute standard about the temperature of the human, so please try to collect the recording of individual temperature in the usual, as a reference for having a fever or not.
- Do not measure the sites of scarred tissue or tissue compromised by skin disorders, because sensing body temperature from sites of scarred tissue or tissue compromised by skin disorders.
- Do not measure the site of forehead temperature if that patients has

trauma on forehead.

- Do not measure if that patient is treated with certain drug therapies.
- Do not immerse the device into water or any other liquid, and not directly sunlight exposure.
- Do not use a mobile or cordless phone near the thermometer when measuring.
- Body temperature may increase in the drug within the effort time limit. Please don't measure.
- In order to ensure the accuracy of measurement data, please don't take measurement of body temperature in strong electromagnetic interference environment (such as microwave, high frequency equipment operation environment).
- Do not forcibly insert the probe in the ear. If you fell abnormal such as pain during measurement, stop using the unit. It may injure the external auditory canal.
- If the ear is cooled such as being exposed to the cold for a long period, wait until the ear is warmed up, then measure. The measured result may be indicated low when you use an ice bag or immediately after coming home from outside in winter.
- Do not use the unit to measurement ear temperature when the external auditory canal is wet such as after swimming or a bath. It may cause measurement inaccurate.
- Do not disassemble, repair, or modify the unit.
- This thermometer is only a perso-

nal device, please do not use with others.

- Do not touch the battery output when measuring.
- The thermometer must be stored according to the technical specifications.
- The materials (ABS) of expect contact with patient had passed the ISO 10993-5 and ISO 10993-10 standard test, no toxicity, allergy and irritation reaction. They are compliance with the MDD requirements. Based on the current science and technology, other potential allergic reactions are unknown.
- The patient is an intended operator. The patient can measure, read data and replace battery under normal circumstances and maintain the device and its accessories according to the user manual.
- The device is not intended for PATIENT transport outside a healthcare facility.

Recommendations

1. Don't use this thermometer for other purposes.
2. It is forbidden to leave the product exposed to any chemical solvent, direct sunshine or high temperature.
3. Don't expose the thermometer under direct sunlight long time so as not to damage the battery.
4. Do not measure while talking on the phone.

5. Please report to MANUFACTURER if any unexpected operation or events occur.

3. Intended use

This thermometer is intended to measure ear canal and forehead temperature at home or hospital, including anyone, such as infants, children and adults.

For the safety reason, children or the baby's temperature must be measured by parent or adults.

4. Temperature measurement mode and range description

The infrared thermometer has the following measurement mode:

1. Forehead temperature measurement mode -- measures the skin surface of human forehead's temperature accurately, take the place of traditional mercury thermometer and electrical thermometer.
2. Ear temperature measurement mode – measures the skin surface of human eardrum's temperature accurately, take the place of traditional mercury thermometer and electrical thermometer.

Normal temperature range for different measuring position

Measuring position	Normal temperature(°C)	Normal temperature(°F)
Anus	36.6-38.0	97.9-100.4
Oral	35.5-37.5	95.9-99.5
Armpit	34.7-37.3	94.5-99.1
Ear	35.5-37.8	95.9-100
Forehead	35.8-38.0	96.4-100.4

Normal forehead temperature range for different ages

Ages	Normal temperature (°C)	Normal temperature (°F)
0-2 years old	36.4-38.0	97.5-100.4
3-10 years old	36.1-37.8	97.0-100.0
11-65 years old	35.9-37.6	96.6-99.7
>65 years old	35.8-37.5	96.4-99.5

NOTE:

The normal temperature and difference between the different body parts is individual. To define yours, measure your temperature for at least 2 weeks at the same ear canal, forehead position and time.

NOTE:

When consulting your physician, communicate that which position the temperature is measured by infrared thermometer MDI231, note the individual's normal temperature range by infrared thermometer MDI231 as additional reference.

NOTE:

Because the forehead temperature is affected obviously by the external environment (eg: environment, air convection and skin tone, etc), we advice that you take the forehead temperature only as reference. When you have a doubt about the measurement result, please use the ear canal temperature to confirm it.

5. Feature

High reliability

This product has passed the manufacturer's internal life and reliability test, time to failure is $\geq 1000\text{h}$.

A wide range of temperature

The measurement range:
 $32.0^{\circ}\text{C}-43.0^{\circ}\text{C}$ ($89.6^{\circ}\text{F} - 109.4^{\circ}\text{F}$).

High accuracy

This product has passed the infrared thermometer performance standards of the European Union and China for measuring clinical requirements, measuring clinical repeatability is no more than $\pm 0.3^{\circ}\text{C}$ (0.5°F).

Humanization design

When the temperature exceeds the range, LCD will display the Lo or Hi prompt.

When operating environment exceeds the designed specifications, LCD will display the Err prompt.

When the thermometer battery po-

wer is insufficient, it has low voltage icon.

Has the hardware self-test function, when hardware malfunction is detected, LCD will display the ErA or ErC prompt.

Power saving function

Start the thermometer without any operation, or no any operation after temperature measurement, the thermometer will shut off and LCD go out with one short beep in $60s \pm 5s$.

Three-color backlight indication function

Has designed the three-color backlight function to indicate the temperature rang:When the correct result is measured, the backlight from LCD will light different color to indicate the temperature range as below:

Temperature Range °C	Temperature Range °F	Indicator Color	Duration of light	Way of Beep
$T < 37.5^{\circ}\text{C}$	$T < 99.5^{\circ}\text{F}$	Green	4 seconds	1 long beep
$37.5^{\circ}\text{C} \leq T < 38.0^{\circ}\text{C}$	$99.5^{\circ}\text{F} \leq T < 100.4^{\circ}\text{F}$	Orange	4 seconds	3 short beep
$38.0^{\circ}\text{C} \leq T$	$100.4^{\circ}\text{F} \leq T$	Red	4 seconds	3 short beep

Memory mode


The thermometer stores the last temperature measurements. The last temperature will be automatically displayed when it is turned on again. The display will show “M” icon.

Measuring indicator light


Measuring indicator light uses a green light to remind during the process of temperature measurement.

6. Overall description

Main component including

LCD display description

7. Operation instruction

Preparation

1. Check battery

Replace the batteries to ensure power supply if there is low voltage icon for the thermometer.

2. Check sensor

If have pollution and spray, please clean it. (The cleaning method see the chapter 8 Care and Cleaning for detailed.) If the sensor's lens is damaged, please stop using.

3. Check thermometer

When you press the Switch measurement button, the system will have self testing of software and hardware. If there are problems, LCD will display "ErC/ErA " symbol. If found, please refer to the chapter 13 Trouble-shooting.

4. In order to make the accurate measuring result, put the thermometer in the measurement environment for 30 minutes.

5. Accuracy of unexpected fluctuations in ambient temperature may decrease the measurement results. When the thermometer at the same measurement position display at different ambient temperature, or measure temperature in front of the air conditioner, it will not be able to obtain accurate results.

6. If you measure forehead temperature, clean forehead and arrange hair, make sure the forehead is naked and clean, in order to ensure the


accuracy of measurement.

Instruction for use

1. Battery installment and replacement

Insert new batteries when the battery symbol appears on the display. Open the battery compartment. Remove the batteries and replace with new batteries, making sure the poles are in the right direction. Slide battery door back until it snaps in place.

NOTE: always change batteries when LCD displays the low battery indicator to make sure a accurate measurement.


Notice:

- Please observe the related national laws of disposing the abandoned battery and don't litter to the garbage can.
- Please take out the battery if the device is not used for long periods of time.
- Please don't put the battery in the fire.

To protect the environment, dispose of empty battery at your retail store or at appropriate collection sites according to national or local regulations.


2. Function settings:

Switch temperature unit between °C and °F:

Remove the battery compartment cover. In temperature measuring state, press the C/F selector, and the selected temperature unit will be displayed on the screen with a “beep” sound. Then you can proceed to temperature measurements.

3. Switch between ear temperature mode and forehead temperature mode:

Ear temperature mode

LCD screen

Forehead temperature mode

Switch between ear temperature mode and forehead temperature mode by removing or installing the forehead temperature measuring cap.

4. Mode selection:

In temperature measuring state, press “1” to select adult mode, child mode or object temperature mode.


5. Enable/disable beep sound:

In temperature measuring state, press “2” to disable or enable the beep sound.

6. Call measurement records:

In temperature measuring state, long press “2” for about 3s to make the beep sound and display the measurement records in the above sequence. Then each time you press “2”, the previous record will be displayed (there are a total of up to 32 temperature value records).


7. Memory clearance:

In temperature measuring state, press and hold “2” and meanwhile press and hold the On/Off key for about 3s, and the screen will be cleared with four “beep” sounds (all records are deleted).

8. Measure ear temperature:

Press the On/Off key to turn on the thermometer and proceed to 2s full-screen display. Then the temperature value of the previous measurement will be displayed for 2s, and then the thermometer will get into the ear temperature mode with a “beep” sound. Insert the temperature measuring probe into the ear canal and press the On/Off key to perform a measurement, and the temperature value will be displayed on the screen with a “beep” sound. Such temperature value is displayed with green green($<37.5^{\circ}\text{C}$), orange ($37.5^{\circ}\text{C} \leq \text{temperature} < 38.0^{\circ}\text{C}$) or red ($\leq 38.0^{\circ}\text{C}$) backlight, as shown below.

NOTE:

Neonate (aged under 1): pull the auricle back and downward and hold it, so that the temperature measuring probe can aim at the eardrum.

Child or adult (aged above 1): pull the auricle back and upward and hold it, so that the ear canal will remain in a straight line and the temperature measuring probe can aim at the eardrum.

NOTE:

1. Keep the probe surface clean, otherwise the measurement can't be accurate.

2. Clean the probe with clean cloth or paper moistened with water or alcohol, and conduct measurement after the water on the surface of the probe all evaporates.

TEMPERATURE TAKING HINTS

- The right ear temperature may differ from the left ear. Therefore, always take the temperature in the same ear.
- The ear must be free from obstruction excess earwax build-up to take an accurate reading.
- External factors may influence ear temperatures, including when an individual has:
 - Been laying on ear or other
 - Had their ears covered
 - Been exposed to very hot or very cold temperatures, or
 - Been recently swimming or bathing
- In these cases, remove the individual from the situation and wait 30 minutes prior to tacking a temperature.
 - Use the untreated ear if prescription ear drops or other ear medication have been placed in the ear canal.
 - Don't use the thermometer outdoors.

9. Measure forehead temperature:
Install the temperature measuring cap to switch to the forehead temperature mode.

Press the On/Off key to turn on the thermometer and proceed to 2s full-screen display. Then the temperature value of the previous measurement will be displayed for 2s, and then the thermometer will get into the forehead temperature mode with a “beep” sound. Aim the temperature measuring probe at the forehead and press the Measure key (if you press and hold the Measure key, the thermometer will keep scanning with several “beep” sounds and then complete the measurement when you release the Measure key) to perform a measurement, and the temperature value will be displayed on the screen with a “beep” sound. Such temperature value is displayed with green ($<37.5^{\circ}\text{C}$), orange ($37.5^{\circ}\text{C} \leq \text{temperature} <38.0^{\circ}\text{C}$) or red ($\leq 38.0^{\circ}\text{C}$) backlight, as shown below.

NOTE:

Make sure the temperature measuring probe aims at and touches the forehead (between the eyebrows) closely.

TEMPERATURE TAKING HINTS

- Before take the temperature, please fix the forehead hair to prevent measured deviation.

- The sweat on head or cosmetic can affect the accuracy of measurement, please maintain the cleanliness of the forehead when measuring.
- It is normal that there may be temperature difference depending on various skin types and color, since different skin type will reflect different voltage of infrared ray.
- Don't use the thermometer outdoors.

NOTE:

1. When measuring, if the environment exceed the 10- 40°C range will be showing Err;
2. If the measuring process is wrongly happened, the backlight can be kind of help to remind that an inaccurate measurement has been taken.

Temperature Range °C	Temperature Range °F	Backlight
T<32.0°C	T<89.6°F	Green
T>43.0°C	T>109.4°F	Red
Duration of light	Way of Beep	LCD symbol
4 seconds	1 long beep	"Lo"
4 seconds	3 short beep	"Hi"

10. Three-color backlight function

The backlight from LCD will light different color to indicate the temperature range as below table:

Temperature Range °C	Temperature Range °F	Indicator Color	Duration of light	Way of Beep
$T < 37.5^{\circ}\text{C}$	$T < 99.5^{\circ}\text{F}$	Green	4 seconds	1 long beep
$37.5^{\circ}\text{C} \leq T < 38.0^{\circ}\text{C}$	$99.5^{\circ}\text{F} \leq T < 100.4^{\circ}\text{F}$	Orange	4 seconds	3 short beep
$38.0^{\circ}\text{C} \leq T$	$100.4^{\circ}\text{F} \leq T$	Red	4 seconds	3 short beep

8. Care and cleaning

The probe tip and lens are the most delicate part of the thermometer. It has to be clean and intact to ensure accurate readings.

If the thermometer is ever accidentally used, please clean the probe and lens as follows:

- Very gently wipe the surface with a cotton swab or soft cloth moistened with alcohol. Only use the thermometer after the alcohol has completely dried out.
- If the lens is damaged, contact the distributor.

Clean the unit body:

- Use a soft, dry cloth to clean the thermometer display and unit body.
- If very dirty, use a soft cloth with alcohol for cleaning.

NOTES:

- Don't use abrasive cleaners.

- Don't use other non-recommended methods to perform disinfection.
- Non-waterproof, don't use the abrasive cleaner to clean the product, don't drop the thermometer in the water or the other liquid.

9. Maintenance

1. We do not authorize any institution or individual to maintain and repair of the product. If you suspect that the products have any questions, please contact the manufacturer or distributor to handle the case.
2. The user must not attempt any repairs to the device or any of its accessories. Please contact the retailer for repair.
3. Opening of the equipment by unauthorized agencies is not allowed and will terminate any claim to warranty.

WARNING: No modification of this equipment is allowed!

10. Calibration

The thermometer is initially calibrated at the time of manufacture. If this thermometer is used according to the use instruction, periodic re-adjustment is not required. If any time your question the accuracy of measurement, please contact distributor or manufacturer, the contact information see last page.

11. Storage

1. Don't put the thermometer under the sunshine, high temperature and moist environment or some-place which maybe get in touch with fire or is vulnerable to vibration.
2. Take out the battery if don't use the device in a long time.

12. Accessories

Only use original accessories. Check that the contents of the delivery are complete.

Quantity	Parts
1pc	MDI231 device
1pc	User manual

13. Trouble-shooting

Troubles or error message	Checklists or situation	Countermeasures or solution
No response/ Automatically reset	The batteries are used up?	Replace new batteries.
	Battery in wrong polarity or type?	Take out the batteries and replace new ones.
	Poor battery contact	Take out batteries and reinsert it correct.
The thermometer show the symbol "Hi" 	Temperature hampered by an air flux.	Please leave the status and wait for 30 minutes to measure. Re-measure according to the manual.
	In the forehead measurement mode: --Temperature readings too close together. -- Measured the other object, such as the sunlight, the air from the fireplace. Hi: Higher than 43.0°C(109.4°F);	
The thermometer show the symbol "Lo" 	The hair and sweat prevent the temperature achievement.	
	Temperature hampered by an air flux.	
	In the forehead measurement mode: -- The measuring distance is too far. -- Measured the other object, such as the air from the air conditioner. Lo: Less than 32.0°C (89.6°F)	

	<p>The ambient temperature is beyond of range of measurement (10°C-40°C/50°F -104.0°F)</p>	<p>Keep the thermometer in the room where temperature is (10°C-40°C/50°F -104.0°F) for 30 minutes</p>
	<p>The hardware is damaged.</p>	<p>Excluding the possibility of temperature allowance first, then send the device to your dealer for repair</p>
	<p>Lower battery, however you can't use it</p>	<p>Replace the new battery.</p>

14. Specifications

Device name	Infrared Thermometer
Model	MDI231
Measurement mode	Forehead and ear measurement modes (Adjusted)
Measurement site	Forehead and ear canal
Power supply	d.c.3V ,1.5 V type AAA (LR 03) batteries
Measuring range:	32.0-43.0°C (89.6°F~109.4°F)

Measuring accuracy: (At laboratory conditions)	$\pm 0.2^{\circ}\text{C}/0.4^{\circ}\text{F}$ during $35.0^{\circ}\text{C}-42.0^{\circ}\text{C}$ ($95.0^{\circ}\text{F}-107.6^{\circ}\text{F}$); other $\pm 0.3^{\circ}\text{C}/0.5^{\circ}\text{F}$
Clinical repeatability:	within $\pm 0.3^{\circ}\text{C}$
Resolution of display	$0.1^{\circ}\text{C}/0.1^{\circ}\text{F}$
Operation condition	$10-40^{\circ}\text{C}(50^{\circ}\text{F}-104^{\circ}\text{F})$, Relative humidity $\leq \text{Rh}85\%$, 70-106kPa
Storage condition	$-20-55^{\circ}\text{C}(-4^{\circ}\text{F}-131.0^{\circ}\text{F})$, Relative humidity $\leq \text{Rh}85\%$, 70-106kPa
Size	140*36*50mm
Weight	about 60g(without batteries)
High body temperature hint	$\geq 38.0^{\circ}\text{C}(100.4^{\circ}\text{F})$
Electric shock	Internally powered ME equipment
Applied part	Type BF applied part, including the whole unit
Mode of operation	Continuous operation
Battery life	≥ 1000 times
Product life	2 years
Software version	V1.1
<p>Note: Not intended to be sterilized. Not for use in an OXYGEN RICH ENVIRONMENT</p>	

*The above specifications are subject to change without prior notice.

15. Standard list

ShenZhen Medek Bio-Medical Co.,Ltd. declares that the MDI231 complies with following applicable standards:

EN 980	Symbols for use in the labeling of medical devices
EN 1041	Information supplied by the manufacturer with medical devices
EN 60601-1	Medical electrical equipment Part 1: General requirements for basic safety and essential performance
EN 60601-1-2	Medical electrical equipment – Part 1-2: General requirements for basic safety and essential performance - Collateral standard: Electromagnetic compatibility - Requirements and tests
EN 60601-1-6	Medical electrical equipment – Part1- 6: General requirements for basic safety and essential performance – Collateral standard: Usability
ISO 80601-2-56	Medical electrical equipment part 2-56: particular requirements for basic safety and essential performance of clinical thermometers for body temperature measurement
EN 60601-1-11	Medical electrical equipment – Part 1-11: General requirements for basic safety and essential performance – Collateral standard: Requirements for medical electrical equipment and medical electrical systems used in home healthcare environment
EN 12470-5	Clinical thermometers – Part 5: Performance of infra-red ear thermometers (with maximum device)
ASTM E1965	Standard Specification for Infrared Thermometer for Intermittent Determination of Patient Temperature
EN 62304	Medical device software - Software life-cycle processes
EN 62366	Medical devices – Application of usability engineering to medical devices

16. Disposal

At the end of the product lifecycle, do not throw this product into the normal household garbage, but bring it to a collection point for the recycling of electronic equipment.

Waste Electrical and Electronic Equipment can have potentially harmful effects on the environment. Incorrect disposal can cause harmful toxins to build up in the air, water and soil and can be harmful to human health.

NOTES:

Handing of battery and wastes method, please act according to the native law to proceed to handle.

Take out the battery if you are not going to use the unit for a long time.

To protect the environment, dispose of empty battery at your retail store or at appropriate collection sites according to national or local regulations.

17. Normalized symbols

	<p>Manufacturer information: The manufacturer is ShenZhen Medek Bio-Medical Co.,Ltd.</p>
	<p>Authorized representative in the European Community. The European representative is: Shanghai International Holding Corp. GmbH (Europe).</p>
	<p>Complies with the European Medical Device Directive 93/42/EEC and amended Directive .2007/47/EC. Notified Body is TUV-SUD.</p>
	<p>Follow operating instructions</p>
	<p>Caution! Consult accompanying documents.</p>
	<p>Type BF applied parts</p>
	<p>Batch code</p>
<p>IP22</p>	<p>IP code of the device: this device's grade of against ingress of solid foreign objects</p>
	<p>Date of manufacture</p>
	<p>Disposal in accordance with Directive 2002/96/EC (WEEE)</p>
	<p>Keep dry</p>

	UP
	Fragile, handle with care
	Stacking layer limit

18. EMC declaration

1. The digital thermometer needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in the ACCOMPANYING DOCUMENTS.

2. Wireless communications equipment such as wireless home network devices, mobile phones, cordless telephones and their base stations, walkie-talkies can effect this equipment and should be kept at least a distance $d=1.0$ m away from the equipment.

Note: As indicated in Table 6 of IEC 60601-1-2 for ME EQUIPMENT, a typical cell phone with a maximum output power of 2W yields $d=1.0$ m at an IMMUNITY LEVEL of 10V/m.

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ShenZhen Medek Bio-Medical Co., Ltd. declara que el MDI231 cumple con los siguientes estándares aplicables: