INVENTUM

maakt 't moment

INFRARED THERMOMETER

infrarood thermometer

TMO430

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English

	normalized symbols	page 27
	introduction	page 28
1	safety instructions	page 29
2	appliance description	page 31
3	body temperature	page 33
4	before first use	page 33
5	how to measure forehead temperature	page 34
6	how to measure ear temperature	page 35
7	temperature taking tips	page 35
8	fever indication	page 36
9	check the memory	page 36
10	date and time setting	page 37
11	voice setting	page 37
12	replacing the battery	page 38
13	specifications	page 39
14	cleaning & maintenance	page 40
15	troubleshooting	page 41
16	important information regarding electromagnetic compatibility (EMC)	page 42
	<i>General terms and conditions of service and warranty</i>	page 52

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normalized symbols

Applied part of type BF: the applied part is the probe tip

Disposal in accordance with directive 2012/19/EU

Refer to Instruction Manual

The name and the address of the manufacturer

The name and the address of the authorised representative in the European Community

Complies with the European Medical Device Directive (93/42/EEC) and amended by directive 2007/47/EC requirements. Notified body TÜV Rheinland (CE0197)

IP22

The first number 2: protected against solid foreign objects of 12,5 mm \varnothing and greater. The second number: protected against vertically water drops when enclosure tilted up to 15 degrees. Vertically falling drops shall no harmful effects when the enclosure is tilted at any angle up to 15 degrees, on either side of the vertical.

Serial number

Exhausted batteries do not belong in household waste. Dispose of the batteries according to the current federal, state and local regulation. As a consumer, you are obligated by law to return spent battery.

NOTES:

1. Under the environment with electrostatic discharge, the unit may malfunction and may require user to reset the unit.
2. Dispose with care; please consult the retailer for details.

Date of manufacture

introduction

The thermometer has been carefully developed for accurate, safe and fast human body temperature measurements in the ear and on the forehead.

The reference body site is core and the mode of operation is adjusted mode. The device is intended for the intermittent monitoring of human body temperature in the home.


It is intended for use on people of all ages except pre-term babies or light-for-date babies. The patient may be an intended operator.

The patient can measure and change batteries under normal circumstances and maintain the device and its accessories according to the user manual.

The thermometer is an ear canal IR thermometer. The forehead temperature is obtained based on the site offset between the ear and the forehead, and the forehead temperature is only used as a reference.

This thermometer is intended for home use by the customer only. The measure site is the ear canal.

Research indicates that the ear is an ideal site for taking body temperature. Because the body temperature is regulated by the hypothalamus, which shares the same blood supply as the tympanic membrane. So please take the ear temperature as much as possible.


1

safety instructions

However, as with any thermometer, proper technique is critical to obtaining accurate temperatures. Therefore, in order to obtain an accurate measurement, please read this manual carefully before using.


- Use of this thermometer is not intended as a substitution for consultation with your physician. Measurement results are for reference only. Contact your physician if you have or suspect any health concerns.
- The infrared sensor must be kept clean, dry, and undamaged at all times to ensure accurate measurements.
- Earwax in the ear canal or lipid on the forehead may cause an inaccurate measurement, so make sure the ear canal or the forehead is clean.
- If you clean the ear or forehead, please wait for 5-10 minutes before measuring.
- Don't touch or blow on the infrared sensor. A polluted infrared sensor may cause inaccuracy.
- Do not use ear/ forehead measurement mode for measurements other than human body temperature.
- Avoid measuring the ear temperature if the ear is inflamed or infected.
- The temperature of the storage area differs a lot from that of the measuring area, please wait for the thermometer temperature to equalize to the room temperature about 30 minutes before use.
- Keep the device out of the reach of children and pets to avoid inhalation or swallowing of small parts.
- Do not allow children to take their temperatures unattended. Children may not be able to use the device according to the instructions in this user manual. It is not a toy.
- Don't throw battery into fire. Battery is likely to explode in fire.
- Don't store the thermometer to temperature extremes (below $-20^{\circ}\text{C}/-4^{\circ}\text{F}$ or over $55^{\circ}\text{C}/131^{\circ}\text{F}$) or humidity extremes (below 15%RH or over 90%RH). Failing to do so may cause inaccuracy.
- Don't use the product in a wet ear hole after swimming or bathing. It is possible to damage the ear hole.


- Don't use the thermometer if the main body is damaged (for example, the infrared sensor is broken). The continuous use of a damaged unit may cause injury, improper results, or serious danger.
- Don't take apart, repair or change any parts of the unit at any time.
- Don't use a mobile phone nearby when the thermometer is working.
- If your thermometer is not going to be used for a long time, it is advised to, remove the battery to prevent possible damage due to chemical leakage. If the battery does leak, remove carefully. Do not allow bare skin to touch leaking fluid.
- This product needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided, and this unit can be affected by portable and mobile RF communications equipment.
- The materials (e.g. ABS) of expect contact with patient had passed the ISO 10993-5 and ISO 10993-10 standards test, no toxicity, allergy and irritation reaction. However, based on the current science and technology, other potential allergic reactions are unknown.
- If you have an allergic reaction to materials, please stop treatment immediately and consult your physician.
- Store the thermometer in a dry, clean place. It is forbidden to leave the product exposed to any chemical solvent, lint, dust, direct sunshine or high temperature.
- When using, you shall not touch the batteries and the patient simultaneously.
- Do not maintain or service the device while the device is in use. If you have any problems with this device, such as setting up, maintaining or using, please contact your customer department of the manufacturer.
- Please report to manufacturer if any unexpected operation or events occur.
- Please avoid mechanical shock. Otherwise the thermometer will be damaged.

Caution: this device should not be used adjacent to or stacked with other equipment and that if adjacent or stacked use is necessary, this device should be observed to verify normal operation in the configuration in which is will be used.

2 appliance description

1. Infrared sensor
2. Memory button [MEM]
3. LCD display
4. Forehead cap
5. Ear measuring button []
6. Forehead measuring button [forehead]
7. Battery compartment
8. TALK button [TALK]
9. Table stand holder

1. Day sign in time mode
2. A.M. and P.M. sign in time mode (“AM”, “PM” will display on LCD when the hour mode is 12hr).
3. Memory position number
4. Indicator of the series of records in memory mode
5. Year in time mode or the series of records in memory mode
6. Year sign
7. Month and day in time mode or memory mode
8. Month sign
9. Night sign in time mode
10. Sleep sign(will display on LCD when sleep mode is “ON”)
11. Low battery
12. Temperature in measuring mode or time in time mode
13. High temperature sign
14. State of measurement sign; it displays on LCD during measurement
15. Sign of temperature measured in ear
16. Sign of temperature measured on forehead
17. Measuring temperature sign
18. SET sign means the setting mode
19. Speaker sign
20. Day sign

3 body temperature

The normal body temperature is a range. A person's normal temperature range tends to decrease with increasing age. The following table shows normal temperature ranges by age:

0 ~2 years	36.4 ~ 38.0 °C	97.5 ~ 100.4 °F
3 ~ 10 years	36.1 ~ 37.8 °C	97.0 ~ 100.0 °F
11 ~ 65 years	35.9 ~ 37.6 °C	96.6 ~ 99.7 °F
> 65 years	35.8 ~ 37.5 °C	96.4 ~ 99.5 °F

The range of normal temperature varies from person to person and can be influenced by many factors such as time of day, level of activity, medications, emotion and so on.

So we recommend that you practice with the thermometer on yourself and family members when you are healthy. This way you know how the thermometer works and you feel more confident of the measurements you take when a family member is ill.


4 before first use

Before you use the appliance for the first time, please do as follows: carefully unpack the appliance and remove all the packaging material. Keep the material (plastic bags and cardboard) out of reach of children. Check the appliance after unpacking for any damage which may have occurred during transport.

1. To achieve accurate measurements, it is very important to check that the infrared sensor is clean before using. To clean the sensor, gently wipe its surface with a cotton swab that is slightly moistened with alcohol and immediately wipe dry with a clean cotton swab. After cleaning, allow at least 20 minutes drying time before taking temperatures.
2. Earwax in the ear canal or lipid on the forehead may cause an inaccurate measurement, so make sure the ear canal or forehead is clean. If you clean the ear or forehead, please wait for 5-10 minutes before measuring.

- If the temperature of the storage area of the thermometer differs a lot from that of the measuring area, please wait 30 minutes before measuring.
- Please open the battery cover and remove the insulating piece for the batteries when you use the thermometer for the first time. The thermometer turns on automatically and will give a beep signal.


Caution:

In waiting mode, the LCD circularly displays the time and the current room temperature.

5 how to measure forehead temperature

- Put on forehead cap before use. Place the thermometer in your hand as shown with your thumb on the **[Forehead Measuring]** button.
 - Do not press the button yet.
 - Make sure the probe is clean before measuring.
 - If the measurement on the forehead is done without the forehead cap, the measurement will be incorrect.
- Place the thermometer gently on the forehead so that the forehead cap is centered between the eyebrows and the hairline. The forehead cap should touch the skin. Press the **[Forehead Measuring]** button to start the measurement.
- Keep the thermometer flat against the forehead until you hear the temperature broadcast or the sound signal. The measurement is complete and the thermometer may be removed from the head.
- Read the temperature on the display.

The innovative infrared technology measures temperature in 2 seconds, whether you take ear or forehead temperature.

6

how to measure ear temperature

1. Remove the forehead cap before use. Place the thermometer in your hand as shown with your finger on the **[Ear measuring]** button.
 - Make sure the probe is clean before measuring.
 - If the measurement is done with the forehead cap, the measurement will be incorrect.
2. Insert the probe very gently and slowly into the ear canal.
 - Hold the thermometer so that the probe faces straight in the direction of the eardrum.

It is important to point the probe window of the thermometer towards the eardrum and its periphery inside the ear.

 - Hold the outer ear and gently pull it toward the rear of the head to straighten the ear canal.
 - Always insert the thermometer into the same ear with the same direction and depth.
3. Press the **[Ear measuring]** button to start the measurement. Keeping the probe flat until you hear the temperature broadcast or the sound signal.
4. Remove the thermometer from the ear canal. The LCD displays the measured temperature. You can measure again after the symbol "000" disappears.

7

temperature taking tips

1. A temperature measurement taken in the right ear may differ from a measurement taken in the left ear. Therefore, always take the temperature in the same ear.
2. External factors may influence ear and forehead temperatures, particularly when an individual has:
 - been lying on one ear or the other
 - had their ears covered
 - been exposed to very hot or very cold temperatures, or

- been recently swimming or bathing.

In these cases, remove the individual from the situation and wait for 30 minutes prior to taking a temperature again.

3. For persons wearing hearing aids or ear plugs, remove the device and wait 30 minutes prior to taking a temperature.
4. Use the untreated ear if ear drops or other ear medications have been placed in the ear canal.


8 fever indication

When the measuring temperature is in the range of 37.6°C ~ 42.9°C (99.6°F ~ 109.3°F), there is a red light under the **[Forehead Measuring]** button for fever indication.

9 check the memory

There are 30 memory positions for recording ear and forehead temperature measurements.


The current measurement is always stored in the last storage space. When all storage spaces have been occupied, the oldest measurement is always deleted from the memory.

Press **[MEM]** button to check the temperature and date like right figure.

Press **[MEM]** button again to check temperature measurement of another memory position.

Press and hold **[MEM]** button to check the memory data fleetly.

If you do not press **[MEM]** button for 5 seconds, the display will go back to waiting mode.

Caution: when replacing the batteries, all memory will be lost.




10 date and time setting

Press and hold **[MEM]** button about 2 seconds to enter into setting mode.

Press **[MEM]** button to increase the parameter which you selected. Press and hold **[MEM]** button to increase the parameter fleetly.

Press **[TALK]** button to switch the settings of hour, minute, year (2001-2099), month, day and sleep mode (ON or OFF).

When the **sleep mode** is active, the thermometer will turn off automatically about 1min or 3min (which is factory default); when the sleep mode is deactivated, the thermometer remains power on. To save power, we recommended to set the sleep mode ON.

11 voice setting


Press **[MEM]** button to turn on the thermometer.

Press and hold **[TALK]** button about 2 seconds to select the talking language you want, it will change to another language or no-talking function [--] automatically.

Keep repeating this step until the right voice setting is selected. The display shows: HL - EN- GE - FR - [--].

CELSIUS/ FAHRENHEIT OPTION

Press and hold **[MEM]** and **[Forehead Measuring]** button to switch between Celsius (°C) and Fahrenheit (°F) when the thermometer in waiting mode, measuring mode or memory mode.



12 replacing the battery

When the low battery indicator “” appears on the LCD, or the thermometer does not function at all, you should replace the batteries as soon as possible.

The thermometer is supplied with two batteries (2xAAA). Remove the battery cover and take out the batteries. Insert new AAA batteries into the battery compartment.

Make sure you are installing the batteries properly, be sure to match the positive and negative ends of the battery to the marking in the battery compartment of the thermometer.

Replace the battery compartment cover and press to close.

Caution:

- To protect the environment, dispose of exhausted batteries according to national or local regulations.
- Keep batteries out of the reach of children.

13 specifications

Power supply	DC 3V (2x AAA)
Body measuring range	32.0°C ~ 42.9°C (89.6°F ~ 109.3 °F)
Laboratory accuracy	± 0.2°C (± 0.4°F) from 35.0°C ~ 42.0°C (95.0°F ~ 107.6°F) ± 0.3°C (± 0.5°F) for other range
Displayed room temperature range	5.0°C ~ 59.9°C (41.0°F ~ 139.8°F)
Room temperature accuracy	± 1°C (± 2°F)
Dimensions	124 mm x 38 mm x 66 mm (LxWxD)
Measuring time	1 - 2 seconds
Measuring interval	10 seconds
Weight	60 g (without batteries)
Forehead operating condition	15.0°C ~ 40.0°C (59.0°F ~ 104.0°F) with a relative humidity of 15% ~ 85% Atmosphere pressure: 700hPa~1060hPa
Ear operating condition	10.0°C ~ 40.0°C (50.0°F ~ 104.0°F) with a relative humidity of 15% ~ 85% Atmosphere pressure: 700hPa~1060hPa
Storage and transportation condition	-20°C ~ 55°C (-4°F ~ 131°F) with a relative humidity of 15% ~ 90% Atmosphere pressure: 700hPa~1060hPa
Service life	2 years
Service life of battery	With new batteries approx. 6 months bases on the use frequency at 5 times a day.
IP classification system	IP22
Classification of ME equipment	Internally powered/ Continuous operation

Caution:

This thermometer has been thoroughly tested and inspected to assure proper performance and operation!

14 cleaning & maintenance

To ensure accurate measurements, keeping the infrared sensor tip clean and free of scratches is very important. Finger prints, earwax or dirt will affect the accuracy of the thermometer.

In order to get an accurate measurement, please clean the probe tip.

Gently wipe its surface with a cotton swab slightly moistened with alcohol and immediately wipe dry with a clean cotton swab.

Use a soft dry cloth to clean the body of the thermometer. Never clean the thermometer with an abrasive cleanser, thinner, benzene or submerge the thermometer into water or other liquids.



After cleaning, place the thermometer in the table stand holder. Store it in a clean, dry place at room temperature. Never expose the thermometer to extreme temperatures, humidity, direct sunlight or shock.

CALIBRATION

The thermometer is initially calibrated at the time of the manufacture. If the thermometer is used according to the instruction, periodic recalibration is not required. If at any time you question the accuracy of the measurement, please contact the retailer immediately.

Don't attempt to modify or reassemble the thermometer.

15 troubleshooting

Troubles	Checklist	Countermeasures
No response/ Automatic reset when pulling out the insulator	Batteries used up?	Change the batteries
	Batteries in wrong polarity?	Take out the batteries, reinsert batteries correctly
	Poor battery contact	
	The measured temperature is lower than 32.0°C/89.6°F or higher than 42.9°C/109.3°F. Please check the operation method.	Follow the steps in the user manual for proper measurement.
	Hardware problem	Contact the manufacturer or point of sale.
The measurement is not accurate or if there is any doubt of the measure result	Please check if the infrared sensor is clean or not	Clean the infrared sensor with a cotton swab according to the user manual.
	Please check that the measuring way is correct or not	Ensure you have read the manual and know how to use the thermometer properly.
	Please check if you have let the thermometer and patient get stabilized in the room for 30 minutes at least	Please keep the thermometer and patient in the measuring room at least 30 minutes before using.
	Are you using the thermometer indoor?	Please take the measurement indoor.
	Please check if you have hold the thermometer in your hand too long and affect the accuracy	Put the thermometer on the table in the room where the measurement is taking place and let it cool down first.

Troubles	Checklists	Countermeasures
	Please check if there is low battery icon on the LCD display	Replace with new batteries

16

Important information regarding electromagnetic compatibility (EMC)

This device should not be used adjacent to or stacked with other equipment and that if adjacent or stacked use is necessary, this device should be observed to verify normal operation in the configuration in which it will be used

Use of accessories other than those specified or provided by the manufacturer of this device could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the device, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

When the operating environment is relatively dry, strong electromagnetic interference usually occurs. At this time, the device may be affected as follows:

- the device stops output;
- the device turns off;
- the device restarts;

The above phenomenon does not affect the basic safety and essential performance of the device, and the user can use it according to the instruction. If you want to avoid the above phenomenon, please use it according to the environment specified in the manual.

Table 1

Guidance and manufacturer's declaration – electromagnetic emissions		
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The device is suitable for use in all establishments, other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Niet toepasbaar	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Niet toepasbaar	

Table 2

Guidance and manufacturer's declaration – electromagnetic immunity			
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	power supply lines: ± 2 kV input/output lines: ± 1 kV	Not applicable	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC61000-4-5	line(s) to line(s): ± 1 kV line(s) to earth: ± 2 kV 100 kHz repetition frequency	Not applicable	Mains power quality should be that of a typical commercial or hospital environment.

Table 2

Guidance and manufacturer's declaration – electromagnetic immunity			
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment -guidance
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0%UT; 0,5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0%UT ; 1 cycle and 70%UT ; 25/30 cycles Single phase: at 0° 0% UT ; 300 cycle	Not applicable	Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.
Power frequency (50Hz/60Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE UT is the a.c. mains voltage prior to application of the test level.			

Table 3**Guidance and manufacturer's declaration – electromagnetic immunity**

The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.

Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 V 0.15MHz to 80MHz 6 V in ISM and amateur radio bands between 0.15 MHz and 80 MHz	Not appli- cable	Portable and mobile RF communications equipment should be used no closer to any part of the device than the recommended separation distance calculated from the equation appropriate for the frequency of the transmitter. Recommended separation distances $d=1,2\sqrt{P}$ 150 KHz to 80 MHz $d=1,2\sqrt{P}$ 80 MHz to 800 MHz $d=2,3\sqrt{P}$ 80 MHz to 2.7 GHz
Radiated RF IEC 61000-4-3	10 V/m, 80 MHz to 2.7 GHz	10 V/m	where, P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer, d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: 

Table 3

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people..

- a Field strengths from fixed transmitters, such as base stations for radio (cellular / cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the device.
- b Over the frequency range 0.15 MHz to 80 MHz, field strengths should be less than 3V/m.

Table 4**Recommended separation distances between portable and mobile RF communications equipment and the device.**

The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)		
	0,15 MHz to 80 MHz $d = 1,2\sqrt{P}$	80 MHz to 800 MHz $d = 1,2\sqrt{P}$	80 MHz to 2,7 GHz $d = 2,3\sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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We do not need to remind you of the importance of service. After all, we develop our products to a standard so that you can enjoy them for many years, without any concerns. If, nevertheless, there is a problem, we believe you are entitled to a solution straight away. Hence our products come with an exchange service, on top of the rights and claims you are entitled to by law. By exchanging a product or part, we save you time, effort and costs.

2-year full manufacturer's warranty

1. Customers enjoy a 2-year full manufacturer's warranty on all Inventum products. Within this period, a faulty product or part will always be exchanged for a new model, free of charge. In order to claim under the 2-year full manufacturer's warranty, you can either return the product to the shop you bought it from or contact the Inventum customer service department via the form at www.inventum.eu/service-aanvraag.
2. The 2-year warranty period starts from the date the product is bought.
3. In order to claim under the warranty, you must produce a copy of the original receipt.
4. The warranty applies only to normal domestic use of Inventum products within the Netherlands.

Breakdowns or faults outside the warranty period

1. Breakdowns or faults in small or large domestic appliances outside the warranty period, can be reported to the customer services department via the form at www.inventum.eu/service-aanvraag or by calling the customer services department.
2. The customer services department may ask you to send the product for inspection or repair. The costs of dispatch will be at your expense.
3. The inspection to establish whether repair is possible involves a charge. You need to grant your permission for this, in advance.
4. In the event of a large domestic appliance, Inventum, at your request, can send out a service engineer. In that case, you will be

charged the call-out costs, as well as parts and labour.

5. In the event of an instruction to repair, the repair costs must be paid in advance. In the event of a repair by a service engineer, the costs of the repair must be settled with the engineer onsite, preferably by means of PIN payment.

Warranty exclusions

1. The following is excluded from the aforesaid warranties:
 - normal wear and tear;
 - improper use or misuse;
 - insufficient maintenance;
 - failure to comply with the operating and maintenance instructions;
 - unprofessional installation or repairs by third parties or the customer himself;
 - non-original parts used by the customer;
 - use for commercial or business purposes;
 - removal of the serial number and/or rating label.
2. In addition, the warranty does not apply to normal consumer goods, such as:
 - dough hooks, baking tins, (carbon) filters, etc.;
 - batteries, bulbs, carbon filters, fat filters etc.;
 - external connection cables;
 - glass accessories and glass parts such as oven doors;
 - and similar items.
3. Transport damage not caused by Inventum is also excluded. Therefore, inspect your new device before starting to use it. If you detect any damage, you must report this to the store where you purchased the product within 5 working days, or to the Inventum customer service department via the form at www.inventum.eu/service-aanvraag. If transport damage is not reported within this period, Inventum does not accept any liability in this respect.
4. The following are excluded from warranty and/or replacement: faults, loss of and damage to the device as a result of an event that is normally insured under the home contents insurance.

Important to know

1. The replacement or repair of a faulty product, or a part thereof, does not lead to an extension of the original warranty period.
2. If a complaint is unfounded, all costs arising from it will be at the customer's expense.
3. Following expiry of the warranty period, all costs of repair or replacement, including administration costs, dispatch and call-out charges, will be charged to the customer.
4. Inventum cannot be held liable for damage caused by external events, unless this liability arises from mandatory statutory provisions.
5. These warranty and service provisions are governed by Dutch law. Disputes will be settled exclusively by the competent Dutch court.

Shenzhen Dongdixin Technology Co., LTD.

Floor 1-2, No. 3 Building, Fanshen Xusheng Industrial Estate,
Xilixiaobaimang, Nanshan District, 518108, Shenzhen China

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Inventum Huishoudelijke Apparaten B.V.
Shanghai International Holding Corp. GmbH (Europe)