DE	Digital TENS/EMS <i>Gebrauchsanweisung</i>	2
EN	Digital TENS/EMS unit <i>Instructions for use</i>	20
FR	Electrostimulateur TENS/EMS numérique <i>Mode d'emploi</i>	37
IT	TENS/EMS digitale <i>Istruzioni per l'uso</i>	55
RU	Тренажер TENS/EMS для мышц <i>Инструкция по применению</i>	73
PL	Cyfrowe urządzenie do elektrostymulacji mięśni i nerwów TENS/EMS <i>Instrukcja obsługi</i>	92
NL	Digitaal TENS/EMS-apparaat <i>Gebruiksaanwijzing</i>	110

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
Read these instructions for use carefully and keep them for later use. Make them accessible to other users and note the information they contain.

Table of contents

1. Included in delivery	21	8.1 Notes on use	27
2. Getting to know your device	21	8.2 Starting use	27
3. Signs and symbols	21	8.3 Program overview	28
4. Intended use	22	8.4 Information regarding the positioning of electrodes	30
4.1 TENS	22	8.5 Customisable programs	31
4.2 EMS/massage	23	8.6 Doctor's function	33
5. Warnings and safety notes	24	8.7 Electric current parameters	33
5.1 Damage	25	9. Cleaning and maintenance	34
5.2 Measures for handling batteries	26	10. Accessories and replacement parts	35
5.3 Notes on electromagnetic compatibility	26	11. What if there are problems?	35
6. Device description	26	12. Technical specifications	35
7. Initial use	27	13. Disposal	36
8. Usage	27	14. Warranty/service	36

1. Included in delivery

Check that the set packaging has not been tampered with and make sure that all contents are present. Before use, ensure that there is no visible damage to the device or accessories and that all packaging material has been removed. If you have any doubts, do not use the device and contact your retailer or the specified Customer Services address.

A	1 x Digital TENS/EMS device (including belt clip)
B	2 x connection cables
C	4 x adhesive electrodes (45 x 45 mm)
D	3 x AAA batteries

2. Getting to know your device

What is the Digital TENS/EMS device and what can it do?

Digital TENS/EMS devices fall into the electrostimulation device category. They provide three basic functions suitable for combined operation:

1. Electrical stimulation of nerve tracts (TENS)
2. Electrical stimulation of muscle tissue (EMS)
3. A massage effect triggered by electrical signals.

The device also features two independent stimulation channels and four self-adhesive electrodes. It offers a wide range of functions for increasing general well-being, pain relief, maintaining physical fitness, relaxation, muscle revitalisation and combating tiredness. For these purposes, you can either choose from pre-set programs or specify your own to suit your individual needs.

The operating principle of electrostimulation devices is based on the imitation of impulses in our bodies that are transferred to nerve and muscle fibres using electrodes via our skin. The electrodes can be applied to many parts of the body; the electrical impulses are completely harmless and virtually painless. In certain applications you will only perceive a slight tingling or vibrating sensation. The electrical impulses that are sent into the tissue influence the transmission of stimulation into nerves, nerve centres and muscle groups in the application area. Electrical muscle stimulation (EMS) is a widespread and generally recognised method, and has been used in sports and rehabilitation medicine for years.

Electrostimulation usually only has an effect after regular applications. With regard to muscles, electrostimulation does not replace regular training. However, it is a useful, supplementary training element.

3. Signs and symbols

The following symbols are used on the device, in these instructions for use, on the packaging and on the type plate:

	WARNING Warning notice indicating a risk of injury or damage to health
	IMPORTANT Safety note indicating possible damage to the device/accessory
	Note Note on important information

	Observe the instructions for use
	Protect from moisture
SN	Serial number
	Isolation of applied parts Type BF Galvanically isolated application part (F stands for “floating”); meets the requirements for leakage currents for type B
	Disposal in accordance with the Waste Electrical and Electronic Equipment EC Directive – WEEE
CE	This product satisfies the requirements of the applicable European and national directives.
	Manufacturer
	The device can emit effective output values above 10 mA, averaged over every 5-second interval
Storage/Transport 	Permissible storage and transport temperature and humidity
Operating 	Permissible operating temperature and humidity
MD	Medical device

	The device may not be used by people with medical implants (e.g. pacemakers), as this may affect their functionality.
	Separate the product and packaging elements and dispose of them in accordance with local regulations.
	Marking to identify the packaging material. A = Material code, B = Material number: 1-7 = Plastics, 20-22 = Paper and cardboard
	Separate the packaging elements and dispose of them in accordance with local regulations.
	Do not dispose of batteries containing harmful substances with household waste
UDI	Unique Device Identifier (UDI) for unique product identification
REF	Item number

4. Intended use

4.1 TENS

4.1.1 Intended use for TENS

The device is intended to treat pain using TENS technology (transcutaneous electrical nerve stimulation).

4.1.2 Target group for TENS

This device is suitable for self-treatment by adults with pain, provided that they have discussed the cause and treatment options with their doctor.

4.1.3 Indications/clinical benefits of TENS

TENS, or transcutaneous electrical nerve stimulation, refers to the electrical stimulation of nerves through the skin. TENS is an effective non-pharmacological method for treating different types of pain

from a variety of causes. It has no side effects if administered correctly. The method has been clinically tested and approved and can be used for simple self-treatment. The pain-relieving or pain-suppressing effect is achieved by inhibiting the transfer of pain to nerve fibres (caused mainly by high-frequency impulses) and by increasing the secretion of endorphins in the body. Their effect on the central nervous system reduces the sensation of pain. The method has been scientifically substantiated and approved as a form of medical treatment. Any symptoms that could be relieved using TENS must be checked by your doctor. Your doctor will also give you instructions on how to carry out a TENS self-treatment regime.

TENS is clinically tested and approved to treat the following complaints:

- Back pain, particularly in the lumbar/cervical spine area.
- Sore joints (e.g. knee, hip and shoulder joints).
- Neuralgia.
- Menstrual cramps in women.
- Pain resulting from injury to the musculoskeletal system.
- Pain caused by circulatory disorders.
- Chronic pain with various causes.

4.2 EMS/massage

4.2.1 Intended use for EMS/massage

Furthermore, the device is intended to allow people to use EMS (electrical muscle stimulation) to strengthen the muscles, for regeneration and for relaxing massage.

4.2.2 Target group for EMS/massage

This device is suitable for self-treatment by adults who want to strengthen, regenerate or massage their muscles. Users should also discuss use of the device with their doctor.

4.2.3 Indications for EMS/massage

In sports and fitness, electrical muscle stimulation (EMS) is used to complement conventional muscle training, to increase the performance of muscle groups and to adjust physical proportions to achieve the desired aesthetic results, amongst other things. There are two different types of EMS application. One is for targeted

strengthening of the muscles (activating application), and the other is to achieve a relaxing, restful effect (relaxing application).

The activating application involves:

- Muscle training to increase endurance and/or
- Muscle training to support strengthening of certain muscles
- Or muscle groups to achieve desired changes in body proportions.

The relaxing application involves:

- Muscle relaxation for easing muscle tension.
- Improving symptoms of muscular fatigue.
- Acceleration of muscle regeneration after intense muscle performance (e.g. after a marathon).

Thanks to integrated massage technology, the device is also able to relieve muscle tension and combat fatigue with a program based on the sensation and effects of a real massage. The positioning suggestions and program tables in these instructions for use allow you to quickly and simply determine the corresponding application (depending on the affected area of the body) and set the device to achieve the desired effects. Thanks to the two separately adjustable channels, the Digital TENS/EMS device offers you the advantage of being able to set the intensity of the impulses independently from each other for two treatment areas on the body, for example to cover both sides of your body or to evenly stimulate larger areas of tissue. The option to individually set the intensity of each channel also enables you to treat two separate areas of the body simultaneously instead of having to treat the individual areas in turn. This saves you time.

5. Warnings and safety notes

Contraindications

To avoid damage to health, we strongly advise against using the Digital TENS/EMS device in the following situations:

- If you have implanted electrical devices (e.g. a pacemaker).
- If you have metal implants.
- If you use an insulin pump.
- If you have a high temperature (e.g. > 39°C).
- If you have a known or acute cardiac arrhythmia or disorders of the heart's impulse and conduction system.
- On acutely or chronically diseased (injured or irritated) skin (e.g. inflamed skin – whether painful or not, reddened skin, rashes, e.g. allergies, burns, bruises, swellings, both open and healing wounds, and post-operative scars where the healing process could be affected).
- If you suffer from a seizure disorder (e.g. epilepsy).
- If you are pregnant.
- If you have cancer.
- After an operation, if strong muscle contractions could affect the healing process.
- If you are connected to a high-frequency surgical device.
- In the case of acute or chronic diseases of the gastrointestinal tract.

General warnings

This device is not a substitute for medical consultation and treatment. Consult your doctor first if you are experiencing any pain or are suffering from an illness. Before using the device, consult your doctor if any of the following apply to you:

- If you suffer from an acute medical condition, in particular if you suspect or have been diagnosed with high blood pressure, a blood coagulation disorder, propensity to thrombo-embolic conditions or recurrent malignant growths.
- If you have any skin conditions.
- If you have unexplained chronic pain in any part of the body.


- Diabetes.
- Any sensory impairment that reduces the feeling of pain (e.g. metabolic disorders).
- If you are receiving medical treatment.
- In the event of complaints linked to the stimulation treatment.
- If you suffer from persistently irritated skin due to long-term stimulation at the same electrode site.

Only use Digital TENS/EMS device:

- On people.
- For the purpose for which it was developed and as specified in these instructions for use. Any improper use may be dangerous.
- For external use.
- With the original accessories supplied, which can be re-ordered. Failure to do this invalidates the warranty.
- This device may be used by children over the age of 8 and by people with reduced physical, sensory or mental skills or a lack of experience or knowledge, provided that they are supervised or have been instructed on how to use the device safely, and are fully aware of the consequent risks of use.
- Children must not play with the device.
- The device is only intended for domestic/private use, not for commercial use.

Do not use the device:

- On the skeletal skull structure, or around the mouth, throat or larynx.
- In the neck/carotid artery area.
- The stimulation should not be applied across or through the head, directly onto the eyes, covering the mouth, to the front of the neck (especially the carotid sinus), or with the electrodes placed on the chest and upper back or across the heart.
- In the genital area.
- The device must never be used near to the heart. The stimulation electrodes must not be placed on any part of the front ribcage (where the ribs and breast-bone are located), especially not on the two large

pectorals. This can increase the risk of ventricular fibrillation and induce cardiac arrest.

General precautions

- Do not use in humid environments (e.g. in the bathroom) or when bathing or showering.
- Do not use after consuming alcohol.
- Always pull gently on electrodes to remove them from the skin to prevent injuries in the unusual case of highly sensitive skin.
- Hold the device away from sources of heat and do not use it in close proximity (approx. 1 m) to shortwave or microwave devices (e.g. mobile phones), as doing so can result in unpleasant current peaks.
- Do not expose the device to direct sunlight or high temperatures.
- Protect the device from dust, dirt and humidity.
- Never immerse the device in water or other liquids.
- The device is suitable for self-treatment.
- For hygiene reasons, the electrodes may only be used on one person.
- If the device does not work properly, or you feel unwell or experience pain, stop using it immediately.
- Switch off the device or the respective channel first before removing or moving the electrodes to prevent unintentional stimulation.
- Do not modify electrodes (e.g. by cutting them). This increases the current density, which is potentially hazardous (max. recommended output value for the electrodes is 9 mA/cm², an effective current density beyond 2 mA/cm² requires increased awareness).
- Do not use whilst asleep, driving a vehicle or operating machinery.
- Do not use whilst undertaking any activity where an unexpected reaction (e.g. strong muscle contractions even at low intensity) could be dangerous.
- Ensure that no metallic objects (e.g. belt buckles or necklaces) come into contact with the electrodes during stimulation. If you are wearing jewellery or have piercings in the area to be treated (e.g. a navel piercing), these must be removed before using the device. Failure to do so could result in spot burns.
- Keep the device away from children to prevent potential hazards.

- Make sure not to confuse the electrode cables including contacts with your headphones or other devices and do not connect the electrodes to other devices.
- Do not use the device whilst using other devices that transmit electrical impulses to your body.
- Do not use in the vicinity of highly flammable substances, gases or explosives.
- Do not use rechargeable batteries. Always use the same battery types.
- During the initial few minutes, use the device while sitting or lying down to minimise the risk of accidental injuries as a consequence of isolated cases of vagal responses (feeling of faintness). If you feel faint, switch off the device immediately, lie down and support the legs in an elevated position (approx. 5–10 minutes).
- Treatment of the skin with moisturising lotions or ointments beforehand is not recommended as this considerably increases the electrode wear and may cause unpleasant current peaks.
- The treatment should be comfortable. If the device does not work properly, or you feel unwell or experience pain, stop using it immediately.
- Keep packaging material away from children (risk of suffocation).

5.1 Damage

- If the device is damaged, do not use it and contact your retailer or the specified Customer Services address.
- To ensure that the device functions effectively, do not drop it or dismantle it.
- Check the device for signs of wear and tear or damage. If there are such signs of wear and tear or damage or if the device was used improperly, it must be returned to the manufacturer or retailer before further use.
- Switch the device off immediately if it is faulty or not working properly.
- Never attempt to open and/or repair the device yourself. Repairs may only be carried out by Customer Services or authorised retailers. Failure to comply with this instruction will void the warranty.

- The manufacturer is not liable for damage resulting from improper or incorrect use.

5.2 Measures for handling batteries

- If your skin or eyes come into contact with battery fluid, wash the affected area with water and seek medical assistance.
- Choking hazard! Small children may swallow and choke on batteries. Therefore, batteries should be stored out of the reach of small children.
- Observe the plus (+) and minus (-) polarity signs.
- If a battery has leaked, put on protective gloves and clean the battery compartment with a dry cloth.
- Protect batteries from excessive heat.
- Risk of explosion! Do not throw batteries into a fire.
- Do not charge or short-circuit the batteries.
- If the device is not to be used for a long period of time, remove the batteries from the battery compartment.
- Use identical or equivalent battery types only.
- Always replace all batteries at the same time.
- Do not use rechargeable batteries.
- Do not disassemble, open or crush the batteries.

5.3 Notes on electromagnetic compatibility


- The device is suitable for use in all environments listed in these instructions for use, including domestic environments.
- The use of the device may be limited in the presence of electromagnetic disturbances. This could result in issues such as error messages or the failure of the display/device.
- Avoid using this device directly next to other devices or stacked on top of other devices, as this could lead to faulty operation. If, however, it is necessary to use the device in the manner stated, this device as well as the other devices must be monitored to ensure they are working properly.
- The use of accessories other than those specified or provided by the manufacturer of this device can lead to an increase in electromagnetic emissions or a decrease in the device's electromagnetic immunity; this can result in faulty operation.

- Failure to comply with the above can impair the performance of the device.

6. Device description

Buttons:

- 1 ON/OFF button 
- 2 E button (Enter)
- 3 Intensity setting buttons (**Ch1 +/-** left, **Ch2 +/-** right)
- 4 Keylock 
- 5 Selection button **UP** and **DOWN**
- 6 M button (menu)

Display (full screen):

- 1 Menu **TENS** / **EMS** / **MASSAGE**
- 2 Low battery
- 3 Keylock
- 4 Timer function (remaining time display) or working time
- 5 Frequency display (Hz), pulse width (μ s) or pause time
- 6 Impulse intensity, channel 1 (**Ch1**)
- 7 Impulse intensity, channel 2 (**Ch2**)
- 8 Operating status display
- 9 Program number

7. Initial use

1. Remove the belt clip from the device (if attached).
2. Press the battery compartment cover on the rear of the device and slide it downwards.
3. Insert the 3 alkaline AAA 1.5V batteries. Make sure that the batteries are inserted the correct way round in accordance with the markings.
4. Close the battery compartment lid again carefully (Fig. 1).
5. Reattach the belt clip, if required.
6. Connect the connection cables with the electrodes (Fig. 2).

-  The electrodes feature clip closures to make connection particularly easy.
7. Guide the connection cable plugs into the socket on the top of the device (Fig. 3).
 8. Do not pull, twist or kink the cables (Fig. 4).

-  Please note that when the battery is replaced or removed, all settings are restored to the factory default settings.

8. Usage

8.1 Notes on use

- The device switches itself off automatically if you do not use it for 2 minutes (automatic switch-off). When the device is switched on again, the LCD screen displays the menu selection and the most recently used menu flashes.
- A brief acoustic signal is output when a valid button is pressed. Two brief acoustic signals are output when an invalid button is pressed.
- You can pause the stimulation at any time by briefly pressing the ON/OFF button . To continue the stimulation, briefly press the ON/OFF button  again and set the desired impulse intensity again.


Fig. 1


Fig. 2

Fig. 3


Fig. 4

8.2 Starting use

Step 1: Look for a suitable program from the program tables (see section “8.3 Program overview”).

Step 2: Place the electrodes on the desired area for treatment (for positioning suggestions see section “8.4 Information regarding the positioning of electrodes”) and connect them to the device.

Step 3: Press the ON/OFF button  to switch on the device.

Step 4: Press the **M** button to navigate through the **TENS** / **EMS** / **MESSAGE** menus and press the **E** button to confirm your selection.

Step 5: Use the **UP/DOWN** setting buttons to select the program number you want and press the **E** button to confirm your selection. At the start of the stimulation treatment, the impulse intensity of **Ch1** and **Ch2** is set to 00 by default. No impulses are sent to the electrodes yet.

Step 6: Use the intensity setting buttons for **Ch1** and **Ch2** to select the impulse intensity you want. The indicator for impulse intensity in the display changes accordingly. If the program is in a pause phase, the intensity cannot be increased.

General information

If you wish to return to the previous selection menu, press the **M** button. By pressing and holding the **E** button, you can skip individual setting steps and start directly with the stimulation treatment.

Keylock

Locks the buttons to avoid them being pressed unintentionally.

1. To activate the keylock, hold down the  button until the symbol is visible in the display (approx. 3 seconds).
2. To deactivate the keylock, hold down the  button until the symbol disappears from the display (approx. 3 seconds).

8.3 Program overview

The Digital TENS/EMS device features a total of over 70 programs:

- 10 TENS programs
- 30 EMS programs
- 10 MASSAGE programs

In all programs you can set the impulse intensity of both channels individually.

You can also set various parameters in TENS programs 8-10 and EMS programs 28-30 to adjust the stimulating effect to the application area.

8.3.1 TENS program table

Progr. no.	Practical areas for application, indications	Phase	Running time (min)
1	Pain relief – acute pain	1	30
2	Pain relief – chronic pain (pulse width modulated)	1	30
3	Endorphin effect (burst)	1	30
4	Sciatica	1	30
5	Atrophy treatment	1	05
		2	15
6	Lumbago	1	20
		2	20
7	Periarthritis	1	15
		2	10

 TENS programs 8–10 can be set individually (see section “8.5 Customisable programs”).

Note: See section 8.4 for the correct electrode position.

8.3.2 EMS program table

Progr. no.	Practical areas for application, indications	Possible electrode positions	Phase	Running time (min)
1	Capillarisation	1–28	1	20
2	Warming up	1–28	1	10
3	Cooling down after training/match	1–28	1	20
4	Maximum power of lower limbs	22, 23, 24, 26, 27	1	05
			2	15
			3	05
5	Resistance of lower limbs	22, 23, 24, 26, 27	1	05
			2	20
			3	05
6	Explosive power of lower limbs	22, 23, 24, 26, 27	1	05
			2	12
			3	05
7	Maximum power of torso and upper limbs	1–20	1	05
			2	15
			3	05
8	Resistance of torso and upper limbs	1–20	1	05
			2	12
			3	05
9	Explosive power of torso and upper limbs	1–20	1	05
			2	12
			3	05
10	Lipolysis of lower limbs	22, 23, 24, 26, 27	1	40
11	Toning of lower limbs	22, 23, 24, 26, 27	1	20
			2	10
12	Strengthening of lower limbs	22, 23, 24, 26, 27	1	20
			2	10
13	Shaping of lower limbs	22, 23, 24, 26, 27	1	20
			2	20

Progr. no.	Practical areas for application, indications	Possible electrode positions	Phase	Running time (min)
14	Increasing size of lower limbs	22, 23, 24, 26, 27	1	15
15	Lipolysis of stomach muscles	18,19, 20	1	40
16	Toning of stomach muscles	18,19, 20	1	20
			2	10
17	Strengthening of stomach muscles	18,19, 20	1	20
			2	10
18	Shaping of stomach muscles	18,19, 20	1	20
			2	20
19	Toning of upper limbs	12–17	1	20
			2	10
20	Strengthening of upper limbs	12–17	1	20
			2	10
21	Shaping of upper limbs	12–17	1	20
			2	20
22	Increasing size of upper limbs	12–17	1	15
23	Toning of hips and thighs	20, 23, 24	1	20
			2	10
24	Strengthening of hips and thighs	20, 23, 24	1	20
			2	10
25	Toning of gluteal muscles	22	1	20
			2	10
26	Strengthening of gluteal muscles	22	1	20
			2	10
27	Shaping of gluteal muscles	22	1	20
			2	20

 EMS programs 28–30 can be set individually (see section “8.5 Customisable programs”).

Note: Refer to section 8.4 for the correct electrode position.

8.3.3 MESSAGE program table


Progr. no.	Practical areas for application, indications	Phase	Running time (min)
1	Muscle strain relief	1	20
2	Relaxing massage	1	25
3	Invigorating massage	1	20
4	Tapping massage	1	15
5	Tapping and jet massage	1	15
6	Vibrating massage	1	20
7	Warm up massage 1	1	20
8	Warm up massage 2	1	10
9	Pressure massage	1	20
10	Kneading and pressure massage	1	20

Note: See section 8.4 for the correct electrode position.

WARNING!

Do not apply the electrodes to the front wall of the chest, i.e. do not massage the large left and right pectoral muscles.

8.4 Information regarding the positioning of electrodes

It is fundamental to the intended success of electrostimulation applications that electrodes are sensibly positioned.

We recommend that you consult your doctor to establish the ideal electrode positions for your intended application area.

The following applies to the selection of electrode positions:

Electrode distance


The greater the distance between electrodes, the larger the stimulated tissue volume. This applies to the area and depth of the tissue volume. At the same time, however, the stimulation intensity of the tissue decreases the further the electrodes are apart. As a result, greater distances between electrodes mean a larger tissue volume is stimulated, but less intensively. Consequently, you must increase the impulse intensity to boost stimulation.

The following guidelines apply to the selection of the electrode distances:

- Sensible distance: approx. 5–15 cm
- At distances below 5 cm, the device primarily stimulates surface structures intensively
- At distances in excess of 15 cm, large areas and deep structures are stimulated very weakly

Relation between electrodes and muscle fibre structures

Adapt the current flow direction to the fibre structure of the muscle according to the muscle layer you would like to treat. If you are targeting superficial muscles, position the electrodes parallel to the fibre structure (A–B / C–D) and if you are targeting deeper layers of tissue, position the electrodes across the fibre structure. You can do this by positioning electrodes as crosses (i.e. diagonally), such as A–D / B–C.

- ❗ As part of pain relief treatment (TENS) using the Digital TENS/EMS device and its 2 separately adjustable channels and 2 adhesive electrodes each, it is advisable either to position the electrodes of a channel so that the area affected by the pain is between the electrodes or to position one electrode directly on the

area affected by the pain and the other electrode at a minimum distance of 2–3 cm.

You may use the electrodes of the second channel to simultaneously treat additional areas affected by pain or use them in conjunction with the electrodes of the first channel to surround the area affected by pain (position electrodes opposite). In this case, we once again recommend positioning electrodes as crosses.

- ❗ Tip for the massage function: always use all 4 electrodes for optimum treatment.

Use the electrodes on clean skin, which is preferably free from hair and grease, in order to prolong the life of the electrodes. If required, clean the skin with water and remove hair prior to treatment.

- ❗ If an electrode should come loose during use, the impulse intensity of both channels is reduced to the lowest level. Apply the electrode again and reset the desired impulse intensity.

8.5 Customisable programs (Applies to TENS 8–10, EMS 28–30)

The TENS 8–10 and EMS 28–30 programs can be set individually according to your needs.

TENS 08 program

TENS 08 is a program that you can also customise. In this program, you can set the impulse frequency to between 1 and 150 Hz and the impulse width to between 50 and 250 μ s.

1. Place the electrodes on the desired area for treatment (for positioning suggestions see section “8.4 Information regarding the positioning of electrodes”) and connect them to the device.
2. Select the TENS 08 program as described in section “8.2 Starting use” (step 3 to step 5).
3. Use the **UP/DOWN** setting buttons to select the impulse frequency you want and press the **E** button to confirm.
4. Use the **UP/DOWN** setting buttons to select your preferred impulse width and press the **E** button to confirm.
5. Use the **UP/DOWN** setting buttons to select the treatment time you want and use the **E** button to confirm.

6. Use the left and right **INTENSITY** setting buttons respectively for **Ch1** and **Ch2** to select the impulse intensity you want.

TENS 09 program

The TENS 09 program is a **burst** program that you can also customise. Various impulse sequences run in this program. Burst programs are suitable for all areas of application to be treated with changing signal patterns (to minimise the level of becoming accustomed to the treatment). In this program you can set an impulse width of between 50 and 250 μ s.

1. Place the electrodes on the desired area for treatment (for positioning suggestions see electrode positions in section 8.4) and connect them to the device.
2. Select the TENS 09 program as described in section “8.2 Starting use” (step 3 to step 5).
3. Use the **UP/DOWN** setting buttons to select your preferred impulse width and press the **E** button to confirm.
4. Use the **UP/DOWN** setting buttons to select the treatment time you want and use the **E** button to confirm.
5. Use the left and right **INTENSITY** setting buttons respectively for **Ch1** and **Ch2** to select the impulse intensity you want.

TENS 10 program

TENS 10 is a program that you can also customise. In this program, you can set the impulse frequency to between 1 and 150 Hz. The impulse width changes automatically during the stimulation treatment.

1. Place the electrodes on the desired area for treatment (for positioning suggestions see electrode positions in section 8.4) and connect them to the device.
2. Select the TENS 10 program as described in section “8.2 Starting use” (step 3 to step 5).
3. Use the **UP/DOWN** setting buttons to select the impulse frequency you want and press the **E** button to confirm.
4. Use the **UP/DOWN** setting buttons to select the treatment time you want and use the **E** button to confirm.
5. Use the left and right **INTENSITY** setting buttons respectively for **Ch1** and **Ch2** to select the impulse intensity you want.

EMS 28 program

EMS 28 is a program that you can also customise. In this program, you can set the impulse frequency to between 1 and 150 Hz and the impulse width to between 80 and 320 μ s.

1. Place the electrodes on the desired area for treatment (for positioning suggestions see electrode positions in section 8.4) and connect them to the device.
2. Select the EMS 28 program as described in section “8.2 Starting use” (step 3 to step 5).
3. Use the **UP/DOWN** setting buttons to select the impulse frequency you want and press the **E** button to confirm.
4. Use the **UP/DOWN** setting buttons to select your preferred impulse width and press the **E** button to confirm.
5. Use the **UP/DOWN** setting buttons to select the treatment time you want and use the **E** button to confirm.
6. Use the left and right **INTENSITY** setting buttons respectively for **Ch1** and **Ch2** to select the impulse intensity you want.

EMS 29 program

EMS 29 is a program that you can also customise. In this program, you can set the impulse frequency to between 1 and 100 Hz. The impulse width changes automatically during the stimulation treatment.

1. Place the electrodes on the desired area for treatment (for positioning suggestions see electrode positions in section 8.4) and connect them to the device.
2. Select the EMS 29 program as described in section “8.2 Starting use” (step 3 to step 5).
3. Use the **UP/DOWN** setting buttons to select the impulse frequency you want and press the **E** button to confirm.
4. Use the **UP/DOWN** setting buttons to select the treatment time you want and use the **E** button to confirm.
5. Use the left and right **INTENSITY** setting buttons respectively for **Ch1** and **Ch2** to select the impulse intensity you want.

EMS 30 program

EMS 30 is a preset program that you can also customise. In this program, you can set the impulse frequency to between 1 and 100 Hz.

You can also set the working time and pause time for this program to between 1 and 30 seconds each.

1. Place the electrodes on the desired area for treatment (for positioning suggestions see electrode positions in section 8.4) and connect them to the device.
2. Select the EMS 30 program as described in section “8.2 Starting use” (step 3 to step 5).
3. Use the **UP/DOWN** setting buttons to select the on time you want and press the **E** button to confirm.
4. Use the **UP/DOWN** setting buttons to select the off time you want and press the **E** button to confirm.
5. Use the **UP/DOWN** setting buttons to select the impulse frequency you want and press the **E** button to confirm.
6. Use the **UP/DOWN** setting buttons to select the treatment time you want and use the **E** button to confirm.
7. Use the left and right **INTENSITY** setting buttons respectively for **Ch1** and **Ch2** to select the impulse intensity you want.

8.6 Doctor's function

The Doctor's function is a special setting to allow you to access your personal program even more easily and directly.

Your individual program settings are instantly recalled and activated when the device is switched on.

You may wish to adjust this individual program following advice from your doctor.

Setting the Doctor's function

- Select your program and the corresponding settings as described in section “8.2 Starting use”.
- At the start of the stimulation treatment, the impulse intensity of **Ch1** and **Ch2** is set to 00 by default. No impulses are sent to the electrodes yet. Before setting the impulse intensity you want using the intensity setting buttons, press and hold the **Ch2** - button for 5 seconds. Storage in the Doctor's function is confirmed by a long acoustic signal.

When you switch on the device again, the program you saved using the Doctor's function is automatically opened directly.

Deleting the Doctor's function

To clear the device again and to re-allow access to other programs, press and hold the **Ch2** - button again for approx. 5 seconds. To do this, the impulse intensity of **Ch1** and **Ch2** must be set to 00. Deletion of the Doctor's function is confirmed by a long acoustic signal.

8.7 Electric current parameters


Electrostimulation devices operate with the following electric current settings, which may affect the stimulation results differently, depending on the setting:

8.7.1 Impulse shape

This describes the time function of the electrical impulse.

It distinguishes between monophasic and biphasic pulse currents. In monophasic pulse currents, the current flows in one direction and in biphasic pulse currents the electrical impulse alternates its direction.

The Digital TENS/EMS device only provides biphasic pulse currents as these relieve muscles, cause little muscle fatigue and provide safer application.


8.7.2 Impulse frequency

The frequency indicates the number of individual impulses per second and is displayed in Hz (Hertz). It can be calculated by determining the cyclic value for the time period.


The relevant frequency determines which types of muscle fibres react best. Slow-reacting fibres react more easily to lower impulse frequencies of up to 15 Hz, whereas fast-reacting fibres only respond from approximately 35 Hz onwards.

Impulses of approx. 45–70 Hz are linked with constant tension in the muscles and quicker fatigue. Higher impulse frequencies are therefore better to use for high-speed strength and maximum power training.

8.7.3 Impulse width

This indicates the duration of an individual impulse in microseconds. The impulse width therefore determines, among other things, the penetration

depth of the electricity, where usually larger muscle mass requires a larger impulse width.

8.7.4 Impulse intensity


Setting the intensity level depends on the individual sensitivity of each user and is determined by a variety of variables, such as site of application, blood supply to the skin, skin thickness and the quality of the electrode contact. The settings should be effective but should never cause an unpleasant sensation, such as pain at the site of application. While a gentle tingling indicates sufficient stimulation energy levels, any setting that causes pain should be avoided.

When using the device for an extended period, you may need to adjust the intensity level, as your muscles may start to adapt to the impulse intensity.

8.7.5 Cycled impulse parameter variation

In many cases it is necessary to cover the overall tissue structure at the site of application by applying several impulse parameters. In the Digital TENS/EMS device, this is achieved by the provided programs, which automatically make a cyclical impulse parameter change. This also prevents individual muscle groups at the site of application being affected by fatigue.

The Digital TENS/EMS device provides sensible default current parameter settings. With this, you can change the impulse intensity at any time during use. In 6 programs you can also set various parameters for stimulation yourself.

9. Cleaning and maintenance

Adhesive electrodes

- To ensure that the adhesive electrodes remain adhesive for as long as possible, clean them carefully with a damp, lint-free cloth or clean the underside of the electrodes under lukewarm, running water and pat dry with a lint-free cloth.

i Before cleaning with water, remove the connection cables from the electrodes.

- Reapply the electrodes to the carrier foil following treatment.

Cleaning the device

- Remove the batteries from the device before cleaning.
- After use, clean the device with a soft, slightly damp cloth. If it is very dirty, you can also moisten the cloth with a mild soapy solution.
- Do not use any chemical or abrasive cleaning agents for cleaning.

! Ensure that no water enters the device.

Reuse of the device

Once it has been properly prepared, the device can be used again. Preparation includes replacement of the treatment electrodes as well as cleaning of the surface of the device using a cloth moistened with a mild soapy solution.

Storage

- Remove the batteries from the device if you will not be using it for a prolonged period of time. Leaking batteries may damage the device.
- Do not make sharp kinks in the connection cables and electrodes.
- Disconnect the connection cables from the electrodes.
- Reapply the electrodes to the carrier foil after use.
- Store the device in a cool, well-ventilated space.
- Never place any heavy objects on the device.

10. Accessories and replacement parts

You can obtain the following replacement parts directly from Customer Services:

<u>Designation</u>	<u>Item number and/or order number</u>
8 x adhesive electrodes (45 x 45 mm)	Item 661.22
4 x adhesive electrodes (50 x 100 mm)	Item 661.21

11. What if there are problems?

The device does not switch on when the ON/OFF button ① is pressed. How to proceed?

- (1) Ensure batteries have been inserted correctly and are in contact with the terminals.
- (2) Replace batteries, if necessary.
- (3) Contact Customer Services.

The electrodes do not stick to the body. How to proceed?

- (1) Clean the adhesive surface of the electrodes using a damp, lint-free cloth. The electrodes must be replaced if they still do not stick securely.
- (2) Clean the skin prior to any application; do not use skincare lotions or oils prior to treatment. Shaving may increase the life of electrodes.

There is no noticeable stimulation. How to proceed?

- (1) Press the ON/OFF button ① to interrupt the program. Check that the connection cables are correctly connected to the electrodes. Ensure that the electrodes are in firm contact with the treatment area.
- (2) Ensure the connection plug is firmly connected to the device.
- (3) Press the ON/OFF button ① to restart the program.
- (4) Check electrode positions and ensure that the adhesive electrodes do not overlap.
- (5) Gradually increase the impulse intensity.
- (6) The batteries are almost empty. Replace the batteries.

The battery symbol is shown. How to proceed?

Replace all batteries.

You are experiencing an unpleasant sensation at the site of the electrodes. How to proceed?

- (1) The electrodes are not positioned correctly. Check their positions and re-position, if necessary.
- (2) The electrodes are worn. This may cause irritated skin, as even distribution of the current across the entire area is no longer guaranteed. For this reason, the electrodes should be replaced.

Skin in the treatment area turns red. How to proceed?

Immediately stop treatment and wait until your skin has returned to its normal condition. If the redness is under the electrode and disappears quickly, there is no risk – this is caused by the locally stimulated, increased blood flow.

However, consult your doctor before you continue treatment if the skin irritation persists and if it is accompanied by an itchy sensation or inflammation. This may be caused by an allergic reaction to the adhesive surface.

12. Technical specifications

Name and model	SEM 43
Type	SEM 43
Output waveform	Biphasic rectangular impulses
Pulse length	50–450 μ s
Pulse frequency	1–150 Hz
Output voltage	Max. 100 V _{pp} (500 ohm)
Output current	Max. 200 mA _{pp} (500 ohm)
Voltage supply	3 x AAA batteries
Treatment time	Adjustable from 5 to 100 minutes
Intensity	Adjustable from 0 to 50
Operating conditions	5°C to 40°C (41°F – 104°F) at a relative humidity of 40–70%

Storage conditions	0 °C – 40 °C (32°F - 104°F) with a relative humidity of 10–90%
Transport conditions	-10°C – 50°C (14°F – 122°F) with a relative humidity of 10–90%
Dimensions	6.6 x 13.5 x 2.9 cm (including belt clip)
Weight	107 g (including belt clip, without batteries), 133 g (including belt clip and batteries)
Height limit for use	3,000 m
Maximum tolerated atmospheric pressure	700–1060 hPa
The serial number is located on the device or in the battery compartment.	

Note: If the device is not used according to the instructions specified, perfect functionality cannot be guaranteed!

We reserve the right to make technical changes to improve and develop the product.

This device conforms with the European standards EN 60601-1 and EN 60601-1-2 (in accordance with IEC 61000-4-2, IEC 61000-4-3, IEC 61000-4-4, IEC 61000-4-5, IEC 61000-4-6, IEC 61000-4-8 and IEC 61000-4-11) and is subject to particular precautions with regard to electromagnetic compatibility. Please note that portable and mobile HF communication systems may interfere with this device.

For more details, please contact our Customer Services at the address indicated.

This device meets the requirements of European Directive 93/42/EEC for medical devices, as well as those of the Medizinproduktegesetz (German Medical Devices Act).

13. Disposal

The empty, completely flat batteries should be disposed of through specially designated collection boxes, recycling points or electronics retailers. You are legally required to dispose of the batteries.

Note: The codes below are printed on batteries containing harmful substances:

Pb = battery contains lead

Cd = battery contains cadmium

Hg = battery contains mercury

For environmental reasons, do not dispose of the device in the household waste at the end of its useful life. Dispose of the device at a suitable collection point in your country. Dispose of the device in accordance with EC Directive – WEEE (Waste Electrical and Electronic Equipment). If you have any questions, please contact the local authorities responsible for waste disposal.

14. Warranty/service

Further information on the guarantee and guarantee conditions can be found in the guarantee leaflet supplied.

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