EN	
	
<div>MILEX[®] a CooperSurgical Brand</div>	
<h2>RING Folding Pessaries</h2>	
For uterine prolapse (Stage I and Stage II) R_xOnly	
AVAILABLE MODELS	
RING without SUPPORT	RING with SUPPORT
<div>[REF]</div> <div>MXPER (Pessary Only) MXKPER (Pessary, TRIMO-SAN™, Jel Jector)</div>	<div>[REF]</div> <div>MXPRS (Pessary Only) MXKPRS (Pessary, TRIMO-SAN™, Jel Jector)</div>
RING with KNOB without Support	RING with KNOB with SUPPORT
<div>[REF]</div> <div>MXPERK (Pessary Only) MXKPERK (Pessary, TRIMO-SAN™, Jel Jector)</div>	<div>[REF]</div> <div>MXPRSK (Pessary Only) MXKPRSK (Pessary, TRIMO-SAN™, Jel Jector)</div>
<div>[CE]0086</div>	
Single-Patient Use	
Wash pessary with mild soap and thoroughly rinse prior to initial use.	
DESCRIPTION	
The Milex® RING Pessary is a medical device made from silicone that is inserted into the vagina to function as a supportive structure of the uterus, bladder and/or rectum.	
WARNINGS	
<ul style="list-style-type: none">Sizes 11 through 13 contain a wire coil - remove before X-ray or MRI. Do not leave pessaries in place for long periods of time, as serious complications can occur which may require surgical intervention. Do not use these pessaries on a patient with a known silicone allergy.	
<div>[CooperSurgical]</div>	

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- Chemicals in various vaginal preparations can interact with the pessary material, resulting in discoloration or deterioration of the pessary. TRIMO-SAN™ does not interact with the pessary material. Use only TRIMO-SAN. Other materials have not been tested for compatibility.

INDICATIONS FOR USE

- The RING without SUPPORT pessary is indicated for support in Stage I and Stage II prolapse.
- The RING with SUPPORT pessary is indicated for support in Stage I and Stage II prolapse complicated by a mild cystocele.
- The RING with KNOB pessary, without SUPPORT is indicated for stress urinary incontinence complicated by Stage I or Stage II uterine prolapse.
- The RING with KNOB with SUPPORT pessary is indicated for stress urinary incontinence complicated by Stage I or Stage II uterine prolapse and mild cystocele.

CONTRAINDICATIONS

- The presence of pelvic infections or lacerations
- A noncompliant patient
- Endometriosis
- Pregnant patient

RECOMMENDED CLEANING INSTRUCTIONS*

- Prepare a cleaning solution by mixing a mild soap (such as Dawn® or equivalent) with tap water using the soap manufacturer's recommended concentration. Prepare this solution in a container large enough to fully submerge the device.
- Soak and Scrub
 - Soak the device in the container of prepared soap solution for a minimum of 5 minutes.
 - Following the 5-minute soak period, scrub the device for a minimum of 15 seconds with a soft-bristled brush, such as a tooth brush and/or pipe brush. Scrub device below water line to prevent aerosolization of contaminants.
 - Following scrub, inspect device for visible soil residue.
- Rinse
 - Remove the device from the soap solution and thoroughly rinse under flowing tap water for a minimum of 30 seconds.
 - Allow the device to dry.

*New pessaries are powdered with food-grade powder.

PATIENT INSTRUCTIONS FOR PESSARY USE

CAUTION

COOPERSURGICAL RECOMMENDS PESSARIES BE INSERTED AND REMOVED BY THE PHYSICIAN OR OTHER HEALTHCARE PROFESSIONAL UNLESS OTHERWISE DIRECTED.

- To ensure the desired correction of your condition, your healthcare professional needs your full cooperation.
- It is essential that your healthcare professional inspect your vagina at frequent intervals for evidence of pressure and/or allergic reaction.
- A gradual increase in the interval of inspection may be considered at the discretion of the treating practitioner.
- It is important to evaluate the effects of the vaginal gel TRIMO-SAN.

REPORT ANY OF THE FOLLOWING SYMPTOMS TO YOUR PHYSICIAN

- Any difficulty in urinating
- Any discomfort
- Any changes in the color or consistency of vaginal discharge
- Any increase in the amount of vaginal discharge or vaginal bleeding
- Any foul odor associated with vaginal discharge
- Vaginal itching
- If the pessary falls out

For medical emergencies and for all medically related advice, consult your healthcare professional.

FOR THE PHYSICIAN / HEALTHCARE PROFESSIONAL

Review these instructions with the patient to establish use regimen.

INSTRUCTIONS

- Wear dry gloves. When necessary, lubricate only the entering end of the pessary with TRIMO-SAN. Hold as illustrated in Figure 1.
- The pessary is folded along the axis of the bigger outer holes (or the inner notches for the RING without SUPPORT) by bringing the small round holes together. The KNOB will be at the top of the arch (see Figure 2). The arch formed points downward as shown with the RING pessary.

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		
Figure 1	Figure 2	Figure 3
3. Direct the pessary past the cervix into the posterior fornix. Allow the pessary to open again into the ring shape after passing the introitus.	4. The index finger is inserted deep into the vagina to turn the pessary approximately 90° so that the KNOB is resting behind the symphysis pubis (see Figure 3). The RING pessaries in this position cannot be folded and pushed out.	5. Ask the patient to sit, stand and bear down slightly. If there is no leakage, and the patient is comfortable with the pessary in position, have her empty her bladder. A properly fitted pessary takes up slack in redundant tissue, holding the uterus higher in the vagina.
6. If the patient can void without difficulty, and the pessary remains in position upon re-examination, and the patient is comfortable with the pessary in place, this is a good indication that the correct size may have been selected. Patient experience may vary.	Note: If the patient is unable to urinate with the pessary in position, remove it and fit her with the next smaller size. Repeat as necessary.	7. Examine the patient while she is in the standing position to ensure the pessary has not shifted position. The patient should not feel the pessary once it is in position. The pessary should not be too loose as it may turn or be expelled and it should not be too tight as it may cause discomfort.
8. The healthcare professional should be able to sweep one finger between the pessary and the vaginal walls. If there is not enough space to do this, the next smaller size should be tried. If excessive space exists, the pessary will not be effective and may rotate or even be expelled.	9. It is sometimes necessary to refit the patient with a different size or type of pessary after a period of time. Do not assume that a replacement will always be the same size as the previous one. Check the fitting to ensure continued patient comfort and relief of symptoms. The useful life of a pessary is limited. Examine frequently for signs of deterioration (such as cracks or breaks in silicone outer surface). A pessary should be replaced if damaged.	10. Ulcerations and erosions frequently occur in cases of complete prolapse due to irritation of the exteriorized cervix or vaginal wall. Whenever possible, reducing the mass and treating the irritation are primary steps before using a pessary. Prolapse reduction may resolve cervical vaginal irritation. Verification of cervical cytology (pap) and/or biopsy as clinically indicated.

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Dawn® علامة تجارية مسجلة للشركة Proctor & Gamble, سينتياي، أوهايو 45202 TRIMO-SAN™ هي علامة تجارية لشركة CooperSurgical, Inc. أما Milex® و CooperSurgical j وميلك هما علامتان تجاريتان مسجلتان لشركة CooperSurgical, Inc. © 2017 CooperSurgical, Inc.