Handleiding / Manuel / Manual

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Manual

Step 1: Preparation

1. Ensure that the pulse oximeter has sufficient power. If necessary, install or replace batteries according to the instructions in the manual.
2. Wash your hands and make sure your fingers are clean and dry before using the pulse oximeter.

Step 2: Placing the Pulse Oximeter

1. Open the pulse oximeter and place your finger in the designated compartment, usually located at the top of the device.
2. Ensure that your finger is fully covered by the sensor and that no external light can interfere with the measurement.

Step 3: Conducting the Measurement

1. Press the power button or the start button on the pulse oximeter to initiate the measurement.
2. Wait patiently while the device performs the measurement. This usually takes a few seconds.
3. Check the pulse oximeter's display to read the results.

Step 4: Reading the Results

1. The pulse oximeter typically displays two important measurements: oxygen saturation (SpO₂) and heart rate (beats per minute).
2. Oxygen saturation indicates the percentage of oxygen in your blood. A normal SpO₂ reading usually falls between 95% and 100%.
3. The heart rate is displayed in beats per minute (bpm).

Step 5: Notes and Monitoring

1. Take notes of measured values and times for future reference if necessary.
2. If using the pulse oximeter regularly, monitor changes in the measurements and consult a healthcare professional if you're concerned about unusual results.

Step 6: Turning Off

1. Turn off the pulse oximeter when you're done with the measurement to conserve battery power.