SILVERGEAR
MASSAGE GUN
RELAX AND RECOVER

User manual

EN

Thank you for purchasing the Silvergear Massage Gun. To let you experience this product to its full potential, please read all instructions carefully before using the product.

MASSAGE HEADS/FUNCTIONS

- A Ball head for large muscle groups like back, arms and thighs
- B U head for the vertebrae
- C Flat head for all body parts
- D Bullet head for body joints

SPECIFICATIONS

1. Power on/off button
2. Speed level LED indicator lights (red)
3. DC charging port
4. AC wall charger

- Level 1 - 2200 rpm
- Level 2 - 2400 rpm
- Level 3 - 2600 rpm
- Level 4 - 2800 rpm
- Level 5 - 3000 rpm
- Level 6 - 3200 rpm

CAUTION

- For adults only.
- This product can only be used on dry, clean skin surfaces of the body. It cannot be used through clothes.
- Only use on soft tissues of the body without causing pain or discomfort. In case of pain, stop immediately. Do not use it on the head or any hard or boned part of the body.
- Bruises may be caused when this device is used with too much pressure. Always keep an eye on the body part you're massaging.
- Keep fingers, hair and vulnerable body parts away from the shaft and back of the massage head, because they might be pinched.
- Do not place any object in the vent of the device.
- Do not immerse in water and do not let water into the vent of the device.
- Charge the device charger only with 24 VDC.
- Do not change or modify the device in any way.

CHARGING

Before first use, make sure you charge the Massage Gun fully for three hours. Only use the supplied AC/DC adaptor. When connected, the charging indicator light will glow red. This will turn green when fully charged.

After that, you can charge the battery at any time at any energy level. The average working time of the battery is two hours, depending on the speed and the applied pressure during use.

USING THE MASSAGE GUN

Unplug the charger before using the device. Attach the massage head of your choice to the massage gun, simply by pushing it in. Push the on/off button. The red LED light turns on as an indicator. Adjust the vibration level by pushing the on/off button. There are six different vibration levels. Press and move it gently on the skin. Do this for approximately 60 seconds per body part. To turn the device off, press the on/off button again until you passed all vibration levels.

MAINTENANCE AND CLEANING

Clean the device with a slightly wet towel and then wipe it dry with a soft cloth. Store it in a dry place or put it in the supplied carton when you travel with the device.

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