

USER INSTRUCTIONS

Portable home use Blackhead remover Device
Deep cleansing
Warm massager



Please read all instructions before operating.

Product Structure



Note: Before using the heating function please put the silicone sleeve on the bottom of the machine

Product Features

Our promise is to bring you a product that is affordable, authentic and result-oriented.

Introduction

The Portable At-Home Blackhead remover Device is a safe, non-invasive system that uses micro-crystalline technology combined with vacuum suction to deeply exfoliate into the epidermis of the skin. These rough particles help loosen the dead skin cells contained within the outermost layer of the epidermis - the stratum corneum. After they are loosened, the vacuum suction pulls the cells from the surface of the skin, revealing a fresh new layer of skin. This technique respects the skin's delicate layers and corrects skin conditions such as wrinkles, fine lines, enlarged pores, flaking skin, acne scars, pigmentation and dry skin. Suitable for all skin types.

Benefits

- Clears pores of impurities
- Removes dead skin cells
- Smooths skin
- Promotes collagen production
- Improves absorption of skin care products
- Improves blood circulation.

Results:

- Visibly clearer skin
- Smoother skin
- Increased elasticity
- Reduced large pores
- Youthful glow
- Reduced scarring
- Improved hydration

Function



Large Round Tip
 For blackheads
 Move around continuously



Small Round Tip
 For sensitive skin and blackheads



Medium Round Tip
 For wrinkles and fine lines



Body Round Tip
 In addition to the face, it can be used in other parts of the body. good for your health

Recommendations

- Never use this device on a daily basis.
- Allow the skin to recover after each treatment for a minimum of 5 days.
- To avoid bruising, do not use excessive pressure, do not stay in one spot and do not use for longer than 10 minutes per treatment.
- Face may be red for 1-24 hours. Avoid sun exposure for the next day.

Prepare Skin for Treatment

- Ensure skin is thoroughly cleansed to remove makeup, dirt or creams.
- Shower, bath, steam face or apply moist towel to face.

How to used

Smear blackhead derived liquid

Drop liquid onto the cotton pad and then apply it for three minutes

Hot application

Use hot compress function to apply the pores for 5minutes

Derived blackhead

Move the smart instrument slowly on your face to clean, then you are sparking without any blackhead.

Wash the face

Wash your face with water after sucking out blackheads

* Nursing face

apply shrink water or use a hydrating mask is better.

Directions

1. Charge the device.
2. Warm massager.
3. Attach desired tip.
4. Press on button and choose desired speed.
5. Move the tip on skin slowly.
6. Turn off the device after use and take off the tip, clean it by cotton swab or wash it.

Specification

| | |
|------------------|--|
| Product Size | 172x35x35mm |
| Working Voltage | 3.3V-4.2V |
| Standby Current | ≤30 μA |
| Battery Capacity | 500mAh |
| Working current | 40 mA |
| Suction Power | 55KPa,60KPa,65KPa |
| Net Weight | 125g |
| Charging Time | 1.5h |
| Working Time | 1.5h for blackhead remover 1h for heating |

WARNING:

1. Read all instructions before use to avoid injury.
2. Prohibit water from seeping into machine housing.
3. Do not use if tips are damaged.
4. Do not use a charger that has a damaged cord or plug.
5. Refrain from use when charging.
6. Always turn off device after use.
7. Discontinue use immediately if discomfort develops.
7. DO NOT use on eyelids or under eyes.