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Product

The Nuvance Codonein Spoon Set is specially designed to remove blackheads, pimples, sebaceous bumps and barley grains hygienically and effectively. It is a safe and professional way to treat blemishes without damaging your skin or causing scars – if used properly.

Total package:

- 4x blackhead removers
- Storage box
- Manual

Safety Tips

- Always clean your hands and tools before and after use.
- Never use the tools on open wounds or inflammations.
- Disinfect the tool with alcohol (min. 70%).
- Don't press too hard – this can cause scarring or damage.
- Cleanse and moisturize your skin afterwards.
- Treat your skin up to 1-2 times a week

How to use?

Preparation

1. Wash your face with a mild cleanser or facewash.
2. Steam your face over a bowl of hot water for 5-10 minutes or use a warm damp towel.
3. This opens up your pores and makes impurities easier to remove.

Use

1. Needle/lancet
 - Use only for closed comedones (white pimples).
 - Carefully pierce the top tip – not deep!
 - Then use a loop to press.
2. Ronde lus (klein of groot)
 - Place the ring over the blackhead.
 - Press down gently and gradually.
 - Wipe away the removed sebum with a clean cloth.
3. Angled loop
 - Ideal for nostrils, chin and edges of the face.
 - Use the same technique as with the straight loop, but at a slight angle.

Post-treatment

1. Cleanse your skin with a mild toner or micellar water.

2. Apply a soothing cream or gel
3. Avoid makeup or sun exposure for at least 4 hours after the treatment.

Maintenance

- Disinfect with alcohol before and after use (min. 70%).
- Store the tools dry in the supplied case.
- Clean the tools thoroughly with warm water and soap from time to time.

Notice

Mobius Lus

This product (or its packaging) is recyclable. Throw it away according to the local recycling regulations.

Triman

This product (or its packaging) is covered by separate collection. Follow the local recycling requirements for correct Waste management.

Green Dot

The Green Dot indicates that the manufacturer contributes financially to a system for the collection, sorting and recycling of Packaging waste. Check the local recycling guidelines for the correct Waste management.

European directives

CE marking on this product indicates that it complies with the relevant European Directives. Learn more about compliance and Technical documentation is available on request available.

Disclaimer

By using this manual, you consent to this disclaimer. TLS Commerce B.V. strives to provide complete, accurate and up-to-date information in the manual and on the accompanying web pages. If you have any questions or have incorrect information, please contact the customer service of TLS Commerce B.V. All information is intended for personal use. TLS Commerce B.V. is not responsible for damage caused by use, incompleteness or inaccuracy of the information in the event of careless use. The information is subject to change without prior notice.

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TLS Commerce B.V. | Nuvance