

MANUAL COLOURFUL LED BEAUTY MASK



Remote control, charging cable, connecting cable, fixed rope



- 1 Led lamp beads
- 2 Bandage connection lock
- 3 Medical grade pads
- 4 VPS material
- 5 The connect of power

- 1 Power indicator light
- 2 Time indicator
- 3 Switch
- 4 Turn on/pause
- 5 Adjust intensity
- 6 Adjust color
- 7 Adjust time



Power & Mask Connector



How to use:

Step1: Plug the USB connector at the end of the wire of your LED MASK.

Step2: Clean your face ,wait till it dries. Applying essence or skin care products。 This makes it easier for light to penetrate your skin.

Step3: Put on your Mask . Make sure that you wear it comfortably.

Step4: Turn on your Mask and you will see tinge of red light covering the whole mask.

Note: the light is at low-level intensity. You may close your eyes if you feel like a little unwell.

Step5: Keep your LED Mask close to your facial skin. You can also change different colors with the remoter.

Note: If you want the mask effect better, we suggest you use professional relational skin care product at the same time , the led light will promote the assimilate of skin care product, Enhances the efficacy of skin care product , then has better effect ,and it will not immediately see the effect ,it need some time to repair you skin,so keep use it and be patient.



470nm-LED Blue light:

Efficiently minimize and heal acne,repair the skin without leaving scars



520nm-LED Green light:

Neutralization,balance skin condition,relieve mental stress and effective dredge lymphoid and edema



590nm-LED Yellow light:

Improve coarse skin and wrinkle redness, fever, ringworm of the skin, increase immunity



630nm-LED Red light:

Whitening pale spot,tender skin and anti-wrinkle,repair damaged skin,smooth the fine wrinkles,shrink pores,hyperplasia of collagen.



LED Purple light:

It is red and blue dual-band light, a combination of two kinds of phototherapy effect,especially in healing and repairing the acne scar



LED Green Blue light:

Enhance the cell energy gradually,promote metabolism



LED Laser light (White light) :

Resolving age spots,improve fine lines and sagging skin.



Very important:

1. A single exposure can not exceed 30 minutes, one day can be used twice, the interval time is more than 6 hours.
2. The following patients are forbidden to use instruments: patients with epilepsy, bleeding tendency, scars, severe skin infections, and other skin diseases.
3. Facial plastic implants (hyaluronic acid and other implants) should not be used within 1 month.
4. Don't drink alcohol within 24 hours after use, do not eat spicy and spicy foods, do not eat seafood and other high-protein foods, avoid exposure (easy tanning)
5. This equipment is not waterproof. Please do not wash it with water and take waterproof measures.

6. Prohibit the use of crowd:

Pregnant women, children, Thyroid disease, Light allergic person, use of drugs can cause allergies.

