Back

105 mm

Cover

105 mm

148 mm

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105 mm

EN - USER INSTRUCTIONS

Your new personal scale has been designed to accurately indicate your weight gain or loss over a period of time and should give you many years of service under normal use.

CARE AND ATTENTION

- Placing your scales on a hard, flat surface will ensure the greatest accuracy and repeatability.
- Always weigh yourself on the same scales placed on the same floor surface, uneven floors may affect the readings.
- Your scales round up or down to the nearest 50g/0.1lb.
- Clean your scale with a damp cloth but do not allow water to get inside the scales. Do not use any chemical cleaning agents. Do not allow your scales to become saturated with water as this can damage the electronics.
Remove the batteries if the scales are not going to be used for a long period of time.
- Treat your scales with care. It is a precision instrument, do not drop or jump on it.

BATTERY CHANGE INDICATOR

When the batteries need to be replaced the following will appear on the screen:

Please open the back board of the scales, then take the batteries carefully out of the batteries compartment.

Do not bend the contacts back while you do this, or they may break. Dispose of the old batteries in a safe place. Do not place on fire.

CAUTION:

Slippery when wet.

148 mm

105 mm

OPERATING INSTRUCTIONS

Step 1.

Put the scales on a hard flat surface (avoiding carpets or any soft surfaces).

Step 2.

Stand on the scales, put your feet evenly on the scale without moving, later after about three seconds your correct weight will be displayed as below picture shown:

Step 3.

When you step off the scale, it will automatically switch off.

NOTE: When the scale Auto-off (the insulation slip has been removed), you just need to operate the scale according to above mentioned STEP 2 and STEP 3.

OVERLOAD INDICATOR

When the scales are overload (over 180kg/396lb/28.3 stone), the following will appear on the screen:

SPECIFICATIONS

- Equipped with 4 high precise strain gauge sensor system
- Tempered safety glass / platform
- Capacity: 180kg (396lb / 28.3 stone)
- Division: 50g (0.1lb)
- Foot-tap switch
- Auto zero resetting & auto power off
- Low power & over-load indicator
- Press "unit" button to convert unit kg/lb/st. The "unit" button is on the backside of the scale (under the battery box with "set/unit" index)

TECHNICAL DATA

Battery (2 x AAA)

Replace the batteries when the "low battery indicator" shows.

148 mm

105 mm

ENVIRONMENTALLY FRIENDLY WASTE PROCESSING

This marking means that this product should not be disposed of with other household waste. In order not to endanger the environment and human health and to promote the reuse of raw materials, this product must be disposed of responsibly. Please dispose used equipment at the appropriate collection points or at the shop where the product was purchased. They will accept this product for environmentally friendly disposal.

WARRANTY

- BCC B.V. gives a 2 year warranty on all defects that are the result of the hidden defects and that make the device unfit for normal use.
- The warranty start at the time of purchase; keep the proof of purchase carefully.
- Within the warranty period, any manufacturing and/or material defects will be rectified by us free of charge, either by repair, replacement of parts or exchange of the device.
- Repair will only be carried out under warranty if it is convincingly demonstrated (by means of your purchase receipt) that the day on which the complaint is submitted falls within the warranty period.
- The warranty expires if the defect is caused by damage as a result of an accident, incorrect use, neglect (e.g. not cleaning properly) or if interventions have been carried out or repairs have been carried out outside the service workshop of BCC B.V. (expect disassembly as stated in the user manual).
- The warranty is also not valid for connection to the wrong voltage, failure to follow the user manual and normal wear and tear of the device.
- BCC B.V. cannot be held liable for material damage or personal accidents as a result of connection in violation of the locally applicable safety regulations and technical standards (e.g. a faulty socket). The guarantee does not give any right to compensation under any circumstances.
- All other claims, including damage, are excluded unless otherwise required by law.
- All other warranty provisions mentioned above will not be used by us.
- Other than the above-mentioned warranty provisions will not be used by us.

SERVICE

For questions, please contact BCC B.V. customer service
Phone number: 020-3348888

What to do if the device does not work?

If your device shows defects, it is best to first contact your nearest BCC B.V. store. Our colleagues will ensure that these are resolved as soon as possible.

Repairs outside the warranty period

Repairs outside the warranty period are always possible.
There are of course costs associated with this.

148 mm