

Usage of the Red Light Therapy Device



Position your Red Light Therapy Device



Plug the Power in



Turn the Power On



Enjoy your Red Light Therapy

- Distance: Apply the light from 0-60 cm to the body.
- Time: Expose each area for 10-20 minutes, depending on distance, repeat on other areas if desired.
- Frequency: Do this 2 times to 14 times a week until results are felt.
- You can turn the Deep Red or NIR lights on or off at will, depending on the desired results.
- Skin rejuvenation, pain relief, hair growth, acne scars, inflammation, sports injuries, neck pain and stiffness, etc.
- Please ask your doctor if you have physiologic contraindications before treatment.

RED LIGHT THERAPY

Draw Power	60W
Input Voltage	DC24V/2.5A
LEDs Qty	40pcs
Wavelength	660/850nm
Product Size	247*147*40mm

Model: PHIRLT40



info@hipower.nl

www.hipower.nl

CE FC    **RoHS** DUTCH DESIGN