80mm

80mm

140mm

[www.optismile.eu](http://www.optismile.eu)

Not for resale

Overschiestraat 63  
1062 XD Amsterdam  
The Netherlands  
info@optismile.eu  
Made in China

Art. No.: 1975 / Batch no. 21B45

 **OptiSmile**  
*For professional Oral Care*

**Universal Oral Cleaner**

*For a healthy oral cavity and fresh breath every day*

EN • Manual NL • Handleiding FR • Mode d'emploi  
DE • Bedienungsanleitung PL • Instrukcja obsługi CZ • Návod k použití  
IT • Manuale SK • Priručka ES • Manual

Thank you for purchasing the OptiSmile Universal Oral Cleaner. This manual will provide instructions for effective and safe use of the device. Please read it carefully before use and store the manual in a safe place for future reference.

### 1. Package content

- Universal Oral Cleaner
- 5 Attachments
- Storage base
- Manual

### 2. Description of parts

1. Dental scaler
2. Tongue cleaner
3. Tartar remover
4. Gum massager
5. Tooth polisher
6. Handle
7. LED light
8. Storage base
9. Control button: on / off / switch mode (vibration & light)

### 3. Specifications

Model:	DW-6037
Material:	ABS / PC / TPE / SS 304
Power:	1 x AAA battery (not included)
Frequency:	15000rpm

### 4. Safety instructions

- For intended use only.
- The product is safe to use by children from the age of 14 and by persons with reduced physical or mental capabilities, or with lack of experience or knowledge, but only under supervision and with taking into account any hazards.
- Do not place, store or use the product near water.
- Do not immerse in water.
- When undergoing treatment for any oral care condition, consult your dentist prior to use.
- Stop using the product and contact your dentist, if you experience any discomfort or pain.
- Not recommended for pregnant women and patients with a pacemaker.
- This appliance is a personal care device and is not intended for use on multiple persons.
- Remove any oral jewellery prior to use of this product.
- Do not use the device if you have sores, blisters or cuts in your mouth.
- Do not allow your child to use the product if it has wobbly teeth, due to adult teeth coming through.

**Note:** When first using, you may feel the vibration is too strong. This is normal for an acoustic dental tool and you will get used to it after a few uses.

### 5. Operation

Press the control button to switch between these 3 modes:

1. LED light
2. LED light & Vibration
3. Vibration

### Function attachments

#### 1. Dental scaler

Remove any superficial stains on your teeth. Work the tip in a 90 degree angle and hold for 10 seconds, while gently moving back and forth.

#### 2. Tongue cleaner

Remove bacterial build-up for a healthier mouth and better breath. Wet and stretch out your tongue to brush it thoroughly. Move back and forth 5-8 times and repeat this 3-5 times for a better effect.

#### 3. Tartar remover

The cone-shaped head wipes off any plaque and removes dental calculus. It will give you clean and healthy teeth. Insert the tip between your teeth and hold it in its place for about 10 seconds.

#### 4. Gum massager

Improve or maintain healthy gums. Hold the dental head in its place for about 10 seconds. Gently move it around and repeat this until you have treated all your gums.

#### 5. Tooth polisher

The round, soft head polishes the tooth surface. It smoothens and whitens your teeth.

### 6. Maintenance

- Clean the device and all its parts after each use.
- Only use water or a mild non-abrasive detergent.
- After cleaning, dry thoroughly with a clean and soft cloth.
- Store the product in a cool and dry place, away from direct sunlight.

**Note:** You can clean the handle with thiomersal or ethyl alcohol liquid with more than 75%.

### 7. Disposal

The crossed out wheeled bin symbol on this product is to remind you that electrical products should not be disposed of with household waste. Please dispose of responsibly at an approved waste or recycling facility. Check with your local authority for recycling advice.

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