

# teesa®

# Sonic toothbrush

**TSA8005  
TSA8006**

**Návod k obsluze**

**CS**

**Bedienungsanleitung**

**DE**

**Owner's manual**

**EN**

**Εγχειρίδιο χρήσης**

**GR**

**Instrukcja obsługi**

**PL**

**Manual de utilizare**

**RO**

**Návod na použitie**

**SK**

## SPECIFICATION

3 work modes for children aged:

4-6: 38 000 strokes/minute

7-9: 41 000 strokes/minute

10 and above: 43 000 strokes/minute

IP rating: IPX7

2 min. auto-off, 30 sec. timer (each mode)

Soft bristles

Noise level: 65 dB

Power / charge indicator

Inductive charger

Usage on one charge: up to 30 days

Battery: 800 mAh

Power supply: 5 V; 1 A

In set: 2x brush head, charging unit



**English**  
**Correct Disposal of This**  
**Product**  
**(Waste Electrical & Electronic Equipment)**



EN

(Applicable in the European Union and other European countries with separate collection systems)

This marking shown on the product or its literature, indicates that it should not be disposed with other household wastes at the end of its working life. To prevent possible harm to the environment or human health from uncontrolled waste disposal, please separate this from other types of wastes and recycle it responsibly to promote the sustainable reuse of material resources. Household users should contact either the retailer where they purchased this product, or their local government office, for details of where and how they can take this item for environmentally safe recycling. Business users should contact their supplier and check the terms and conditions of the purchase contract. This product should not be mixed with other commercial wastes for disposal.

Made in China for LECHPOL ELECTRONICS Sp. z o.o.  
Sp.k., ul. Garwolińska 1, 08-400 Miętne.

teesa®

[www.teesa.pl](http://www.teesa.pl)