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## *Body scales*

PRODUKTNUMMER: 106472/106473/106474  
MODELL: DBPW001/DBPW002/DBPW003

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## INSTRUCTIONS

### General

ATTENTION! Observe the safety and assembly instructions in order to avoid the risk of injury or damaging the product.

**Only for domestic use and not for commercial use!**

Important: Read these instructions carefully and thoroughly. Keep these instructions to read through at a later date. If you give the product to anyone else at some point in the future, please ensure you also pass on this manual.

Keep small parts out of the reach of children. Pay attention especially when unpacking it. Keep plastic bags and other packaging away from children reach. Danger of Suffocation!

Check the delivery for completeness. Later complaints cannot be accepted.

Check all elements and parts for damages. Despite careful checks since may happen that even the best product takes damage during the transport. In this case, do not build up of your product. Defective parts can endanger your health.

Never modify the product! Changes will void the warranty and the product may be unsafe or even dangerous.

### Intended Use

The scales is designed for the weighing of people in the home. It is not suitable for use in clinics, medical practices, etc.. The maximum carrying capacity of 180 kg must not be exceeded.

Use the product only for its designed purpose of use. Manufacturer will not assume any responsibility in case of damages cause of an unappropriated use. Any modification applied by customer on the original product will have a negative impact on its safety, cause danger and warranty will be cancelled.

## Safety instructions

This device is not intended for use by persons (adults and children) with physical, sensory or mental impairments or lacking experience and/or knowledge of use unless they are supervised by a person responsible for their safety or have received instructions from that person on how to use the applicator.

Children should be supervised to make sure they do not play with the device.

Do not submerge the scales in water or other liquids.

Place the scales on a stable, even surface.

Do not overload the scales. Pay attention to the "Technical Data".

Protect the scales from sharp jolts and vibrations.

Do not place the scales in the vicinity of devices that emit electromagnetic radiation, such as mobile phones, wireless devices, and the like. This can lead to incorrect display indicators or malfunctions.

Do not place any objects on the scales when you are – for example – storing them in a cupboard. The sensitive measuring electronics in the scales could be damaged by heavy objects.

Do not attempt to repair the device yourself in the event of a malfunction. Repairs should only be carried out by authorized service centers.

Do not use the appliance in a location where extreme temperatures or humidity may occur. Keep the scales away from liquids.

## Risk of injury

1. Do not step onto the outer edge of any one side of the scales. Risk of tipping!
2. Do not step onto the scales with wet feet or when its surface is damp. Risk of slipping!
3. The weighing surface is made of glass. If you fall down, there is a risk of the glass breaking and thus a risk of injury from glass splinters.

## BATTERY SAFETY INSTRUCTIONS

- Do not disassemble batteries!
- Remove discharged batteries from the device immediately!
- Increased risk of leakage, avoid contact with skin, eyes and mucous membranes! If battery acid comes in contact with any of these parts, rinse the affected area with copious amounts of fresh water and seek medical attention immediately!
- If a battery has been swallowed seek medical attention immediately!
- Insert the batteries correctly, observing the polarity!
- Always keep the battery compartment well closed!
- Remove the batteries from the device if it is not going to be used for an extended period!
- Keep batteries out of children's reach!
- Do not attempt to recharge these batteries! There is a danger of explosion!
- Do not short circuit! There is a danger of explosion!
- Do not throw into a fire! There is a danger of explosion!
- Keep unused batteries in their packaging away from metal objects in order to prevent short circuiting!

## SPECIFICATIONS

### Model: DBKW001/ DBKW002/ DBKW003

- Maximum weight capacity: 180kg
- Minimum measurable weight: 2kg
- Graduation: 0.1kg
- User selectable units: kg/lb/st
- Low-battery indicator: LO
- Overload indicator: EEEE
- Operating temperature 5°C ~ 40°C

## OPERATION

1. Get one CR2032 Lithium put in the battery seat of the back of scale
2. Choose the weight unit and move to the right position.
3. The scale must put on the smooth hard floor.
4. The scales automatically switch on when you step onto the weighing surface with a little force.
5. Always stand in the center of the platform, thus enabling the scale to take your weight exactly.
6. Within seconds, personal weight is established electronically and displayed
7. Once you step off the scale after weighing, it will switch itself off automatically after few seconds.

IMPORTANT: REMOVE THE PLASTIC UNDER THE BATTERY

## ERROR DISPLAY

- When the screen shows "LO", it means that the power is low, please change the batteries as soon as possible.
- When the scales are overloaded, "EEEE" will appear on the display.

## ADVICE FOR CARE

If the scale surface is dirty, please clean the scale with soft cloth, do not use any chemical detergent and allow water to get inside the scale.

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## Disposal:

Once item end-of-life, please dispose valuable raw materials for a proper recycling. If you are not sure how to handle it in a correct way please contact your local disposal or recycling centers for advice.

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This marking indicates that this product should not be disposed with other household wastes throughout the EU. To prevent possible harm to the environment or human health from uncontrolled waste disposal, recycle it responsibly to promote the sustainable reuse of material resources. To return your used device, please use the return and collection systems or contact the retailer where the product was purchased. They can take this product for environmental safe recycling.

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Copyright by
Deuba GmbH & Co. KG