# MEDISANA<sup>®</sup>

home of wellness

- DE** Shiatsu-Massage-Sitzauflage <sup>MCS</sup>
- GB** Shiatsu massage seat cover <sup>MCS</sup>
- FR** Matelas de massage Shiatsu pour siège <sup>MCS</sup>
- IT** Cuscino a sedile per massaggio shiatsu <sup>MCS</sup>
- ES** Respaldo de masaje de Shiatsu <sup>MCS</sup>
- PT** Coxim com efeito de massagem Shiatsu <sup>MCS</sup>
- NL** Shiatsu-massagekussen <sup>MCS</sup>
- FI** Shiatsu-hierontaistuinalusta <sup>MCS</sup>
- SE** Shiatsu massagesits <sup>MCS</sup>
- GR** Μαξιλάρι για μασάζ Σιάτσου <sup>MCS</sup>

Art. 88932

**Gebrauchsanweisung**  
**Instruction Manual**  
**Mode d'emploi**  
**Istruzioni per l'uso**  
**Instrucciones de manejo**

Bitte sorgfältig lesen!  
Please read carefully!  
A lire attentivement s.v.p.!  
Da leggere con attenzione!  
¡Por favor lea con cuidado!

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## GB

- ① 4 three-dimensional vibrating twin massage heads (Shiatsu massage)
- ② Fastening straps (rear)
- ③ Plug with cable
- ④ 2 Massage motors (Vibration massage)
- ⑤ Remote control pocket
- ⑥ Control unit
- ⑦ On/Off button with LED (switches the appliance on/off)
- ⑧ DEMO button with LED (performs all massage functions)
- ⑨ HEAT button with LED (switches the heating on/off)
- ⑩ Shiatsu massage button with 2 LEDs (3 massage zones to choose from: whole back, upper back or lower back)
- ⑪ Vibration massage button with 4 LEDs (on/off, 3 intensity levels to choose from: L low, M medium or H high)
- ⑫ Timer button with 3 LEDs (3 times to choose from: 5/10/15 minutes)

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## IMPORTANT INFORMATION! RETAIN FOR FUTURE USE!

Read the instruction manual carefully before using this device, especially the safety instructions, and keep the instruction manual for future use. Should you give this device to another person, it is vital that you also pass on these instructions for use.

### *Explanation of symbols*

This instruction manual belongs to this device. It contains important information about starting up and operation. Read the instruction manual thoroughly. Non-observance of these instructions can result in serious injury or damage to the device.

#### **WARNING**

These warning notes must be observed to prevent any injury to the user.

#### **CAUTION**

These notes must be observed to prevent any damage to the device.

#### **NOTE**

These notes give you useful additional information on the installation or operation.

Protection category II

Lot number

Manufacturer

**before using the unit**

- Always carefully inspect the cable, the control unit and the Shiatsu massage seat cover for damage before use. Do not attempt to use a defective unit.
- Do not use the unit if there are visible signs of damage on the unit or the cable components, the unit does not work properly or the seat or the control unit have been dropped or become damp. In order to prevent hazards, always send the unit in to the service centre.
- The zip fastener on the massage device must not be opened. It is only fitted for reasons relating to the manufacturing of the device.

**operating the unit**

- The Shiatsu massage seat cover must only be used for its intended purpose as described in the instruction manual.
- Using the Shiatsu massage seat cover for any other purpose invalidates the warranty.
- Never leave the unit unattended if it is connected to the power supply.
- The maximum operating time for one session is 15 minutes. Allow it to cool off for at least 15 minutes before using it again.
- Do not allow the Shiatsu massage seat cover to come into contact with pointed or sharp objects.
- Never put the unit down next to an electric stove or other sources of heat.
- Never use the Shiatsu massage seat cover if it is folded up or crumpled.
- Never use the device under covers, such as blankets or pillows.
- Operate the device in an upright position only, never place it flat on the floor.
- Make sure that all moving parts can move freely at all times.
- Never sit or stand on the moving parts of the device with all your weight and do not place objects on the device.

**The device is only intended for household use and must not be used for commercial or medical purposes.  
If you have health concerns, consult your doctor before using the Shiatsu massage seat cover.**

**The massage seat cover must only be used in enclosed spaces!**

**Do not use the massage seat cover in damp rooms (e.g. when bathing or showering).**

**maintenance and cleaning**

- The unit is maintenance free.
- You may only clean the device yourself. If the device malfunctions please do not repair it yourself as this leads to the automatic expiry of all guarantee rights. Ask your specialist dealer and have repairs carried out by the service.
- Do not use the unit if there are visible signs of damage on the unit or the cable components, the unit does not work properly or the massage seat cover or the control unit have been dropped or become damp. In order to prevent hazards, always send the unit in to the service centre.

# 1 Safety instructions / 2 Useful Information

**GB**

- Do not immerse the device in water or other fluids.
- If, nonetheless, liquid should penetrate the device, unplug the mains cable immediately.

**Do not wash!**

**Do not dry clean!**

## 2 Useful Information

### Thank you very much

Thank you for your confidence in us and congratulations! You have acquired a **MEDISANA** quality product with this Shiatsu massage seat cover **MCS**. This unit is designed for massaging the back and upper thigh areas.

We recommend that you carefully read the following information on use and care of the product so that you experience its benefits and enjoy using your **MEDISANA** Shiatsu massage seat cover **MCS** for many years.

### 2.1 Items supplied and packaging

Please check first of all that the device is complete and is not damaged in any way. In case of doubt, do not use it and contact your dealer or your service centre.

The following parts are included:

- 1 **MEDISANA** Shiatsu massage seat cover **MCS**
- 1 Control unit
- 1 Instruction manual

The packaging can be reused or recycled. Please dispose properly of any packaging material no longer required.

If you notice any transport damage during unpacking, please contact your dealer without delay.

### WARNING

**Please ensure that the polythene packing is kept away from the reach of children! Risk of suffocation!**

## GB 2 Useful Information / 3 Operating

MEDISANA<sup>®</sup>  
home of wellness

### 2.2 What is Shiatsu massage?

Shiatsu is a form of finger pressure massage and one of the most important massaging techniques along with manual lymphatic drainage and reflexology. This form of physical therapy was developed in Japan but stems from traditional Chinese massage. Gentle stroking and beneficial pressure application techniques are used to awaken life forces and encourage them to flow through the body. The object of the massage is not just to ease tension, but also to optimise the flow of life force and harmonise body and soul.

### 2.3 How does the MEDISANA Shiatsu massage seat cover work?

In the **MEDISANA** Shiatsu massage seat cover **MCS** you have a modern device for genuine Shiatsu style back massages. The **MEDISANA** Shiatsu massage seat cover **MCS** has four three-dimensional rotating twin massage heads for Shiatsu massage of the three zones "Upper Back", "Lower Back" and "Full Back". The device is equipped with an additional 3-level vibration function and a heat function in the massage heads. The unique combination of intensive shiatsu and vibration massage with radiated heat makes the massage seat cover an effective and relaxing treatment. Their settings can be comfortably adjusted by means of the easy-to-handle control unit **6**.

## 3 Operating

### 3.1 Usage

The use of the shiatsu massage seat cover will improve your wellbeing after a stressful day. The massage effect is particularly useful for treating tense muscles and tired tissue. The unit provides pleasant relaxation for overcoming everyday stress. The massaging can improve your well-being, and is also useful in the sport and fitness areas.

### 3.2 Operation

#### Start-up

- Place the massage seat cover on a high chair or armchair and attach it to this by means of the fasteners **2** on the back. Ensure that the seat cover is supported over its entire area and has a sufficiently high backrest.
- Insert plug **3** into the mains socket. Ensure that the power supply is easily accessible.
- Take a seat on the shiatsu massage seat cover and familiarise yourself with the functions.
- Sit comfortably and upright. Ensure you sit centered and such that the massage heads **1** can move to the left and right of your spinal column.

#### Operating the unit

The **MCS** shiatsu seat cover has three functions for an effective massage: fingertip massage, vibration massage and heat. The heat function can only be used together with the shiatsu massage, during which the vibration massage can also be activated. Do not use the unit for more than 15 minutes at a time. After approx. 15 minutes, the unit will switch off automatically. Always allow the unit to cool down fully before you use it again.

### Buttons on the control unit

Use the buttons on control pad **6** to operate the massage seat cover:

- Switch on the device by pressing On/Off button **7**. The LED illuminates. The device starts automatically with the 3D massage on the full back setting. The LED shows the current massage region.
- Select the region you wish to massage with the Shiatsu massage button **10**.
  1. *Entire back*  
The Shiatsu massage heads **1** move from the "Park" position at the base of the backrest, up to the shoulder area, massaging the back in a circular motion whilst doing so. Once the uppermost position has been reached, the direction of rotation of the massaging heads is reversed. The heads then travel down to the base of the backrest. Once the heads reach this point, the direction of rotation is reversed again and then travel up the back again.
  2. *Lower back*  
The Shiatsu massage heads **1** travel up and down the lower half of the back in a circular motion. The reversal of the direction of rotation takes place when the respective end positions are reached.
  3. *Upper back*  
The Shiatsu massage heads **1** travel up and down the upper half of the back in a circular motion. The reversal of the direction of rotation takes place when the respective end positions are reached.
  4. *Stop*
- You can enhance the shiatsu massage at any time with a vibration massage for your thighs by pressing the vibration massage button **11**.
- The vibration massage in the upper thigh area can be operated separately by pressing vibration massage button **11**. The LED next to "L" illuminates and both vibration massage motors **4** are switched on. They begin with a gentle massage. Press button **11** again and the vibration level increases (medium intensity massage). The LED next to "M" illuminates. Press button **11** a third time and the vibration level increases (highest intensity massage). The LED next to "H" illuminates. After button **11** is pressed a fourth time, the vibration massage is switched off and the LED extinguishes.
- The heat function can be activated at any time with the shiatsu massage, during which the vibration massage can also be activated. Press the heat button **9** to activate the function; the LED on the control unit and the massage heads light up red. The heat is radiated through the massage heads **1**. Press the button **9** again to de-activate the heat function; the LED and the red light of the massage heads go out.
- *Timer function:*  
The automatic switch-off function is activated with selection of the massage region and is preset to 15 minutes. You can set the automatic switch-off function to 5, 10 or 15 minutes by pressing the TIMER button **12**.

Please note that changing the timer setting makes the timer start to count again from the beginning.

To avoid overheating the device, select the maximum massage duration at the beginning of the massage and do not change it again during the massage.

Do not use the massage device for longer than 15 minutes and then allow it to cool. Using the device for longer can overstimulate the muscles and cause tension instead of relaxation.

The massage can be finished at any time by pressing the On/off **7** button.

- *Demo function:*  
Press the DEMO button **8**, the device starts automatically with the 3D massage on the different back setting with different functions. The LEDs show the current massage regions and functions. The device will stop once all functions worked.
- *To switch off the device:*  
Press the On/Off button **7**. The power LED begins to flash and the massage heads travel back down to the "Park" position. When this position is reached, the LED extinguishes. If you press the On/Off button **7** when one of the massage functions is switched on, the vibration motors just stop immediately. First the massage heads travel downwards to the "Park" position. The LED flashes meanwhile and then extinguishes.
- After each use, switch the device off using the On/Off button **7** and remove the plug **3** from the mains socket. Always store control unit **6** in its bag **5**.

#### 4.1 Cleaning and maintenance

- Before cleaning the device, make sure that it is turned off and the mains plug is removed from the socket. Allow the unit to cool down.
- Only clean the Shiatsu massage seat cover with a soft cloth slightly dampened with a mild detergent (e.g. soap suds). Under no circumstances should you use aggressive detergents, strong brushes, solvents, petrol, thinners or alcohol. Rub the seat dry with a clean dry cloth.
- When cleaning, never immerse the unit in water and make sure that no water gets into the unit.
- Wait until the device is completely dry before using it again.
- Undo any kinks in the cable.
- You should ideally store the device in its original packaging and keep it in a clean, dry place.

#### 4.2 Disposal

This product must not be disposed of together with domestic waste. All users are obliged to hand in all electrical or electronic devices, regardless of whether or not they contain toxic substances, at a municipal or commercial collection point so that they can be disposed of in an environmentally acceptable manner. Consult your municipal authority or your dealer for information about disposal.

#### 4.3 Technical Data

Name and model	: <b>MEDISANA</b> Shiatsu massage seat cover <b>MCS</b>
Voltage supply	: 220 - 240 V~ 50 - 60 Hz
Power consumption	: approx. 60 W
Operating time	: max. 15 minutes
Timer	: after approx. 5, 10 and 15 minutes
Operating conditions	: only use in dry rooms
Storage conditions	: in a cool, dry place
Dimensions (L x W x H)	: approx. 105.5 x 42 x 9.5 cm
Weight	: approx. 5.6 kg
Power cable length	: approx. 2.1 m
Article No.	: 88932
EAN Code	: 40 15588 88932 5

**In accordance with our policy of continual product improvement, we reserve the right to make technical and optical changes without notice.**

**GB 5 Warranty****MEDISANA**  
home of life**Warranty and  
repair terms**

Please contact your dealer or the service centre in case of a claim under the warranty. If you have to return the unit, please enclose a copy of your receipt and state what the defect is.

The following warranty terms apply:

1. The warranty period for **MEDISANA** products is three years from date of purchase. In case of a warranty claim, the date of purchase has to be proven by means of the sales receipt or invoice.
2. Defects in material or workmanship will be removed free of charge within the warranty period.
3. Repairs under warranty do not extend the warranty period either for the unit or for the replacement parts.
4. The following is excluded under the warranty:
  - a. All damage which has arisen due to improper treatment, e.g. non-observance of the user instructions.
  - b. All damage which is due to repairs or tampering by the customer or unauthorised third parties.
  - c. Damage which has arisen during transport from the manufacturer to the consumer or during transport to the service centre.
  - d. Accessories which are subject to normal wear and tear.
5. Liability for direct or indirect consequential losses caused by the unit are excluded even if the damage to the unit is accepted as a warranty claim.

 **MEDISANA AG**  
Jagenbergstraße 19  
41468 NEUSS  
GERMANY  
eMail: [info@medisana.de](mailto:info@medisana.de)  
Internet: [www.medisana.de](http://www.medisana.de)

The service centre address is shown on the attached leaflet.

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