

2-IN-1 Skin Tag Removal Kit

Skin Tag Removal Kit works by stopping the skin tag's blood supply and results can be seen within just a few days of applying a band. The bands are suitable for small to medium sized skin tags. The Skin Tag Removal Kit designed to be used on skin tags on the face and body but we do not recommend using it around the sensitive eye area.

Skin Tag Removal Kit Contents :

2 × Band applicator cone	1 × Band remover
40 × Removal bands	10 × Cleansing swabs
36 × Skin repair patches	1 × Instruction manual

Instrucitons for use :

- 1) Clean the skin tag and the surrounding area with one of the provided cleansing swabs.
- 2) Push one of the bands up and onto the top section of the cone.
- 3) Place the cone over the skin tag and level with the surrounding skin area.
- 4) Place the band remover onto the cone and push lightly until the band pops off the end.
- 5) The band will now be held securely in place where it wil remain until the skin tag has dropped off.

Cleansing swab ingredients :

Swabs are satrurated with isopropyl alcohol 70%.

