## duct features

Strong environmental temperature adaptability

New control probe structure to ensure more accurate measurement

Two measurement modes of body temperature and object temperature, one bottom mode switching operation is convenient and practical;

Fever alarm prompt, backlight colour prompt;

Automatically save 32 sets of temperature measurement values;

Large size LCD screen high brightness three colour backlight, Clear and soft display;

You can select Two temperatures units: Celcius and Fahrenheit

Automatic shutdown saves power and energy;

Simple operation mode.

## cations

se read this manual carefully before use. It is forbidden to immerse the head of the temperature gun in any liquid, and it is not permitted to use it in an excessively high or low-temperature environment for an extended period. It is not allowed to collide, fall, mix with sharp objects, and dismantle by yourself. Do not use in a robust electromagnetic interference environment. Put the forehead temperature out of the reach of children. The measurement results can't replace the doctor's diagnosis.

## duct performance, major structural components

Measuring distance:	3cm~ 5cm	
Measuring range:	For the human body, 34°C~42.9°C (93.2°F~109.2°F)	For object surface 0°C~100°C (32°F~212°F)
Measuring accuracy	Object surface temperature +/- 1°C/1.8°F	Human body temperature measurement within (35-42)      above (35-42) +/- 0.2°C/0.4°F      +/-0.3°C/0.54 °F
Operating ambient temperature:	0°C~42°C (32°F~107.6°F)	
Power supply:	DC1.5V * 2 (2PC AAA BATTERY)	
Battery status:	Low battery warning	
Memory function:	Memorize 32 measurements	
Display unit:	Celcius and Fahrenheit	
Automatic shutdown:	Automatic shutdown with no action in 30 seconds	
Product measurement:	104g 149*77*43 mm	

Infrared detector  
 LCD  
 Memory switch  
 °C/°F  
 Set  
 Memory switch  
 Battery

## play icons

Battery indication		Display	Power is low. Please replace the battery
		Display and flashing	The battery is too low to work replace immediately
		No Display	Fully charged
Sound		Display	Buzzer sound on
		No display	Buzzer sound off
Temperature measurement mode		Body temperature mode	
		Object surface temperature mode	
Display unit	°C °F	Celsius temperature Fahrenheit temperature	
Display value		Measured temperature information	
Memory	M	Display value of memory value and group number	

## sets of body temperature memory data queries

When turning on the device, press the "memories" button, and the LCD will display the first memory value; press again to display the next record, a total of 32 groups of memory data.

## Function Key Description

Function key	Function description
SET	Short press can switch between human temperature and object temperature modus, long press set in order to enter setting menu
°C/°F	Switch between Celsius and Fahrenheit
Memories	Review previous measured temperatures

product provides prompt tone switch, temperature compensation function  
ng, the comparison table of setting menu is as follows:

press the set key to enter the setting mode. F1 is the prompt tone switch, and  
rompt tone is turned on or off by the C/F key. Press the set key again, F2 is the  
erature compensation, increase the temperature compensation through the C/F  
and reduce the temperature compensation with the memory query key. Press  
set key again to exit and save the setting mode.

## **ety instructions**

r the battery is installed, press the power on button for 1 second, after the  
ce is turned on. Confirm that the measurement mode above the screen is body  
erature mode (non object temperature mode). Aim the forehead temperature  
probe at the middle of the forehead and  
it vertical. The distance is about 3-5cm. Press the measurement trigger to start  
measurement. After about 1 second, the screen will display the measured value.  
n the device is not used, it will shut down automatically after 30 seconds.

Before and after use, please keep the sensor and probe cavity clean.  
Please place the forehead temperature gun in a temperature stable  
environment for use. When ambient temperature changes greatly (e.g. from  
indoor to outdoor), please place for about thirty minutes after measurement.  
Do not begin to measure body temperature immediately after measuring the  
temperature of an object at very high or very low temperature. Please place  
it for ten minutes.

When the measurement object comes from a place that is quite different  
from the measured ambient temperature, it should stay in the test  
environment for at least five minutes.

Try not to measure in the forehead blowing, drenching water, sweating,  
coating, smear cosmetics, etc. Do not take body temperature within 30  
minutes after exercise, bath, meal. temperature. When the working  
environment temperature is lower than 0 C the screen will display Err 1,  
backlight red, when the ambient temperature is higher than 42 C, the screen  
will display Err 2, backlight red;

## **age conditions, methods**

ucts must be kept clean and kept in a dry place. Don't put the electronic  
read temperature gun in place with electric shock. Do not store the electronic  
read temperature gun in an extreme temperature environment above 55C or  
w -10C and humidity above 85%. The transportation process should prevent  
heavy fall and so on, can not be transported in the extreme temperature  
ronment above 55°C or below -10 C, humidity above 85%.

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