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## Flowee spijkermat & kussen

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# Directions for the Nail Mat & Pillow

The trendy nail mat Flowee is a modern version of the bed of nails from India. Your blood circulation and nerve endings are stimulated due to the 'sharp points' on the mat. You also produce the 'feelgood' hormones endorphin and oxytocin. This results in a feeling of relaxation, pain relief and euphoria. Flowee creates a good flow of energy. This is the foundation of our good health: your body and spirit remain more in balance.

Use the Flowee to relax a stiff neck and shoulder muscles or to ease headache and/or muscular pain. The Flowee is also ideal for inner peace. You can even combine it with Yoga and Mindfulness. Have you got tired or cold feet? Then the mat can be very invigorating. Use the Flowee whenever it suits you.

## Usage

The 3 most common positions are the back, the neck and the feet. These are simple to apply and also have the best effect. You can use the Flowee on practically every part of your body. You can discover this for yourself once you have more experience. The mat can be used at whatever time of the day suits you!

### Use for the back

The lying back position is the most widely used position on the Flowee. There are numerous nerves situated in the back and these influence our entire body. This lying down position is also the most effective for relaxation. Lay the mat on your bed, sofa or floor. It is also possible to use a pillow for under your head to lie more comfortably. Lie down carefully on your back on the mat and spread your weight evenly. You will feel some pain during the first few minutes but that will quickly disappear. Breathe calmly and relax.

If the mat feels too uncomfortable during the initial stages, then you can wear a thin t-shirt or place a thin sheet on the mat. Even if your pain threshold is low you can still use the Flowee with ease. It is important to listen to the signals your body is sending. Your body gradually adapts and in between layer is no longer necessary.

Most users experience a feeling of peace and relaxation during their first sessions. After a few minutes you experience a warm feeling or tingling, this is the result of heightened blood circulation. Increase the use of the Flowee gradually. For the first sessions 10 to 15 minutes, increasing it as you see fit. 15 to 20 minutes is sufficient for the effect, but you may use it longer. Falling asleep on the mat due to the deep relaxation is not uncommon. You are able to sleep more deeply due to the Flowee so you can get a better night's sleep and wake up more refreshed.

### Use for the feet

Tired feet bothering you? Place the mat on the floor, sit down on a chair or sofa and put your feet on the mat (with or without socks). Press your feet gently on the sharp points and experience the invigorating tingles from the Flowee. You can decide yourself what feels nice and pleasant. You can keep your feet still, or move them up and down for a massage effect. If you like you can even stand up on it.

Due to stimulation of the nerve endings and the improved blood circulation the tiredness quickly disappears from your feet.

### Use for the neck and shoulders

Headache or a stiff neck or shoulder muscles? Use the Flowee pillow and lay your head down on it for 15 minutes. Lying with your head down on the mat itself can also feel nice and relaxed. Tip! For people with little hair a thin layer of material in between can feel more comfortable.

The improved circulation makes the muscles relax, besides which the endorphin that is produced has a pain killing effect.

### **Tips to extend the life of the Flowee:**

- You can clean the nail mat by hand washing it at a maximum of 30°; (remove the filling)
- Do not use in the tumble dryer;
- Keep it out of direct sunlight and store in a dry place;
- Do not place any heavy objects on the mat;
- Air the mat regularly.

### **Important!**

Use of the Flowee is **not** recommended in the following cases:

- Problems with blood clotting, heart trouble and/or in combination with the use of blood thinning medication;
- On problem skin areas where the Flowee comes in contact. Avoid contact with, for example large birthmarks and/or skin infections. Be careful while turning because of the sharp points;
- During pregnancy: consult with your doctor/midwife before using the Flowee.

If in doubt, always consult a doctor before using the nail mat.

### **Extra tips for using the Flowee:**

- If you have little baggage space and want to take your Flowee on holiday with you: remove the filling from the mat.
- Use an alarm if you have to keep an appointment shortly after you use the Flowee. There is a chance that you will fall asleep.

**Wishing you much relaxation and vitality with the Flowee!**