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# HANDLEIDING

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**MELVY.**

# HANDBOOK

Congratulations on your brand new Aceso Infrared Lamp. This powerful device is named after the goddess of the healing process and delivers over 140mw/cm<sup>2</sup> of light to the body in two different wavelengths: 660nm and 850nm. After careful study, it has been clinically proven that these wavelengths have a positive effect on the human body.

Read this manual carefully to learn everything you need to know to safely and effectively enjoy the benefits of infrared therapy.

## PACKAGE CONTENTS

Infrared lamp Aceso	1x
Protective goggles	1x
Power cord	1x

## WHAT IS INFRARED THERAPY

Infrared therapy uses specific wavelengths of natural light to improve the health and function of cells without harmful UV rays, excess heat, or negative effects. Light therapy is one of the earliest recorded forms of medicine. Sun therapy was first used by the Egyptians, and forms of light therapy were also applied by the ancient Greeks, Chinese, and other cultures.

There is no doubt that light has biological effects. In fact, your body needs light to become healthier. The visible red light of 660nm promotes skin health and collagen production. The invisible near-infrared light of 850nm penetrates deeper into tissue, organs, and joints and promotes muscle recovery and reduces joint pain.

## BENEFITS OF ACESO

- ✓ Reduces stress.
- ✓ Reduces fatigue.
- ✓ Helps with depression.
- ✓ Reduces muscle and joint pain.
- ✓ Stimulates blood circulation.
- ✓ Increases collagen production and metabolism.
- ✓ Improves your immune system.
- ✓ Stimulates better sleep.
- ✓ Improves tissue repair.
- ✓ Reduces scarring.
- ✓ Durable, strong LED lighting (660nm, 850nm).
- ✓ Light and compact size.
- ✓ Free protective glasses.
- ✓ Available in different colours.

## DISPLAY FUNCTIONS

The Aceso infrared lamp features a display with four functions. More info for each function are given below.

### PROGRAM

Choose from first mode with R/NIR+ (660nm and 850nm). Second mode with R+ (660nm) or the third NIR+ (850nm).

### SMART TIMER

Infrared lamp Aceso features a smart timer that switches off the device after 5, 10, 15, 20, 25 or 30 minutes.

### RECOVERY+

Pulsating light output of 10Hz, 20Hz or 30Hz. We recommend using eye protection with Recovery+.

### BRIGHTNESS

Adjust the brightness of the Infrared Lamp Aceso to a strength of 25%, 50%, 75% or 100%.

## OPERATING INSTRUCTIONS

- STEP 1: Connect the power cable to the device.
- STEP 2: Plug the adapter into a wall socket.
- STEP 3: Use the smart timer to select the duration.
- STEP 4: Choose Recovery+ (optional).
- STEP 5: Choose your programme: R/NIR+, R+ or NIR+.
- STEP 6: Choose the light intensity: 25%, 50%, 75% or 100%.
- STEP 7: Position yourself at a distance of 20-30cm.

Make sure your skin is exposed for best results. In addition, we recommend using Aceso for about 15-20 minutes per session.

Results may vary depending on goals and consistency of treatment. Consistent daily use for 8-12 weeks can yield systemic benefits such as improved skin health, increased energy, hormonal health and sleep quality. Daily use of the device maximises response and healing rates.

## WARNING

Always consult your doctor before using the product in case of medical conditions. Also contact your doctor if you are sensitive to light or take certain medications, are pregnant, have cancerous lesions or have recently used steroids before using the device.

This product is not intended to diagnose, cure or prevent diseases and therefore no guarantee can be made about the ability of this product to improve physical, skin or mental condition.

Stop use immediately and contact your healthcare provider if you experience discomfort or have any concerns.

**Are you content? Then share your experience through a review!**

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**Any questions? Send us an e-mail. We will be happy to help you!**

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