**DE Tageslichtlampe LT 460**

**GB Therapy Light LT 460**

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**IT Lampada a luce diurna LT 460**

**ES Lámpara de luz diurna LT 460**

**PT Lâmpada de luz diurna LT 460**

**NL Daglichtlamp LT 460**

**FI Kirkasvalovalaisin LT 460**

**SE Dagsljuslampa LT 460**

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**Zero-Plus International Limited**  
Room 1004, 10/F., Join-In Hang Sing  
Centre, 71-75 Container Port Road, Kwai  
Chung, New Territories, Hong Kong

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**EC**

**REP**

I-Plus Company Ltd.  
4B Christchurch House, Beaufort  
Court, Medway City Estate, Ro-  
chester, Kent. ME2 4FZ, UK

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## IMPORTANT INFORMATION! RETAIN FOR FUTURE USE!

Read the instruction manual carefully before using this device, especially the safety instructions, and keep the instruction manual for future use. Should you give this device to another person, it is vital that you also pass on these instructions for use.

### *Explanation of symbols*

This instruction manual belongs to this device. It contains important information about starting up and operation. Read the instruction manual thoroughly. Non-observance of these instructions can result in serious injury or damage to the device.

### **WARNING**

These warning notes must be observed to prevent any injury to the user.

### **CAUTION**

These notes must be observed to prevent any damage to the device.

### **NOTE**

These notes give you useful additional information on the installation and operation.

**IP20**

Ingress protection  
rating

Class II equipment

**LOT**

LOT number

Manufacturer

Date of  
Manufacture

Keep  
dry

## **Intended use**

The **MEDISANA** Therapy Light **LT 460** is intended for use as a light bath respectively therapy light to treat seasonal affective disorder (SAD), also known as the “winter blues”.

The use may support energy boost and stabilization of cardiac frequency as may help to experience a more “awake” feeling.

The device is intended only for the treatment of humans.

## **Contraindications**

Before using the appliance, consult your doctor if you:

- suffer(ed) from severe depression,
- had a recent eye surgery or a diagnosed eye condition, which cause avoidance of bright light
- are on specific medication (i.e. analgesics, antidepressants, psychotropic drugs or malaria tablets);
- high blood pressure / vitreous haemorrhage
- suffer from skin sensitivity to light, such as systemic lupus erythematosus

## **Possible side effects**

You may experience the following side effects:

- headaches
- eyestrain
- sleeping problems
- nausea

Measures to avoid, minimize or alleviate these side effects:

- Increase the treatment distance and do not look into the light directly
- Avoid using the appliance during the late evening or before sleep
- Stop using the appliance for a few days to make the side effects disappear and then try again.

## **Safety information**

- The appliance needs special precaution regarding EMC (electromagnetic compatibility) and needs to be installed and put into service according to the EMC information provided in the chapter “Normative reference” (see appendix).
- Use the device only according to its intended purpose as specified in the instruction manual. The warranty will be invalidated if the device is used for purposes other than those for which it is intended.
- The device is not determined for commercial use, but only for domestic use.
- The device needs no calibration, inspection or maintenance.
- Before use make sure, that all packaging material has been removed and that the device or parts of the device are not showing any visible damage. If in doubt, contact the customer service.
- Do not use the device outside or in damp rooms. Do not put the device in water or other fluids.
- Never take hold of the device with wet hands, if it is connected to the mains. Never let any humidity get on the device. Only use it if completely dry.
- Do not cover the appliance during operation or while it is cooling down.
- Only use the device in dry locations, where a solid stand is given and where the unit cannot fall. There is a danger of hurts or fractures to anyone or anything below if it falls.
- Do not operate the device in the vicinity of anaesthetic mixtures, that are inflammable in air or oxygen.
- Never touch the appliance and another person simultaneously.
- Protect the appliance from strong impacts.

- Before connecting the device to your power supply, please ensure that the supply voltage stated on the rating plate is compatible with your mains supply.
- The mains cable of the device need to be laid in that way, that the risk of stumbling or strangulation is completely excluded.
- Keep the cord away from hot objects and open flames.
- Always disconnect the mains plug and let the device cool down before you touch it!
- Always pull on the plug itself when disconnecting the device from the power supply. Never pull on the power cord.
- To completely disconnect the unit from the power source unplug the adaptor plug from the power source.
- If the mains cable, the adaptor or the device is damaged, contact the customer service centre.
- The adaptor contains a transformer. Do not cut off the adaptor to replace it with another plug, as this causes a hazardous situation.
- Do not attempt to repair the device yourself in the event of a malfunction or damages, as this would invalidate the warranty. Ask your specialized dealer. Repairs should only be carried out by authorised service centres to avoid hazards.
- Never use a damaged device.
- This appliance does not contain parts, which can be exchanged or repaired by the user.
- The appliance is protected against unauthorized use. Do not modify this appliance without authorization of the manufacturer.
- If this appliance is modified, appropriate inspection and testing must be conducted to ensure continued safe use of the equipment.
- This device is not designed to be used by persons (including children) with limited physical, sensory or mental abilities, or by persons with insufficient experience and/or knowledge, unless under observation by a person responsible for their safety, or unless they have been instructed in the use of the device.

- Children must be supervised to ensure that they do not play with the device.
- Do not leave the device unattended when it is in operation.
- Only store the device in locations, which are out of reach of children.
- If you have health concerns of any kind, contact your general practitioner.
- People with eye- or retinal diseases, with optic nerve impairments as well as diabetics or people with similar health problems need to check with their doctor beforehand, if the use of this device is uncritical for them.

**Thank you very much**

for your confidence in us and congratulations!

With the **Therapy Light LT 460** you have purchased a top quality product from **MEDISANA**. In order to achieve the desired effect with your Therapy Light **LT 460** from **MEDISANA** in the long term, we recommend that you read the following information on its use and main-tenance carefully.

**2.1 Scope of delivery and packaging**

Please check first of all that the device is complete and is not damaged in any way. In case of doubt, do not use the appliance and contact your supplier or your service centre.

The following parts are included:

- 1 **MEDISANA** Therapy Light **LT 460**
- 1 mains adaptor
- 1 stand
- 1 instruction manual
- 1 storage bag

The packaging can be reused or recycled. Please dispose properly of any packaging material no longer required. If you notice any transport damage during unpacking, please contact your supplier without delay.

### **WARNING**

**Please ensure that the polythene packing is kept away from the reach of children! Risk of suffocation!**

**2.2 How does the therapy light LT 460 work?**

In autumn and winter - the seasons with less sunshine hours and short days - the lacking sunlight influences the chemical processes in the brain. The hypothalamus function will be constrained and it is presumed, that the production of the hormone melatonin is being influenced. It may cause unwanted tiredness.

The **MEDISANA** Therapy Light **LT 460** provides bright electrical light in order to simulate the missing sunlight in autumn and winter. This is generated without the harmful UV radiation of the sunlight. The light therapy helps to level the production of the hormone melatonin and relieve the symptoms of winter depression (SAD).

### 3.1 Commissioning

- Remove the stand ⑤ from the stand storage cavity ④.
- Plug the stand ⑤ in one of the holes on the back of the device. Either use one of the holes ③ to put up the device vertically. Or plug the stand as it is shown in the picture (see „*Device and controls*“) to put up the device horizontally.
- Place the device approx. 20 to 30 cm away from your eyes. The light should reach your face sideways.
- Place the device on a stable, level surface (e.g. on a table). Make sure that the device will not drop or fall down.
- Connect the power plug to the socket for power plug ⑥ on the back of the device and to a mains socket.

#### NOTES

- **Set up the device in a place near a mains socket.**

### 3.2 Application

1. To switch on the device, press the On-/Off-button ② once. The lamps will emit a light with *high intensity*.
2. Press the On-/Off-button ② again to lower the intensity of the light - *low intensity*,
3. Press the On-/Off-button ② again to switch off the device.

As long as you enjoy the soothing light of the device, you may read a book, phone or do your work. Always keep in mind to keep the correct distance between the device and your eyes (= 20 to 30 cm).

#### Recommended use regarding the duration of application:

| High intensity |                      |                     |
|----------------|----------------------|---------------------|
| Distance       | Application duration | Application area    |
| 30 cm          | approx. 2 hours      | 400 cm <sup>2</sup> |
| 20 cm          | approx. 1 hour       | 370 cm <sup>2</sup> |

| Low intensity |                      |                     |
|---------------|----------------------|---------------------|
| Distance      | Application duration | Application area    |
| 20 cm         | approx. 2 hours      | 370 cm <sup>2</sup> |

- The closer you are to the panel, the shorter you will have to use the device to achieve the desired effect.
- Do not look directly into the light or the light screen ❶ for a longer period of time.
- During the first applications of the device, you may experience eye fatigue or headache. This is normal and should soon pass. If the symptoms persist, lower the intensity and/or the duration of the application, until you are feeling well again.

#### **Application over longer periods**

Repeat the treatment for at least 7 consecutive days during the months with less sunlight. Depending on your individual requirement, you may extend this period.

## **4 Miscellaneous**

### **4.1 Cleaning and main- tenance**

- Before cleaning the appliance, switch it off, pull the power plug out of the socket and let the device cool down completely.
- Never touch the device with wet hands.
- Clean the appliance using a soft, lint-free cloth, if necessary dampened with mild, warm soapy water and well wrung out. Do not allow water to penetrate the appliance.
- Never immerse the device in water or other liquids.
- When cleaning the device, never use aggressive cleaning agents, scouring agents or solvents; neither should you use scouring pads, steel wool or tough brushes. This may damage the surfaces.

#### 4.2 Storage

- If you do not use the device, keep it in a dry place out of reach of children.
- Let the device cool down completely before storing it.
- If you do not use the device, disconnect it from the mains supply.
- Pay attention to the storage conditions in the chapter *4.5 Technical specifications*.

#### 4.3 Disposal

This product must not be disposed together with the domestic waste. All users are obliged to hand in all electrical or electronic devices, regardless of whether or not they contain toxic substances, at a municipal or commercial collection point so that they can be disposed of in an environmentally acceptable manner. Consult your municipal authority or your dealer for information about disposal.

#### 4.4 Trouble-shooting

| Problem          | Possible cause          | Solution                                                                      |
|------------------|-------------------------|-------------------------------------------------------------------------------|
| No light         | Device is switched off. | Use On-/Off button <b>2</b> to switch it on.                                  |
|                  | No power connection     | Check mains adaptor and power cable.                                          |
|                  | No power connection     | The power cord may be defective. Contact your retailer or the service center. |
| No or weak light | LEDs at end of lifespan | Contact your retailer or the service center.                                  |
|                  | LEDs defective          | Contact your retailer or the service center.                                  |

|                                                   |                                                |                                                                                          |
|---------------------------------------------------|------------------------------------------------|------------------------------------------------------------------------------------------|
| <b>4.5<br/>Technical<br/>specifi-<br/>cations</b> | Name and model                                 | : <b>MEDISANA Therapy Light LT 460</b>                                                   |
|                                                   | Power supply                                   | : Input: 100-240 V~ 50/60 Hz<br>Output: 12 V, 1 A                                        |
|                                                   | Adapter model                                  | : CS12F120100FGF / CS12F120100FBF                                                        |
|                                                   | Power                                          | : 12 W                                                                                   |
|                                                   | Lighting elements                              | : LED                                                                                    |
|                                                   | Illuminance (Distance)                         | : > 10,000 lux (10 cm); approx. 5.000 lux (20 cm)                                        |
|                                                   | Dimensions                                     | : 26 x 15 x 1,7 cm                                                                       |
|                                                   | Weight                                         | : 550 g                                                                                  |
|                                                   | Operating conditions                           | : 0 °C to 40 °C, 30 to 90 % relative humidity,<br>atmospheric pressure: 700-1060hPa      |
|                                                   | Storage and transport conditions               | : - 10 °C to + 60 °C, 30 to 90 % relative humidity,<br>atmospheric pressure: 700-1060hPa |
|                                                   | Maximum output of radiation / spectral radiant | : 14 Wm <sup>-2</sup> sr <sup>-1</sup>                                                   |
|                                                   | Article No.                                    | : 45220                                                                                  |
|                                                   | EAN number                                     | : 40 15588 45220 8                                                                       |

**CE 0123**

**In accordance with our policy of continual product improvement, we reserve the right for technical and design changes without notice.**

The current version of this instruction manual can be found under  
[www.medisana.com](http://www.medisana.com)

**4.6  
Guidelines /  
Standards** This device is certified in accordance with EC Guidelines and carries the CE symbol (conformity symbol) "CE 0123". The specifications of EU Guideline "93/42/EEC" for medical products are met.

**Electromagnetic compatibility:** Please check the appendix for details on this measurement data.

### Warranty and repair terms

Please contact your supplier or the service centre in case of a claim under the warranty. If you have to return the unit, please enclose a copy of your receipt and state what the defect is.

The following warranty terms apply:

1. The warranty period for **MEDISANA** products is three years from date of purchase. In case of a warranty claim, the date of purchase has to be proven by means of the sales receipt or invoice.
2. Defects in material or workmanship will be corrected free of charge within the warranty period.
3. Repairs under warranty do not extend the warranty period either for the unit or for the replacement parts.
4. The following is excluded under the warranty:
  - a. All damage which has arisen due to improper treatment, e.g. non-observance of the user instructions.
  - b. All damage which is due to repairs or tampering by the customer or unauthorised third parties.
  - c. Damage which has arisen during transport from the manufacturer to the consumer or during transport to the service centre.
  - d. Accessories which are subject to normal wear and tear.
5. Liability for direct or indirect consequential losses caused by the unit are excluded even if the damage to the unit is accepted as a warranty claim.

#### **Zero-Plus International Limited**

Room 1004, 10/F., Join-In Hang Sing  
Centre, 71-75 Container Port Road, Kwai  
Chung, New Territories, Hong Kong

imported & distributed by  
**MEDISANA AG**  
Jagenbergstraße 19  
41468 NEUSS  
GERMANY

|           |            |
|-----------|------------|
| <b>EC</b> | <b>REP</b> |
|-----------|------------|

I-Plus Company Ltd.  
4B Christchurch House, Beaufort  
Court, Medway City Estate,  
Rochester, Kent. ME2 4FZ, UK

The service centre address is shown on the attached leaflet.

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| Guidance and manufacturer's declaration – electromagnetic immunity                                                                                                                                                       |                                                                                                                                                                                                                                   |                                                                                                                                                                                                                                   |                                                                                                                                                                                                                                                                                                                      |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| The model LT 460 (LND018915V2 ) is intended for use in the electromagnetic environment specified below. The customer or the user of the model LT 460 (LND018915V2) should ensure that it is used in such an environment. |                                                                                                                                                                                                                                   |                                                                                                                                                                                                                                   |                                                                                                                                                                                                                                                                                                                      |
| Immunity test                                                                                                                                                                                                            | IEC 60601 test level                                                                                                                                                                                                              | Compliance level                                                                                                                                                                                                                  | Electromagnetic environment – guidance                                                                                                                                                                                                                                                                               |
| Electrostatic discharge (ESD)<br>IEC 61000-4-2                                                                                                                                                                           | ±6 kV contact<br>±8 kV air                                                                                                                                                                                                        | ±6 kV contact<br>±8 kV air                                                                                                                                                                                                        | Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.                                                                                                                                                                       |
| Electrical fast transient/burst<br>IEC 61000-4-4                                                                                                                                                                         | ±2 kV for power supply lines                                                                                                                                                                                                      | ±2 kV for power supply lines                                                                                                                                                                                                      | Mains power quality should be that of a typical commercial or hospital environment.                                                                                                                                                                                                                                  |
| Surge<br>IEC 61000-4-5                                                                                                                                                                                                   | ±1 kV differential mode                                                                                                                                                                                                           | ±1 kV differential mode                                                                                                                                                                                                           | Mains power quality should be that of a typical commercial or hospital environment.                                                                                                                                                                                                                                  |
| Voltage dips, short interruptions and voltage variations on power supply input lines<br>IEC 61000-4-11                                                                                                                   | <5 % $U_T$<br>(>95 % dip in $U_T$ )<br>for 0,5 cycle<br><br>40 % $U_T$<br>(60 % dip in $U_T$ )<br>for 5 cycles<br><br>70 % $U_T$<br>(30 % dip in $U_T$ )<br>for 25 cycles<br><br><5 % $U_T$<br>(>95 % dip in $U_T$ )<br>for 5 sec | <5 % $U_T$<br>(>95 % dip in $U_T$ )<br>for 0,5 cycle<br><br>40 % $U_T$<br>(60 % dip in $U_T$ )<br>for 5 cycles<br><br>70 % $U_T$<br>(30 % dip in $U_T$ )<br>for 25 cycles<br><br><5 % $U_T$<br>(>95 % dip in $U_T$ )<br>for 5 sec | Mains power quality should be that of a typical commercial or hospital environment. If the user of the model LT 460 (LND018915V2) requires continued operation during power mains interruptions, it is recommended that the model LT 460 (LND018915V2) be powered from an uninterruptible power supply or a battery. |
| Power frequency (50/60 Hz) magnetic field<br>IEC 61000-4-8                                                                                                                                                               | 3 A/m                                                                                                                                                                                                                             | 3 A/m                                                                                                                                                                                                                             | Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.                                                                                                                                                                            |
| NOTE $U_T$ is the a.c. mains voltage prior to application of the test level.                                                                                                                                             |                                                                                                                                                                                                                                   |                                                                                                                                                                                                                                   |                                                                                                                                                                                                                                                                                                                      |

| Guidance and manufacturer's declaration – electromagnetic immunity                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |                          |                  |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------|------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| The model LT 460 (LND018915V2 ) is intended for use in the electromagnetic environment specified below. The customer or the user of the model LT 460 (LND018915V2) should ensure that it is used in such an environment.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |                          |                  |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| Immunity test                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    | IEC 60601 test level     | Compliance level | Electromagnetic environment – guidance                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| Conducted RF IEC 61000-4-6                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       | 3 Vrms 150 kHz to 80 MHz | 3 Vrms           | Portable and mobile RF communications equipment should be used no closer to any part of the model LND018915V2, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.<br><br>Recommended separation distance                                                                                                                                                                                                                                                                                                                                                           |
| Radiated RF IEC 61000-4-3                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        | 3 V/m 80 MHz to 2,5 GHz  | 3 V/m            | $d=1.167\sqrt{P}$<br><br>$d=1.167\sqrt{P}$ 80MHz to 800MHz<br><br>$d=2.333 \sqrt{P}$ 800MHz to 2.5 GHz<br>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, <sup>a</sup> should be less than the compliance level in each frequency range. <sup>b</sup> Interference may occur in the vicinity of equipment marked with the symbol above.  |
| NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.<br>NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |                          |                  |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| <sup>a</sup> Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the model LT 460 (LND018915V2) is used exceeds the applicable RF compliance level above, the model LT 460 (LND018915V2) should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the model LT 460 (LND018915V2). |                          |                  |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| <sup>b</sup> Over the frequency range 150 kHz to 80 MHz, field strengths should be less than [V1] V/m.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |                          |                  |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |

| Guidance and manufacturer's declaration – electromagnetic immunity                                                                                                                                                       |            |                                                                                                                                                                                                     |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| The model LT 460 (LND018915V2 ) is intended for use in the electromagnetic environment specified below. The customer or the user of the model LT 460 (LND018915V2) should ensure that it is used in such an environment. |            |                                                                                                                                                                                                     |
| Emissions test                                                                                                                                                                                                           | Compliance | Electromagnetic environment – guidance                                                                                                                                                              |
| RF emissions<br>CISPR 11                                                                                                                                                                                                 | Group 1    | The model LT 460 (LND018915V2) uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. |
| RF emissions<br>CISPR 11                                                                                                                                                                                                 | Class B    |                                                                                                                                                                                                     |
| Harmonic emissions<br>IEC 61000-3-2                                                                                                                                                                                      | Class C    |                                                                                                                                                                                                     |
| Voltage fluctuations / flicker emissions<br>IEC 61000-3-3                                                                                                                                                                | complies   |                                                                                                                                                                                                     |

| Recommended separation distances between portable and mobile RF communications equipment and the LT 460 (LND018915V2)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |                                                             |                                               |                                             |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------|-----------------------------------------------|---------------------------------------------|
| The model LT 460 (LND018915V2) is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Model LT 460 (LND018915V2) can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the model LT 460 (LND018915V2) as recommended below, according to the maximum output power of the communications equipment.                                                                                                                        |                                                             |                                               |                                             |
| Rated maximum output power of transmitter W                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     | Separation distance according to frequency of transmitter m |                                               |                                             |
|                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 | 150 kHz to 80 MHz                                           | 80 MHz to 800 MHz                             | 800 MHz to 2,5 GHz                          |
|                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 | $d = \left[ \frac{3,5}{V_1} \right] \sqrt{P}$               | $d = \left[ \frac{3,5}{E_1} \right] \sqrt{P}$ | $d = \left[ \frac{7}{E_1} \right] \sqrt{P}$ |
| 0.01                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            | 0.117                                                       | 0.117                                         | 0.233                                       |
| 0.1                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             | 0.369                                                       | 0.369                                         | 0.738                                       |
| 1                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               | 1.167                                                       | 1.167                                         | 2.333                                       |
| 10                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              | 3.690                                                       | 3.690                                         | 7.377                                       |
| 100                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             | 11.67                                                       | 11.67                                         | 23.33                                       |
| For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer. <b>NOTE 1</b> At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies. <b>NOTE 2</b> These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people. |                                                             |                                               |                                             |

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