## 3-in-1 Ultrasonic Cavitation Massager InnovaGoods® - Wellness Beauty

### 3-in-1 Ultrasonic Cavitation Massager InnovaGoods®

### 3-in-1 Ultrasonic Cavitation Massager InnovaGoods®

## InnovaGoods®

### 3-in-1 Ultrasonic Cavitation Massager InnovaGoods®

### 3-in-1 Ultrasonic Cavitation Massager InnovaGoods®

### 3-in-1 Ultrasonic Cavitation Massager InnovaGoods®

### User's Manual

### 3-in-1 Ultrasonic Cavitation Massager InnovaGoods®

### 3-in-1 Ultrasonic Cavitation Massager InnovaGoods®

### 3-in-1 Ultrasonic Cavitation Massager InnovaGoods®

## ENGLISH

### 3-in-1 Ultrasonic Cavitation Massager InnovaGoods®

- Charger: 110-240 V / 50-60 Hz
- Input: DC 15 V / 800 mAh
- 3 massage functions

### 3-in-1 Ultrasonic Cavitation Massager InnovaGoods®

### 3-in-1 Ultrasonic Cavitation Massager InnovaGoods®

### 3-in-1 Ultrasonic Cavitation Massager InnovaGoods®

### 3-in-1 Ultrasonic Cavitation Massager InnovaGoods®

- Clean the skin and apply moisturising cream/gel on it evenly; keep the skin constantly hydrated; do not use directly on dry skin.
- Plug the 3-in-1 Ultrasonic Cavitation Massager InnovaGoods® into the mains electricity and press the power button.
- Ultrasound function: press the «SONIC» button to activate the function and change the intensity (high: H – low: L).
- Infrared function: press the «INFRARED» button to activate the infrared heat function; the red light will illuminate on the head.
- Electrostimulation function: clean the skin (you can apply a conducting gel, not included); connect the cable to the device and the patches; stick the patches to the skin; press the «EMS» button and choose from the modes: Tapping, Massage, Knead, Scraping or Slimming; select the desired intensity.
- Use in sessions of between 5 and 15 mins; do not exceed 30 mins each day.
- In the event of overheating, the Ultrasonic Cavitation Massager InnovaGoods® will switch off automatically.
- Products with a gel texture are recommended for skincare since they are better conductors.

### 3-in-1 Ultrasonic Cavitation Massager InnovaGoods®

### 3-in-1 Ultrasonic Cavitation Massager InnovaGoods®

### 3-in-1 Ultrasonic Cavitation Massager InnovaGoods®

Make sure that the voltage of the power source corresponds to that of the appliance. Unplug the appliance whenever you are not using it. Do not pull on the cable and ensure it does not tangle. Keep the cable away from hot and sharp surfaces. Before cleaning, make sure that appliance is disconnected from the power supply and has cooled completely. Clean with a damp cloth. Do not use abrasive products for cleaning. Do not immerse the appliance or cable in liquid. Do not handle it with wet hands Do not use near water or

damp areas. Suitable for domestic use only. Use the device for its intended purpose only. Store the device in a dry, dust-free place. Do not leave the appliance unsupervised when it is turned on. This appliance is not a toy. Keep out of reach of children. Do not modify the appliance or any of its components. If the appliance, the cable or the pieces looked damaged or do not work correctly, do not use them or attempt to repair them; they should only be repaired by a qualified technician. This device is not intended for use by children or persons with reduced physical, sensory or mental capabilities, unless they are supervised by an adult responsible for their safety. This device is not household waste. Please follow the environmental regulations in place for its recycling Discard the product in accordance with the European Directive 2002/96/EC. Enquire about your local recycling legislation for electrical and electronic products and follow these steps. Do not dispose of this product with household waste. Following these instructions will help prevent harmful consequences to personal and environmental safety.

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