KÖNIG
HEALTH CARE

HC-DT10

MANUAL (p. 2)
DIGITAL THERMOMETER

ANLEITUNG (p. 4)
DIGITALTHERMOMETER

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ENGLISH

How to use:

1. Before using, please disinfect the probe with alcohol first.
2. Press the on/off button, Bi sound will be heard and display „88.88°C for about 2 seconds.
3. After that, the thermometer will display the last temperature memory record for about 2 seconds, and then display L °C (Lo°C), signal °C flashing, that means the thermometer is waiting for use. If the temperature measured is below or above the display range the thermometer will display either the „L“ or the „H“ symbol. The display range for this thermometer is 32.00°C to 43.00°C.
4. Use the thermometer to measure body temperature. The value of temperature displayed will arise gradually and °C is flashing.
5. The measurement process usually takes about 1 min (in oral). After measurement, °C will stop flashing and the alarm signal will sound for about 5 seconds, at the same time, the thermometer will display the user's body temperature.
6. This thermometer can automatically shut off within 10 minutes after measurement. It is recommended to turn off the power manually right after measurement in order to save energy.

Take temperature:

When you measure your body temperature, please remain to measure at the same place for some time.

1. To take temperature under tongue, before measurement please close your mouth for 5 minutes, then insert the sensor part under the tongue, keeping nestled and close your mouth. LCD reading is final reading.
2. To take temperature in the armpit, before measurement please completely wipe the sweat there away, then close the armpit for about 5 minutes and insert the sensor part against the center of the armpit closely. Put both left and right waist intersection, making the thermometer in contact with the user's body closely. Add 0.5 °C to LCD reading to obtain final reading.
3. To take temperature in the rectum, inserts the sensor part into anus gently. Minus 0.5 °C to LCD reading to obtain final reading.

Cautions:

1. Do not bend or drop the thermometer.
2. Please disinfect the probe with alcohol, do not boil the probe.
3. Do not put the thermometer in contact with chemical reagent or diluents.
4. Avoid exposing the thermometer under direct sunlight.
5. Don't disassemble the thermometer except replacing the battery.
6. This thermometer has already been calibrated while making. If use according to the operation instruction, does not need to calibrate again.

Battery replacement:

1. The battery needs to be replaced by a new battery when the LCD displays  or .
2. Please open the battery cover, take out the old battery with a toothpick.
3. Place a new 1.5V button battery (LR/SR-41).
4. Replace the battery cover.

: Degree of protection against electric shock is Type B applied part.

: The disposal of battery and device shall comply with the local environment requirements.

The device comply with the standard of IEC 60601-1:1998+ A1:1991+ A2:1995 and IEC60601-1-2:2001

Declaration of conformity

We,

Nedis B.V.

De Tweeling 28

5215MC 's-Hertogenbosch

The Netherlands

Tel.: 0031 73 599 1055

Email: info@nedis.com

Declare that product:

Name: KONIG

Model: HC-DT10

Description: Digital thermometers

Is in conformity with the following standards

EMC: EN 60601-1-2:2001

LVD: EN 60601-1:1998+A1:1991+A2:1995

Following the provisions of the 93/42/EEC directive.

's-Hertogenbosch, 2007, 04 December

Mrs. J. Gilad

Purchase Director

Warranty:

No guarantee or liability can be accepted for any changes and modifications of the product or damage caused due to incorrect use of this product.

General:

Designs and specifications are subject to change without notice.

All logos brands and product names are trademarks or registered trademarks of their respective holders and are hereby recognized as such.

Attention:

This product is marked with this symbol. It means that used electrical and electronic products should not be mixed with general household waste. There is a separate collections system for these products.

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Copyright ©