


# Myofascial Physiotherapy Device User Manual



- 4. Configuration**
- Lower on/off button
  - Device level indicator (red LED) (off)
  - High speed level indicator (red LED) (off)
  - Removable, detachable silver metal head when removing the vibration head, the silver metal head will be out.
  - DC charging pin
  - AC/DC charger
  - 7.4V machine with 8.4V1ade, 12V machine with 12.6V1ade, USB with data cable (customer provided 5V2A charging head)

**Maintenance, cleaning and storage:**  
Clean the device with a slightly wet towel, and then wipe it dry with a soft cloth for storage. Do not put it into the water and enter when you are traveling.

**Warnings and reminders**  
Do not immerse it into water. Keep it away from liquids or heat sources. Keep the vent free from dust and scraps.  
Do not remove the screws or try to disassemble the device.  
Do not continuously use it for more than one hour. Please stop the device for 30 minutes and then use it again.  
Unplug the charger from the battery after charging or before the use.  
It is only for adults. If you have any use, do not use this device. Please consult with your doctor before using this device.

**Charging**  
Before using it for the first time, charge it fully for 4-6 hours.  
When using, connect the device to the charger provided in the charging part of the battery, and then plug the AC side cable into the socket.  
The LED on the battery flashes regularly and displays the battery level, indicating that the battery is being charged.  
When the battery is fully charged, the LED on the device flashes and displays the battery level, indicating that the battery is being charged.  
When all LEDs stop flashing and one of them indicates that the battery is fully charged.  
The battery can be charged at any time and any level.  
It is not recommended to use up the battery to reach LED off.  
The average operating time of the device is two hours, depending on the speed and the applied pressure during the use.  
Do not use the machine when it is charged.

**Safety Instructions**  
**Danger: In order to reduce the risk of electric shock, fire and personal injury, this product must be used in accordance with the following instructions.**


This is only for adults only.  
It can only be used on dry, clean skin surfaces of the body, and cannot be used through clothes. Press and move it gently on the skin. The duration on each part is about 60 seconds.  
Use this device only on soft tissues of the body without causing pain or discomfort. Do not use it on head or any hard or bony part of the body.  
Only use the vibrating head that can produce the best results.  
Bruise may be caused whether under moderate pressure or normal pressure. Please always check the area that receives physiotherapy.  
Stop using it immediately once there is any sign of pain or discomfort.  
Keep fingers, hair or other parts of the body away from the shaft and back of the vibrating head, because they might be pinched.  
Do not place any objects in the vent of the device.  
Do not immerse it in water or allow water to enter the vent of the device.  
Do not drop the device or use it improperly.  
Carefully check the device and battery before each use.  
Do not change or modify the device in any way.  
Never operate or change the device in an unattended way.



**Front position**  
Shoulder girdle: 20-30cm  
Cervical vertebrae: 15-20cm  
Thoracic vertebrae: 15-20cm  
Lumbar vertebrae: 15-20cm  
Sacrum: 15-20cm  
Pelvis: 15-20cm  
Anatomical diagram of the spine and ribcage.




**Back position**  
Cervical vertebrae: 15-20cm  
Thoracic vertebrae: 15-20cm  
Lumbar vertebrae: 15-20cm  
Sacrum: 15-20cm  
Pelvis: 15-20cm  
Anatomical diagram of the back and neck.



**Lower limb position**  
Hip: 15-20cm  
Knee: 15-20cm  
Ankle: 15-20cm  
Foot: 15-20cm  
Anatomical diagram of the lower limb.

Indications and Contraindications	Indications	Contraindications
<p>The deep penetration of the device can reduce pain from joint sprains, bursitis, muscle cramps, neuritis and other musculoskeletal disorders in a short period of time. Almost all soft tissue diseases are associated with trauma, excessive muscle use or abnormal posture, and most of the scar tissues involved in these diseases can be treated with this device. Avoid excessive stress and prolonged physiotherapy when using the device. Also avoid sensitive areas such as head, face and cervical vertebrae. Stop physiotherapy if there is inflammation, swelling or increased pain.</p> <p>This is very likely to reduce the types of contraindications mentioned above, and practical experience has shown that in many cases, it is also a physiotherapy recommendation to combine vibration training into physiotherapy programs. This must be done by, or in the company of, a doctor, expert or physical therapist.</p>	<ul style="list-style-type: none"> <li>&gt; Pain and cramps due to muscle injury.</li> <li>&gt; Aneurysm, bleeding, use of blood thinner</li> <li>&gt; Heart disease, with a pacemaker or defibrillator</li> <li>&gt; Help the flow of edematous fluid in the swollen area</li> <li>&gt; Relax thickened connective tissues and area</li> <li>&gt; Within 90 days of intra-articular fixation, the device should not be used within 3 inches of the incision.</li> <li>&gt; Reduce the accumulation of lactic acid in muscles</li> <li>&gt; Increase joint mobility</li> <li>&gt; Eliminate muscle fatigue</li> </ul>	<ul style="list-style-type: none"> <li>&gt; Sensitive area: head, face, cervical vertebrae, spine</li> <li>&gt; Positions close to the bones with less muscle coverage: the back of the foot, the back of the hand</li> </ul>



**Foot position**  
Anatomical diagram of the foot.

Use Method

Use Method