

	100*80mm
	150g

## WEARABLE HEATING PAD USER MANUAL

3 VIBRATION  
MODES



CORDLESS  
HEATING PAD  
FOR PERIOD PAIN

### Product Parameters

<b>Name:</b>	Wearable Cordless Heating Pad
<b>Size:</b>	200*92*25mm
<b>Use Scope:</b>	Waist Abdomen
<b>Warm Gear:</b>	50/55/60°C
<b>Color:</b>	White / Pink / Black
<b>Charging Ways:</b>	Type-C Charge
<b>Material:</b>	ABS+Deer Velvet Fabrics
<b>Charging Voltage:</b>	5V

### Functional Features

Machines using graphene heating principle. Three temperatures are adjustable. 122°F (50°C) for daily warming, 131°F (55°C) for period nursing, 140°F (60°C) for cramps and pain. Also can with the vibration massage.

1



- **Vibration Switch:**  
Long press for 2s to turn on the vibration switch, short press to switch mode, long press for 2s to turn off.
- **Warm Switch:**  
Long press for 2s to open the temperature switch, short press to switch mode, long press for 2s to close.

2

### Safety Instructions

- 1 Please do not use it while sleeping.
- 2 The ban on near the fire, and all kinds of corrosive liquid used in an adverse environment.
- 3 Please avoid using it in the same position for a long time. Attention to confirm the skin condition, if there is any abnormal immediately remove it.

### Disabled People

- 1 Pregnant women, cervical tumor patients disabled.
- 2 Gout onset period or not completely healed with gout are forbidden to use.
- 3 Gravity disease and heart disease are forbidden to use.
- 4 The skin wound is still unhealed or patients with skin abnormalities are forbidden to use.



3