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User Manual

Product introduction

The instrument uses photon beauty IPL technology to irradiate the skin directly with intense pulsed light to produce photochemical effect. On the one hand, it can regenerate collagen fibers and elastic fibers in the skin, restore skin elasticity, improve blood microcirculation, and remove or reduce wrinkles; on the other hand, intense pulsed light can penetrate the skin, and be selectively absorbed by pigment groups and capillaries in tissues, heat energy can coagulate the blood without damaging normal skin tissues, and the pigment molecules is heated so as to eliminate the expansion of capillaries, remove color spots, and remove excess hair.

Product summary

Function Introduction

1. The display effect is 999999 times.
2. Insert the adapter in the default display screen flashes three times and goes off. Press and hold the power on key for 2 seconds to start the machine. The refrigeration is turned on by default, and the refrigeration indicator is also on, and the gear 1 indicator is on for a long time.
3. Insert the adapter in the default display screen flashes three times and goes off. Long press the power on key 2-3s turn on and turn off.
4. After the machine is turned on, the default gear is gear 1. Short press the power on key to enter the shift mode. Gear 1 is gear 1, gear 5 is gear 5.
5. After turning on the machine, touch the skin, press the lighting key for a short time, the light tube will flash, and long press for 1.5 seconds to turn on or off the automatic lighting mode, and the corresponding icon will be on or off, and the indicator light will be on or off.
6. After the machine is turned on, the head sensor function can detect the hair gear automatic light flashing, press the lighting button to light up, and in the automatic mode, the light will come out when touching the skin.
7. Power supply: DC12V / 3a 3C. Refrigeration output: 2.8-3.2v hair output: 12v-0.19a-6000rpm.
8. Output voltage of large capacitor: 1: 270, 2: 290, 3: 310, 4: 330, 5: 350.
9. Double-click the flash button to turn the refrigeration mode on or off.

Product Error Report

Wrong Labeling	Warning
ERR1	Cooling Wind Blows Not, Turn
ERR2	Warning High temperature at outlet of hair
ERR3	Power Input and Machine abnormality
ERR4	Insufficient power input voltage (less than 12V)

Matters Needing Attention

1. Persons with skin sensation impairment should not use it.
2. Always clean face, axilla and other skin diseases, allergic constitution can not be used.
3. Pregnant women, Simultaneous treatment period, tattoo period, can not use.
4. Those who have been exposed to sunlight 2 weeks cannot be used.
5. Dark skin color, light skin color and tattoo can not be used.

Instructions

- Preparation before use:
1. Please remove the hair on the surface of the skin and do a good job of skin care treatment, but before shaving, you must wet the hair using soap or shower gel. Shaving is not the hair in the last, otherwise "dry shaving" will easily damage the skin on the skin surface and cause allergy.
 2. Clean the skin and make sure there is no obvious hair on the surface. It is recommended to use a depilatory cream before using the photon laser removal, it is strictly forbidden to apply conditioning agents, skin oil and other products. Otherwise, there will be no light.

Product Parameters

Serial Number	Test Item	Test Item, general use test, basic performance, adjustable transformer, power factor, hair removal
1	Test Condition	10 °C – 50 °C, humidity below 80%
	The Power Adapter	Input: 100-240V-50 / 60Hz 1.2A Output: 12V-3A 3C certification
	Japan Voltage	12V
	Japan Current	3A
2	Output output power	Light mode: 36W
	Temperature protection	When the temperature is higher than 80 °C, it will stop working and restart when it is less than 50 °C.

Treatment stage cycle reference table

1-4 weeks	5-6 weeks	7-8 weeks	More than 8 weeks
More than 75% hair falls out naturally	Two times a week Loss of new hair and reduced growth	Once a week The hair becomes thinner and less	1-2 times a week According to the state of hair growth, selective use that when it grows use it again and it no longer grows

One course of hair removal is irradiation for 3 times. After the first exposure to the depilated part, the second and third times should be done, do not irradiate the same place for 3 times continuously.

The higher the gear, the better the effect. It is recommended to use 1-2 gear in the early stage, and 5 gear can be selected after skin adaptation in the later stage. The interval time of treatment depends on the growth rate of your personal hair, and the situation of different parts of the body will also be different.

Cleaning and preservation

1. After use, please turn off the power supply of the device, unplug the plug and let the instrument cool down.
2. Wet the cleaning soft cloth with a little water, and then use it to clean the following parts:
 - Filter
 - Skin contact surface
3. Store the instrument in a dust free dry place with a temperature between 0°C and 60°C.

Product maintenance

1. Please turn off the power supply of the device immediately after use, unplug the plug and let the instrument cool down.
2. Wet the cleaning soft cloth with a little water, and then use it to clean the following parts:
 - Skin contact part
3. Store the instrument in a dust free dry place with a temperature between 0°C and 60°C.

Operation Precautions

1. During operation, the eyes shall not flash directly at the light mouth or directly flash to the eyes.
2. After removal of skin after operation is a normal phenomenon, ice compress or cold cooling can be used.
3. When cleaning the dirt on the glass in the light outlet, it is forbidden to spray alcohol directly on the light outlet.
4. When the instrument is in operation, no gel like cosmetics are allowed on the skin, such as gel, so as not to damage the instrument.
5. Try not to eat sensitive food (such as spicy, heat, resistant, etc.) during operation.
6. The patients with darker skin may have delayed skin reaction, and the energy should not be increased because the immediate reaction is not obvious.
7. Pay attention to the moisturizing and sunscreen of skin after shaving.
8. In the process of use, please do not touch the air inlet and cooling hole, so as not to affect the heat dissipation.
9. The product does not have waterproof function, please do not use it in bathrooms, shower and other places that may easily get water.
10. If the shell and window of the product are damaged or soaked in water, please send it back to the manufacturer for maintenance to determine whether it can be repaired and reused. If the product is immersed in water under charged condition, please unplug the power supply and contact the instrument to ensure your safety.
11. Avoid storage in places with strong magnetic field, without moisture and static or high temperature.
12. Do not strongly impact or knock the instrument.
13. It is forbidden to use the product for purposes other than the specified functions.

Use Complete

1. The skin may have slight redness, swelling, tingling or warm feeling. These reactions are completely harmless and soon disappear.
2. After using the instrument for 3 hours, you can safely apply some lotion, cream, fragrance agent, moisturizing product or cosmetics in the course of treatment. If the skin is highly sensitive after the end of the treatment, please wait for these reactions to disappear and then apply any product on the skin. If the skin is itchy after applying the product, please wash it with water.