



Mini massage stick

LCD model (8 modes and 19 gears)

1. It can be used for about 100 minutes after one hour of charging. The greater the strength, the greater the power consumption. If it cannot be used during charging, it can be used after charging. 2. The bottom button of the foot is button on the right side of the stick. Please power on and press the bottom button to increase the force. Decrease the force when it is to switch the mode.

Heat usage

The temperature intensity is divided into 10 gears (increased by 0.10 gears). When using for the first time, please do not set a lower strength. Please adjust the strength after a heating to low strength. When the circulating power of the battery gradually decreases, the current intensity will also gradually weaken. If you continue to use it, please charge the full.

(0-mode)

When the intensity is set to 0, it will automatically over-charge every 10 seconds.

(19 gear adjustment)

If the intensity is gear 1-10, it is 0.10 gears. If the intensity is gear 11-19, it will automatically increase by 0.10 gears. "gear 0 will automatically shut down within 10 seconds."

Please not a lower intensity when using it for the first time, and adjust after adapting to the low intensity. There are 8 modes in total. For details, please refer to the following massage mode guide.

Technical parameter

Model number	BR-B010	Color	Black/Blue
Voltage	3.7V	Size	82*62*16
How to use	4-5 (1) (gear) 19 gears	Weight	31
Frequency	1-100HZ	Material	PC+ABS
Material	ABS+PC	Power	200mA

Matters needing attention

- When charging, do not unplug the foot, just plug in the charging wire for charging.
- Two clips must be installed properly when using.

ITEM NO: BR-B010

MADE IN CHINA

