Technische Specificaties

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NorahLux Lumnia

Feel the power and know what makes the difference!

This portable Red Light Therapy device is the most affordable in its class. This makes it an ease to bring along during travels. Despite her size she is as powerful per cm² as her bigger brothers! This means that the Lumnia is applicable for all therapeutic purposes as described on our website.

Certificate

All NorahLux lights are CE and RoHS certified

Instructions

This model is plug & play. Connect the power supply to the wall socket and the other end into the backside of the light. From this point on the light is ready for use. Turn the light on by pushing the buttons RED and NIR. They can operate independently to get red light only, near-infrared light only or in combination. In most cases the combination of red and nir will give the best results. The remaining buttons give you an option to automatically turn off the light after a set time as indicated. Most times the 10 or 15 minutes button will suffice. You can find more information and answers to questions on the FAQ page on our website.

Distance and Time

With experience we gained during the years using Red Light Therapy it has become very easy to explain the best use for results. Optimal results are made with sessions at a distance of 10 to 15cm and a time of 10 to 15 minutes per treated area, for instance the face area.

Distance to the light: 10cm to 15cm

Duration of session: 10 to 15 minutes

Repeat: If possible, 1x daily

Video instructions

Instructions and other data is presented on video on our Youtube channel NorahLux Rood Licht Therapie. We will try to add English versions soon! You can go to the video by scanning the QR-code.

Safety

Contrary to most conventional treatments Red Light Therapy is safe, non-invasive, is free of chemicals, medicine and UV radiation and is not associated with harmful side effects.

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