Usage of the Red Light Therapy Device







Plug the Power in



Turn the Power On



Enjoy your Red Light Therapy

- Distance: Apply the light from 0-60 cm to the body.
- Time: Expose each area for 10-20 minutes, depending on distance, repeat on other areas if desired.
- Frequency: Do this 2 times to 14 times a week until results are felt.
- You can turn the Deep Red or NIR lights on or off at will, depending on the desired results.
- Skin rejuvenation, pain relief, hair growth, acne scars, inflammation, sports injuries, neck pain and stiffness, etc.
- Please ask your doctor if you have physiologic contraindications before treatment.



| Adress: | Markenhaven 28 |
|--------------|------------------|
| City: | Amersfoort |
| Zipcode: | 3826 AC |
| Country: | Netherlands |
| Phonenumber: | +31 85 3033 501 |
| Email: | info@phipower.nl |
| Website: | www.phipower.nl |



INSTRUCTION MANUAL RED LIGHT THERAPY

(Timer Control Design)





Technical specification

| Model | PHIRLT40 | Beam Angle | 30 degrees |
|---------------|----------------------|---------------------|------------------|
| LED Power | 200W | Shell Color | White |
| Draw Power | 60W±5% | Product Size | 247*147*40mm |
| Input Voltage | 100-240V AC,50/60 Hz | Working Temperature | –20 °C to 60 °C |
| LEDs Qty | 40pcs | Life Span | 50,000 hours |
| LED type | 3W/5W Mixed | Timer Control | 10/15/20 minutes |
| Wavelength | 660nm/850nm | | |

Deep Red + Near Infrared (NIR) Light

This type of light will go the deepest into your skin. It will help to address the clarity and tone of your skin as well as the deeper bone and muscle. Great for elastin and collagen production, this type of light helps to heal skin with wounds and scars

Near Infrared (NIR) Light - 850nm

- · Pain relief
- · Speeds up wound healing
- · Increases circulation
- · Healing of muscle and ligament
- · Boosts enzymes to improve body
- · Boosts mitochondria function in cells
- · Increases cell energy
- · Alleviates allergy symptoms
- · Speeds up healing of burns
- · Skin repair

Deep Red Light - 660nm

- Boosts liver function
- · Beneficial for lungs
- Promotes circulation Pain relief
- · Beneficial for bones
- · Beneficial for heart
- · Promotes wound healing
- Beneficial for kidneys Boosts liver function

- Skin repair
- Reduces scar tissue
- Enhances the senses
- Improves nutritional metabolism
- Regenerates blood
- Promotes muscle healing

Take it anywhere with you





