Interdental Power Flosser

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EN

Thank you for purchasing the Interdental Power Flosser. Please read the operating instructions and safety warnings carefully before using the device. Keep this manual for later consultation.

Intended use

This interdental cleaning device is exclusively for cleaning the interdental spaces with the nylon picks recommended by OptiSmile. Any other use or modification of the appliance is considered non-intended use and involves considerable risks. OptiSmile assumes no liability or responsibility for any damage arising from improper use or failure to observe these operating instructions.

Content

- 1x Interdental Power Flosser
- 1x Holder
- 60x Replacement nylon picks

Power source

- 2 x 1.5V AA batteries (not included)

Inserting/ changing batteries

1. Carefully pull off the handset.
2. Take out (if inserted) the old batteries and insert two 1.5V type AA batteries conform the polarity marking.
3. Close the battery compartment and make sure the sealing ring is in position. Replace the batteries when the effectiveness decreases.

Insert/change replacement nylon picks

1. If a nylon pick is pinned, pull the pick out.
2. Insert a new nylon pick in the device.

Switching ON/OFF

1. Turn the device ON by pressing and holding the ON/OFF button (LED lights up).
2. The device turns off when the ON/OFF button is released.

Cleaning interdental spaces

Due to the length of the nylon picks it's sufficient to only clean from the front of your teeth. Cleaning of the interdental spaces does not replace daily teeth brushing.

1. Insert the nylon pick into the interdental space. Hold the device vertically during cleaning.
2. Turn ON the device and move the device up and down along the interdental space. Continue with each interdental space.
3. When you're finished cleaning, rinse off the nylon pick. Only rinse the nylon pick since the device is not waterproof.

Cleaning & storage

- Clean the device with a damp cloth and a mild cleaning agent. Make sure no liquids enter the inferior housing.
- Dry all parts before using or storing, store the device in the holder.
- If the device isn't used for a long time, remove the batteries and store in a place where it is protected from dust, humidity and direct sunlight.

Warnings

- Before using the device, check on visible damages. DO NOT use when it has been dropped or damaged.
- Do not use force when using the device in order to avoid injuries in the oral cavity.
- Keep device and packing materials out of reach of babies and toddlers.
- The device may be used by children from 12 years old and people with reduced physical sensory or mental capabilities, lack of experience or knowledge if they have been instructed how to use the device or when they are under supervision.
- Do not submerge in water. The device is only splash proof.
- Do not expose to extreme temperatures, high humidity, heat sources, hot surfaces and open fires.
- Do not modify the device and do not open the housing.
- The replacement nylon picks should not be inhaled or swallowed.
- Only use original Optismile replacement picks. Purchase the picks on www.optismile.eu

Battery safety instructions

- Never throw batteries into a fire.
- Never short-circuit batteries.
- Do not recharge the batteries.
- Check the condition of the batteries regularly. Leaking batteries can cause damage to the device.
- Use protective gloves in case of leaking batteries. This prevents possible burns to the skin.
- Store batteries in a safe place, out of reach of children.

Disposal of the appliance

Never dispose of the appliance in your normal domestic waste. This product is subject to the provisions of the European Directive 2012/19/EU.

Dispose of the appliance via an approved waste disposal company or your municipal waste disposal facility. Please comply with all applicable regulations. Please contact your waste disposal center if you are in any doubt.

Your local community or municipal authorities can provide information on how to dispose of the worn-out product.

Disposal of the batteries

Used batteries must not be disposed of in the household waste.

Consumers are legally obliged to dispose of (rechargeable) batteries at a collection point in their community/city district or at a retail store. The purpose of this requirement is to ensure that batteries are disposed of in an environmentally friendly manner. Only dispose of batteries when they are fully discharged.

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